

CORMAY MICROALBUMIN CONTROL



CONTROL SERUM FOR DETERMINATION OF ALBUMIN CONCENTRATION IN URINE AND CEREBROSPINAL FLUID

II GENERATION

CORMAY MICROALBUMIN CONTROL 1 x 2 ml Cat. No 4-461

Control serum CORMAY MICROALBUMIN CONTROL is prepared from human serum diluted with saline. It contains preservatives to inhibit microbial growth. This control is designed for routine quality control when measuring albumin in urine and in cerebrospinal fluid.

The exact concentrations of albumin are printed on the label.

PREPARATION

1. Control serum is ready for use.
2. Before every use mix by gently inverting the vial.

STABILITY AND STORAGE

1. Tightly closed control serum is stable at 2-8°C until the expiry date printed on the package.
2. Once opened, control serum is stable for 4 weeks when stored at 2-8°C.
3. **Protect from light! Do not freeze!**

NOTES

- Product for in vitro diagnostic use only.
- Control serum must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.
- Do not use after expiry date.
- **Control serum's value is lot-specific.**
- Control serum is based on human serum that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Control serum contains sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 04. 2017.

MANUFACTURER

PZ CORMAY S.A.
ul. Wiosenna 22,
05-092 Łomianki, POLAND
tel.: +48 (0) 81 749 44 00
fax: +48 (0) 81 749 44 34
<http://www.pzcormay.pl>