PRESTIGE 24i LQ MG

Nr kat. 4-229, 4-429, 4-462 (PL)

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(EN)

Danger

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

SPECIMEN ^{6,7}

Serum, heparinized plasma free from hemolysis, 24-hours urine.

Recommended anticoagulants: heparine lithium, sodium or ammonium salt.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain approximately 3 times the magnesium concentration found in normal serum.

Urine preparation: acidify urine with some drops of concentrated hydrochloride acid to pH 1.0. Then dilute 1 part of acidified urine with 4 parts of distilled water. Multiply the result by 5. Mix well samples before analysis.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.

Avoid foaming.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- **Biolis 24i Premium:** When performing assays in the analyzer there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, it is recommended to **precede the MG test by wash cycle** with BIOLIS WASHING SOLUTION (Cat. No 3-126, 3-127) according to the guidelines contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

- **Biolis 30i:** When performing assays in the analyzer, there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of magnesium concentration should be performed **in the separate order** (follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER).

REFERENCE VALUES

serum / plasma	mg/dl	mmol/l
newborn 2 – 4 d	1.5 – 2.2	0.62 – 0.91
children 5 mo – 6 y	1.7 – 2.3	0.70 – 0.95
6 – 12 y	1.7 – 2.1	0.70 – 0.86
12 – 20 y	1.7 – 2.2	0.70 – 0.91
adults	1.6 – 2.6	0.66 – 1.07
24-hours urine	mg/24h	mmol/24h
	72.9 – 145.8	3 – 5

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum, CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every week (Prestige 24i, Biolis 24i Premium) or every 2 weeks (Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**
0.02 mg/dl (0.008 mmol/l)
- **LoD (Limit of Detection):**
0.03 mg/dl (0.012 mmol/l)
- **LoQ (Limit of Quantitation):**
0.07 mg/dl (0.029 mmol/l)
- **Linearity:**
up to 5.0 mg/dl (2.06 mmol/l)

Specificity / Interferences

Haemoglobin up to 0.313 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl, triglycerides up to 1000 mg/dl and calcium up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	2.02	0.02	1.04
level 2	4.31	0.04	0.92
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	2.0	0.08	3.7
level 2	4.1	0.11	2.7

Method comparison

A comparison between magnesium values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 59 serum samples gave following results:

$$y = 0.9927x + 0.0361 \text{ mg/dl};$$

$$R = 0.975 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	22	MG
Auto Rerun SW		
ON		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	0.14	5.0
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
		OFF

- Prestige 24i Premium, Biolis 24i Premium

Item No.	22	Item Name	MG	Optical	
Data information					
Units	mg/dl				
Decimals	2				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	END method				
Main Wave Length	505 nm				
Sub Wave Length	700 nm				
Method	Xylidyl blue				
Correlation					
Y=	Slope	Intercept			
	1	0			

Item No.	22	Item Name	MG	Optical
Aspiration				
Kind	Single			
Vol.	Kind	Vol.	Add	Units
	Sample	3	5	µl
	Reagent 1	250	10	µl
Data Process				
Read	Main	Start	End	
	Sub	9	10	
Abs.Limit				
Low	-0.1	High	2	
Correction value				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point				
1				
Span				
3				
Third mixing				
OFF				

Item No.	22	Item Name	MG	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	1.6	2.6	1.6	2.6
Urine				
Plasma	1.6	2.6	1.6	2.6
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	22	Item Name	MG	Optical	
Auto Rerun SW					
ON					
Auto Rerun Condition (Absorbance)					
	Lower	OFF			
	Higher	OFF			
Auto Rerun Range (Conc.)					
	First Dil	Low Value	Dil	High Value	Dil
Serum		0.06		5.0	
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Prozone)					
OFF					
Dilution					
100:Dil2					

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