PRESTIGE 24i HDL DIRECT

Nr kat. 4-179, 4-479 (PL)

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LITERATURA

- Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl:1 4-13.
- Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85:1234-41.
- Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
- Camps, J, Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol, Clinical Chemistry, 1999; 45:685-688

Data wydania: 04.2021.

PRESTIGE 24i HDL DIRECT

Cat. No 4-179, 4-479

(EN)

Stability on board of the analyser at 2-10°C: Prestige 24i - 12 weeks, Biolis 24i Premium - 12 weeks.

INTENDED USE

Diagnostic kit for determination of HDL-cholesterol concentration (direct method) used in the automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL). The principle role of HDL in lipid metabolism is the uptake and transport of cholesterol from peripheral tissues to the liver. Low HDL cholesterol (HDL-C) levels are strongly associated with an increased risk of coronary artery disease. The HDL-C determination is used to identify high-risk patients.




METHOD PRINCIPLE

The assay is a homogeneous method for directly measuring HDL-cholesterol concentration in serum or plasma, without any off-line pretreatment or centrifugation steps.

Accelerator selective detergent methodology.

During the first phase, LDL, VLDL particles and Chylomicrons generate free non-HDL cholesterol, which through an enzymatic reaction, produce hydrogen peroxide. The generated peroxide is consumed by a peroxidase reaction with DSBmT yielding a colourless product.

During the second phase, specific detergent solubilises HDL-Cholesterol. In conjunction with cholesterol oxidase (CO) and cholesterol esterase (CE) action, peroxidase and 4-AAP develop a coloured reaction which is proportional to HDL-Cholesterol concentration.

REAGENTS

Package

	Cat. No 4-179 (24-TRAY)	Cat. No 4-479 (36-TRAY)
1-Reagent	4 x 40 ml	6 x 23 ml
2-Reagent	4 x 15 ml	6 x 9 ml

The reagents are stable up to the kit expiry date printed on the package when stored at 2-8°C.

– URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁴

serum/ plasma	40 – 60 mg/dl
	1.04 – 1.55 mmol/l

As HDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control, it is recommended to use, with each batch of samples the CORMAY SERUM HN (Cat. No 5-172), CORMAY SERUM HP (Cat. No 5-173) and CORMAY LIPID CONTROL LEVEL 1 (Cat. No 5-179), CORMAY LIPID CONTROL LEVEL 2 (Cat. No 5-180).

For the calibration of automatic analysers the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended. The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium) with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or manual procedure is used.

- **LoB (Limit of Blank):**
0 mg/dl (0.0 mmol/l)
- **LoD (Limit of Detection):**
0.5 mg/dl (0.013 mmol/l)
- **LoQ (Limit of Quantitation):**
4.0 mg/dl (0.10 mmol/l)
- **Linearity:**
up to 155 mg/dl (4.01 mmol/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Bilirubin conjugated up to 60 mg/dl, bilirubin total up to 60 mg/dl, haemoglobin up to 1 g/dl, ascorbic acid up to 100 mg/dl, Intralipid up to 1800 mg/dl, triglycerides up to 2000 mg/dl and gamma-globulins up to 5000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	73.2	0.43	0.58
level 2	30.4	0.22	0.72
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	80.9	3.55	4.4
level 2	32.7	1.46	4.5

Concentrations in the test

1-Reagent

Buffer	
Cholesterol oxidase (<i>E.coli</i>)	< 1000 U/l
Peroxidase (horseradish)	< 1300 ppg U/l
N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT)	< 1 mM
Accelerator	< 1 mM
Preservative	< 0.06 %
Ascorbic acid oxidase (<i>Curcubita</i> spp.)	< 3000 U/l

2-Reagent

Buffer	
Cholesterol esterase (<i>Pseudomonas</i> spp.)	< 1500 U/l
4-aminoantipyrine (4-AAP)	< 1 mM
Detergent	< 2 %
Preservative	< 0.06 %

Warnings and notes

- Protect from direct sunlight and avoid contamination!

SPECIMEN

Serum, heparinized or EDTA plasma.

Anticoagulants containing citrate should not be used.

Blood should be collected only if the patient has been fasting for 12 -14 hours.

Serum: Collect whole blood by venepuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).

Plasma: Specimens may be collected in EDTA or lithium or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).

Serum and plasma should not remain at 15-30°C longer than 14 hours. If assays are not completed within 14 hours, serum or plasma should be stored at 2 - 8°C for up to 1 week. If specimens need to be stored for more than 1 week, they may be preserved at less than -20°C for up to 3 months. Samples may be frozen once.

Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT – LIPASE II GEN, HDL DIRECT – TG, HDL DIRECT – TG MONO, HDL-DIRECT – UA. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- **Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT – TG mono, HDL DIRECT – LIPASE, HDL DIRECT – UA, HDL DIRECT

Method comparison

A comparison between HDL cholesterol values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 58 serum samples gave following results:

$$y = 0.9994 x + 1.4363 \text{ mg/dl;}$$

$$R = 0.986 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl:1 4-13.
2. Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85:1234-41.
3. Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
5. Camps, J. Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol, Clinical Chemistry, 1999; 45:685-688.

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- Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl:1 4-13.
- Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85:1234-41.
- Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
- Camps, J. Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol. Clinical Chemistry, 1999; 45: 685-688.

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Item name	6	HDL-D
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	1.8	168
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range		
OFF		

▪ Prestige 24i Premium, Biolis 24i Premium

Item No.	6	Item Name	HDL-D	Optical
Data information				
Units	mg/dl			
Decimals	1			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	600nm			
Sub Wave Length				
Method	Direct			
Correlation				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	6	Item Name	HDL-D	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		4	5	µl
Reagent 1		180	10	µl
Reagent 2		60	10	µl
Data Process				
Read	Start	End	Abs.Limit	
Main	51	52	Low	-0.2
Sub	29	31	High	2
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	6	Item Name	HDL-D	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	40	60	40	60
Urine				
Plasma	40	60	40	60
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	6	Item Name	HDL-D	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			1.1	200
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
ON				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:DiI2				

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