PRESTIGE 24i MYOGLOBIN

Nr kat. **4-262, 4-481** (PL)

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Cat. No **4-262, 4-481** (EN)

INTENDED USE

Diagnostic kit for determination of myoglobin concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Myoglobin (Mb) is a hemo-protein present in cardiac and skeletal muscle cells and is released into blood circulation when these cells are damaged. The determination of serum Mb level is useful in the diagnosis of myocardial infarction, muscular dystrophy, myositis and myopathy, and also for the assessment of treatment and disease prognosis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between Mb in a sample and anti-Mb antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of Mb in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-262 (24-TRAY)	Cat. No 4-481 (36-TRAY)
1-Reagent	1 x 20.5 ml	2 x 13.5 ml
2-Reagent	1 x 8.5 ml	2 x 5.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 8 weeks (Hitachi 911).

Concentrations in the test

suspension of latex particles sensitized with anti-Mb (rabbit) antibodies (pH 7.3) 0.12 w/v%
 glycine buffer solution (pH 9.0)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 1-Reagent put on basic position in reagent tray.
 2-Reagent put on start position in reagent tray.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ³

serum, plasma	< 70 ng/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL II (Cat. No 4-290).

For the calibration of automatic analyzers systems the CORMAY MYOGLOBIN CALBRATORS kit (Cat. No 4-279) is recommended.

The calibration curve should be prepared every 8 weeks (Hitachi 911), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzer Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 2 ng/ml
- LoD (Limit of Detection):** 4 ng/ml
- LoQ (Limit of Quantitation):** 10 ng/ml
- Linearity:** up to 925 ng/ml

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.96 g/dl, ascorbate up to 500 mg/l, bilirubin up to 62 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	90	1.03	1.15
level 2	296	2.23	0.76
Reproducibility (day to day) n = 80	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	95.7	2.89	3.0
level 2	306.4	5.94	1.9

Method comparison

A comparison between myoglobin values determined at **Biolis 30i** (y) and at **BS-800** (x) using 69 serum samples, gave following results:

$$y = 0.992x - 3.5653 \text{ ng/ml}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between myoglobin values determined at **Biolis 30i** (y) and at **BS-800** (x) using 31 plasma samples, gave following results:

$$y = 0.9779x - 5.3264 \text{ ng/ml}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	43	MYO
Auto Rerun SW		
	OFF	
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Prozone Range		
	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	43	Item Name	MYO	Optical
Data information				
Units	ng/ml			
Decimals	0			
Analysis				
Type	RATE method			
Main Wave Length	570nm			
Sub Wave Length	800nm			
Method	Immunoturbidimetry			
Calibration				
Type	Spline1			
Std sample conc.				
Blank	0	#1 *	#2 *	
#3	*	#4 *	#5	
#6				
Correlation				
Y=	Slope	1	X+	Intercept
				0

Item No.	43	Item Name	MYO	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Main	Start	End	
	Sub	36	52	
Vol.				
	Kind	Vol.	Add	Units
	Sample	6	5	µl
	Reagent 1	160	10	µl
	Reagent 2	50	10	µl
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Correction value				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	99			
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	43	Item Name	MYO	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	70	0	70
Urine				
Plasma	0	70	0	70
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	43	Item Name	MYO	Optical
Auto Rerun SW				
	OFF			
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
	OFF			
Dilution				
	99:Di1			

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