PRESTIGE 24i LIPOPROTEIN (a)

Nr kat. 4-263, 4-482 (PL)

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Cat. No **4-263, 4-482** (EN)

INTENDED USE

Diagnostic kit for determination of lipoprotein (a) concentration intended to use automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lipoprotein(a) is a complex, cholesterol-carrying particle in the blood related to LDL which, when present at high levels, may be associated with the development of atherosclerosis and coronary heart disease, independent of LDL cholesterol and apoB. The structural component of Lp(a) distinguishing it from LDL is apolipoprotein(a), a large protein attached by disulfide bonding to the apoB-100 component of LDL. The similarity of the apo(a) sequence to those of plasminogen and hepatocyte growth factor suggests that the role of lipoprotein(a) in promoting the development of atherosclerosis may come from its capability to:

- 1) interfere in the breakdown of blood clots,
- 2) stimulate atherosclerotic cell proliferation.

Lp(a) levels is largely due to hereditary factors and considered to be useful in assessment of atherosclerosis risks.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between Lp(a) in a sample and anti-Lp(a) antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (700 nm), with the magnitude of the change being proportional to the quantity of Lp(a) in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-263 (24-TRAY)	Cat. No 4-482 (36-TRAY)
1-Reagent	1 x 38 ml	2 x 23 ml
2-Reagent	1 x 20 ml	2 x 12.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package.

Concentrations in the test

suspension of latex particles sensitized with anti-Lp(a) antibodies (rabbit) (pH 7.3) 0.4 w/v%
 glycine buffer solution (pH 9.0)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C.
- Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Application guide for Biolis 30i analyser is available on request.

Actions required:

- **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPOPROTEIN (a) – TG MONO. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

REFERENCE VALUES ⁴

serum, plasma	< 30 mg/dl
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY Lp(a) CONTROL N (Cat. No 4-492) and CORMAY Lp(a) CONTROL P (Cat. No 4-493)

For the calibration of automatic analysers systems the CORMAY Lp(a) CALIBRATORS kit (Cat. No 4-281) is recommended.

The calibration curve should be prepared every time the test is performed, with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers TBA-30R and Biolis 24i Premium. Results may vary if a different instrument is used.

- **Sensitivity:** 9.3 mg/dl.
- **Linearity:** up to 80 mg/dl.

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 6000 mg/dl, intralipid up to 500 mg/dl, RF up to 500 IU/ml do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	20.9	0.1	0.33
level 2	43.5	0.3	0.74

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 66 serum samples gave following results:

$$y = 1.108x - 1.44 \text{ mg/dl;}$$

$$R = 0.989 \quad (R - \text{correlation coefficient})$$

A comparison between CORMAY reagent (y) and commercially available assay (x) using 57 plasma samples gave following results:

$$y = 1.079x - 0.16 \text{ mg/dl;}$$

$$R = 0.990 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Prestige 24i Premium, Biolis 24i Premium

Item No.	44	Item Name	LpA	Optical
Data information				
Units	mg/dl			
Decimals	1			
Calibration				
Type	Spline1			
Std sample conc.				
Blank	0	#1	*	#2 *
#3	*	#4	*	#5 *
#6				
Analysis				
Type	END method			
Main Wave Length	700nm			
Sub Wave Length				
Method	Immuno			
Correlation				
	Slope	Intercept		
Y=	1	X+	0	
Item No.	44	Item Name	LpA	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample	5	5		µl
Reagent 1	160	10		µl
Reagent 2	80	10		µl
Data Process				
Read	Start	End	Abs.Limit	
Main	41	42	Low	-3
Sub	30	31	High	3
Correction value				
Blank correction		0.6735		
End Point Limit		2		
Linear Check (%)		0		
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Item No.	44	Item Name	LpA	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	30	0	30
Urine				
Plasma	0	30	0	30
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Item No.	44	Item Name	LpA	Optical
Auto Rerun SW				
OFF				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Absorbance Range				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
Prozone Range				
OFF				
Dilution				
99:Dil1				