PRESTIGE 24i RF

Nr kat. **4-269, 4-488** (PL)

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Cat. No. **4-269, 4-488** (EN)

INTENDED USE

Diagnostic kit for determination of rheumatoid factor levels used in automatic analysers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Rheumatoid factor (RF) is an autoantibody against human IgG commonly seen in sera at a high concentration in some conditions, particularly in patients with rheumatoid arthritis. The measurement of RF value is useful in evaluating the diagnosis, effects of therapy and prognosis of RA, systemic lupus erythematosus, chronic hepatopathy, etc.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between RF in a sample and denatured human IgG which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of RF in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-269 (24-TRAY)	Cat. No 4-488 (36-TRAY)
1-Reagent	2 x 40 ml	3 x 23 ml
2-Reagent	2 x 15 ml	3 x 9 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks.

Concentrations in the test

suspension of latex particles sensitized with denatured human IgG (pH 7.3) 0.17 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

In analysers Prestige 24i and Biolis 24i Premium for reagent blank 0.9% NaCl is recommended.

In analysers Biolis 30i for reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in analyser, there is a probability of cross-contamination affecting the tests results: RF – CALCIUM ARSENAZO, RF - CALCIUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in analyser, there is a probability of cross-contamination affecting the tests results: RF – CALCIUM ARSENAZO, RF – CALCIUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁴

serum, plasma	< 18 IU/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL I (Cat. No 4-288)

For the calibration of automatic analysers systems the CORMAY RF CALBRATORS kit (Cat. No 4-277) is recommended.

The calibration curve should be prepared every 4 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 0.0 IU/ml

- LoD (Limit of Detection):** 1.3 IU/ml
- LoQ (Limit of Quantitation):** 11 IU/ml
- Linearity:** up to 113 IU/ml

For higher concentration of RF dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.0 g/dl, bilirubin up to 66 mg/dl, intralipid up to 20 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[IU/ml]	[IU/ml]	[%]
level 1	18.4	0.35	1.88
level 2	37.2	0.31	0.84
Reproducibility (day to day) n = 80	Mean	SD	CV
	[IU/ml]	[IU/ml]	[%]
level 1	17.5	0.87	5.0
level 2	38.9	0.65	1.7

Method comparison

A comparison between RF values determined at **Biolis 30i** (y) and at **BECKMANN COULTER AU680** (x) using 66 serum samples gave following results:
 $y = 0.9762x + 0.7766$ IU/ml;
 $R = 0.997$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Galvin J.P. et al.: Particle enhanced photometric immunoassay systems., Clin Lab. Assay (Pap.Annu.Clin.Lab.Assays Conf.), 4th,73 (1983).
- Singer J.M. et al.: The latex fixation test. I. Application to the serologic diagnosis of rheumatoid arthritis, Amer. J. Med., 888 (1956)
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 960, (2006).
- Internal reference range studies.

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Item name	50	RF
Auto Rerun SW		
ON		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower OFF	Higher OFF
Auto Rerun Range (Result)		
	ON Lower	ON Higher
Serum	11.0	120
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
Prozone Range	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	50	Item Name	RF	Optical
Data information				
Units	IU/ml			
Decimals	1			
Analysis				
Type	END method			
Main Wave Length	660nm			
Sub Wave Length				
Method	Immuno			
Calibration				
Type	Spline1			
Std sample conc.				
Blank	0	#1 *	#2 *	*
#3	*	#4 *	#5 *	*
#6				
Correlation				
Slope	Intercept			
Y= 1	X+ 0			

Item No.	50	Item Name	RF	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	8	5	µl	
Reagent 1	240	10	µl	
Reagent 2	80	10	µl	
Data Process				
Read	Start	End		
Main	50	52		
Sub	36	37		
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Correction value				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
First	Start	End	Limit (%)	
Second				Low

Item No.	50	Item Name	RF	Optical
Normal Range				
	Male	Female		
	Low High	Low High		
Serum	0 18	0 18		
Urine				
Plasma	0 18	0 18		
CSF				
Dialysis				
Other				
Panic Range				
	Male	Female		
	Low High	Low High		
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	50	Item Name	RF	Optical	
Auto Rerun SW					
ON					
Auto Rerun Condition (Absorbance)					
Lower	OFF				
Higher	OFF				
Auto Rerun Range (Conc.)					
	First Dil	Low Value	Dil	High Value	Dil
Serum		11.3		120	
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Prozone)					
OFF					
Dilution					
99:Dil1					

PRESTIGE 24i RF

• **Biolis 30i**

Item no	50	Item name	RF	Specimen	SERUM/ PLASMA	OPTICAL												
Data information				Aspiration volume														
UNITS		IU/mL		TYPE														
DECIMALS		1		Double														
Analysis				<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>8</td> <td>240</td> <td>80</td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>				SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	8	240	80	BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	8	240	80															
BOTTLE (ml)																		
METHOD		END method		FIRST DIL.														
Main Wave Length		660 nm																
Sub Wave Length																		
CORRELATION (Y= AX + B)				Data processing read														
A =		1		START														
B =		0		END														
Blank value				<table border="1"> <tr> <td>MAIN</td> <td>51</td> <td>53</td> </tr> <tr> <td>SUB</td> <td>35</td> <td>36</td> </tr> </table>			MAIN	51	53	SUB	35	36						
MAIN	51	53																
SUB	35	36																
<table border="1"> <tr> <td>• WATER</td> <td>° REAGENT</td> </tr> </table>				• WATER	° REAGENT	ABS LIMIT												
• WATER	° REAGENT																	
				-3 TO 3														
Calibration				Collection value														
TYPE		Spline		END POINT														
STABILITY				2.5														
				LINEARITY CHECK (%)														
				0														
				Prozone check														
				<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> </table>				START	END	LIMIT (%)	FIRST				SECOND			
	START	END	LIMIT (%)															
FIRST																		
SECOND																		
				MINIMUM ABS.														
				MEAN														
				VARIATE														
				°HIGH														
				•LOW														

Item No	50	Item Name	RF	Specimen	SERUM/PLASMA	OPTICAL
Reference intervals						
MALE				FEMALE		
LOW	HIGH	LOW	HIGH			
0	18	0	18			
Panic range						
MALE				FEMALE		
LOW	HIGH	LOW	HIGH			
Decision limit						
		MALE		FEMALE		
Reaction check						
		°ON		•OFF		
CHECK						
LOW						
HIGH						
VL CHECK				VH CHECK		
°ON		•OFF		°ON		•OFF
Auto rerun						
		•ON		°OFF		
Auto rerun range (conc.)						
Re	Value	Dil.	Re	Value	Dil.	
	11			113		
Auto rerun condition (abs.)						
					DIL.	
LOWER	°ON		•OFF			
HIGH	°ON		•OFF			
Auto rerun condition (prozone)						
		°ON		•OFF		
SAMPLE VOL.						
Dilution						
		•DIL 1		° DIL 2		