PRESTIGE 24i LDL DIRECT

Nr kat. **4-180, 4-490** (PL)

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Cat. No. **4-180, 4-490** (EN)

INTENDED USE

Diagnostic kit for determination of LDL-cholesterol concentration (direct method) used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL).

LDL are synthesised in the liver by the action of various lipolytic enzymes on triglyceride rich VLDL. LDL-cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL-cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

METHOD PRINCIPLE

The assay is a homogenous method for directly measuring LDL-cholesterol concentrations in serum or plasma, without the need for any off-line pretreatment or centrifugation steps. Liquid selective detergent method.

The method is in a two reagents format and depends on the properties of a unique detergent. This detergent (1-REAGENT) solubilizes only the non LDL particles (HDL, VLDL, CM). The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (2-REAGENT) solubilizes the remaining LDL particles and chromogenic coupler allows for color formulation. The enzyme reaction with LDL-cholesterol in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

REAGENTS

Package

| | Cat. No 4-180 (24-TRAY) | Cat. No 4-490 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-REAGENT | 4 x 37.5 ml | 6 x 21.5 ml |
| 2-REAGENT | 4 x 14.5 ml | 6 x 8.5 ml |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 12 weeks, Biolis 24i Premium - 12 weeks.

Concentrations in the test 1-REAGENT

| | |
|---|----------------|
| Buffer | |
| Detergent 1 | < 1.0 % |
| Cholesterol esterase (<i>Pseudomonas</i> spp.) | < 1500 U/l |
| Cholesterol oxidase (<i>Cellulomonas</i> spp.) | < 1500 U/l |
| Peroxidase (horseradish) | < 1300 ppg U/l |
| 4-aminoantipyrine | < 0.1 % |
| Ascorbic Acid Oxidase (<i>Curcubita</i> spp.) | < 3000 U/l |
| Preservative | |
| 2-REAGENT | |
| Buffer | |
| Detergent 2 | < 1.0 % |
| N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT) | < 1.0 mM |
| Preservative | |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

Warning

H317 May cause an allergic skin reaction.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water

SPECIMEN

Serum, sodium heparinized or EDTA plasma. Anticoagulants containing citrate should not be used. Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours). Plasma: Specimens may be collected in EDTA or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours). If not analysed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens need to be stored for more than 5 days, they may be preserved at -80°C. Samples may be frozen once. Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.
1-REAGENT put on basic position in reagent tray.
2-REAGENT put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LDL DIRECT – LIPASE, LDL DIRECT – LIPASE II GEN, LDL DIRECT – TG, LDL DIRECT – TG MONO, LDL DIRECT – UA. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LDL DIRECT – TG, LDL DIRECT – LIPASE, LDL DIRECT – UA, LDL DIRECT – URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES

| NCEP* Classification: | Adults | |
|--------------------------|-----------|-------------|
| | mg/dl | mmol/l |
| Optional | < 100 | < 2.59 |
| Near optional | < 130 | < 3.37 |
| Bordeline high | 130 - 159 | 3.37 - 4.12 |
| High | 160 - 189 | 4.14 - 4.90 |
| Very high | ≥ 190 | ≥ 4.92 |

* National Cholesterol Education Program

As LDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control, it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172), CORMAY SERUM HP (Cat. No 5-173) and CORMAY LIPID CONTROL LEVEL 1 (Cat. No 5-179), CORMAY LIPID CONTROL LEVEL 2 (Cat. No 5-180).

For the calibration of automatic analysers the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended. The calibration curve should be prepared every 12 weeks (Prestige 24i) or every 12 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.6 mg/dl (0,016 mmol/l)
- LoD (Limit of Detection):**
0.7 mg/dl (0,018 mmol/l)
- LoQ (Limit of Quantitation):**
4,5 mg/dl (0,12 mmol/l)
- Linearity:**
up to 700 mg/dl (18,13 mmol/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Triglycerides up to 1293 mg/dl, bilirubin conjugated up to 20 mg/dl, bilirubin total up to 20 mg/dl, haemoglobin up to 0.5 g/dl, ascorbic acid up to 500 mg/l and gamma-globulins up to 5000 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--|-----------------|---------------|-----------|
| level 1 | 140.5 | 2.24 | 1.60 |
| level 2 | 58.2 | 1.02 | 1.75 |
| Reproducibility (day to day) n = 80 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
| level 1 | 145.1 | 7.3 | 5.0 |
| level 2 | 66.5 | 3.33 | 5.0 |

Method comparison

A comparison between LDL cholesterol values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 54 serum samples gave following results:
y = 1.1066 x - 4.6206 mg/dl;
R = 0.992 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | | |
|----------------------|-------|-----|-----------|
| Prozone Check | | | |
| | Start | End | Limit (%) |
| First | | | |
| Second | | | Low |
| Third | | | Low |

| | | | |
|--|-------|--------|-----|
| Item name | 25 | LDL-D | |
| Auto Rerun SW | | | |
| ON | | | |
| Auto Rerun Condition (Absorbance) | | | |
| Absorbance Range | | Lower | OFF |
| | | Higher | OFF |
| Auto Rerun Range (Result) | | | |
| | ON | ON | |
| | Lower | Higher | |
| Serum | 4.2 | 199 | |
| Urine | | | |
| Plasma | | | |
| CSF | | | |
| Dialysis | | | |
| Other | | | |
| Prozone Range | | | |
| OFF | | | |

• Prestige 24i Premium, Biolis 24i Premium

| | | | | |
|-------------------------|------------|-----------|-------|---------|
| Item No. | 25 | Item Name | LDL-D | Optical |
| Data information | | | | |
| Units | mg/dl | | | |
| Decimals | 1 | | | |
| Calibration | | | | |
| Type | Linear1 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 |
| #3 | | #4 | | #5 |
| #6 | | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 600nm | | | |
| Sub Wave Length | 700nm | | | |
| Method | Direct | | | |
| Correlation | | | | |
| Slope | | Intercept | | |
| Y= | 1 | X+ | 0 | |

| | | | | |
|----------------------|--------|-----------|-------|---------|
| Item No. | 25 | Item Name | LDL-D | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Vol. | | | | |
| | Kind | Vol. | Add | Units |
| Sample | 3 | 5 | | µl |
| Reagent 1 | 210 | 10 | | µl |
| Reagent 2 | 70 | 10 | | µl |
| Data Process | | | | |
| Read | Start | End | | |
| | Main | 51 | 52 | |
| | Sub | 29 | 31 | |
| Abs.Limit Low ~ High | | | | |
| | -0.2 | 2 | | |

| | | | | |
|-------------------------|-------|-----|-----------|--|
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | 1 | | | |
| Span | 3 | | | |
| Third mixing | | | | |
| OFF | | | | |
| Correction value | | | | |
| Blank correction | | | | |
| End Point Limit | | | | |
| 2 | | | | |
| Linear Check (%) | | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | Low | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 25 | Item Name | LDL-D | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | 0 | 100 | 0 | 100 |
| Urine | | | | |
| Plasma | 0 | 100 | 0 | 100 |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | |
|--|-----------|-----------|-------|---------|
| Item No. | 25 | Item Name | LDL-D | Optical |
| Auto Rerun SW | | | | |
| ON | | | | |
| Auto Rerun Condition (Absorbance) | | | | |
| | | Lower | OFF | |
| | | Higher | OFF | |
| Auto Rerun Range (Conc.) | | | | |
| | First Dil | Low | High | |
| | | Re | Value | Dil |
| Serum | | | 4.2 | 700 |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Auto Rerun Condition (Prozone) | | | | |
| OFF | | | | |
| Dilution | | | | |
| 100:Di2 | | | | |