## PRESTIGE 24i UIBC

Nr kat. **4-186, 4-494** (PL)

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### LITERATURA

- Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- Wick M, Pingerra W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- Guder Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- Perrotta, G., Iron and Iron-Binding Capacity, In: Pesce, A.J., Kaplan, L.A. eds., Methods in Clinical Chemistry, C.V. Mosby, St. Louis, 1258 - 1261,1987
- Tietz NW (ed). Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1701-1703; 1821t (1999).
- Burtis CA, Bruns DE. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 7th ed. St. Louis, p. 306. 2014
- Schreiber WE. „Iron and Porphyrin Metabolism” in Pesce, A.J., Kaplan, L.A. editors. Clinical Chemistry: Theory, Analysis, Correlation. Mosby Inc., an affiliate of Elsevier Inc., St. Louis, 755-770, 2010

Data wydania: 04.2021.

## PRESTIGE 24i UIBC

Cat. No **4-186, 4-494** (EN)

### INTENDED USE

Diagnostic kit for determination of unsaturated iron binding capacity intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.


### INTRODUCTION

The total iron content of the body is about 3 to 3.5 g. Of this amount about 2.5 g contained in erythrocytes or their precursors in the bone marrow. Plasma contains only about 2.5 mg of iron. Iron is transported as Fe (III) bound to the plasma protein apotransferrin. The apotransferrin-Fe (III) complex is called transferrin. Normally only about one third of the iron binding sites of transferrin are occupied by Fe (III). The additional amount of iron that can be bound is the unsaturated (or latent) iron-binding capacity (UIBC). The sum of the serum iron and UIBC represents total iron binding capacity (TIBC). TIBC is a measurement for the maximum iron concentration that transferrin can bind.

Serum UIBC levels vary in disorders of iron metabolism where iron capacities are often increased in iron deficiency and decreased in chronic inflammatory disorders, malignancies or in course of hemochromatosis.

### METHOD PRINCIPLE

Direct, colorimetric method with ferrozine:

In an alkaline environment known ferrous ion concentration incubated with serum, binds specifically with transferrin at unsaturated iron binding sites. Remaining unbound ferrous ions are measured with the chromogen reaction.

The difference between the amount of excess iron and the total amount added to the serum is equivalent to the quantity bound to transferrin. This is the UIBC (unsaturated iron binding capacity) of the sample.

### REAGENTS

#### Package

	Cat. No 4-186 (24-TRAY)	Cat. No 4-494 (36-TRAY)
1-Reagent	1 x 40 ml	1 x 23 ml
2-Reagent	1 x 12 ml	1 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C Prestige 24i, Biolis 24i Premium – 8 weeks.

### Concentrations in the test

<b>1-Reagent</b>	
buffer (pH 8.4)	0.25 mol/l
ammonium iron (II) sulfate	20 µmol/l
thiourea	90 mmol/l
detergent	0.1 %
sodium azide	<0.1 %

### 2-Reagent

sodium ascorbate	150 mmol/l
sodium chloride	75 mmol/l
3-(2-pyridyl)-5,6-bis(2-[5-furyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) preservatives	≥ 10 mmol/l 0.3 %

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagents.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.
- Products contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.
- The negative UIBC value will occur when the patient's serum iron level exceeds the binding capacity of the transferrin.
- For the diagnostic purpose UIBC determination should be performed at the same time with iron determination. Obtained result have to be interpreted in relation to the result of iron concentration and the percentage saturation of transferrin with iron ions<sup>1</sup>.
- 1-Reagent contains thiourea. May produce an allergic reaction (EUH208).
- 2-Reagent contains 1-[1,3-Bis (hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis (hydroxymethyl) urea. May produce an allergic reaction (EUH208).

### SPECIMEN

Serum, heparin plasma.

Separate serum/plasma at the latest 2 h after blood collection to avoid hemolysis. Samples should be taken in the morning from patients, since iron values decrease during the course of the day.

Discard contaminated specimens.

Anticoagulants such as EDTA, oxalate and citrate must not be used, as they bind iron ions and prevent reaction with chromogen<sup>4</sup>.

Serum can be stored up to 3 days at 20-25°C, 7 days at 4-8°C or up to one month at -20°C. Plasma can be stored up to 7 days at 4-8°C or up to month at -20°C<sup>4</sup>.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyzer there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of unsaturated iron binding capacity should be performed **in the separate order** (follow the recommendations contained in the instruction: 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER).
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results. To avoid this effect follow the recommendations contained in the instruction: 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES 5,6

The reference values were calculated from the serum iron (SI) and TIBC ranges reported in literature, in accordance with mathematic formula:

$$\text{UIBC} = \text{TIBC} - \text{SI}$$

Reference values for UIBC are listed in table below:

serum / plasma	µg/dl	µmol/l
Females	80 – 375	14 – 67
Males	75 – 360	13 – 64

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) is recommended.

The calibration curve should be prepared every 7 days (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):**  
4.0 µg/dl (0.72 µmol/l)
- LoD (Limit of Detection):**  
8.0 µg/dl (1.43 µmol/l)
- LoQ (Limit of Quantitation):**  
40 µg/dl (7.16 µmol/l)
- Linearity:**  
up to 533 µg/dl (95.41 µmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin interfere even in small amounts, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl, copper up to 3.5 mg/dl and zinc up to 15 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	92.00	1.79	1.94
level 2	124.00	2.77	2.23
Reproducibility (day to day) n = 80	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	99.0	8.30	8.4
level 2	139.9	7.63	5.5

#### Method comparison

A comparison between UIBC values determined at **Biolis 30i** (y) and at **BS-400** (x) using 62 serum samples, gave following results:

$$y = 0.9357x + 8.8541 \text{ µg/dl ;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- Wick M, Pingerra W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- Guder Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- Perrotta, G., Iron and Iron-Binding Capacity, In: Pesce, A.J., Kaplan, L.A. eds. Methods in Clinical Chemistry, C.V. Mosby, St. Louis, 1258 - 1261,1987
- Tietz NW (ed). Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1701-1703; 1821t (1999).
- Burtis CA, Bruns DE. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 7th ed. St. Louis, p. 306. 2014
- Schreiber WE. „Iron and Porphyrin Metabolism” in Pesce, A.J., Kaplan, L.A. editors. Clinical Chemistry: Theory, Analysis, Correlation. Mosby Inc., an affiliate of Elsevier Inc., St. Louis, 755-770, 2010

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- Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- Wick M, Pingerra W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- Guder Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- Perrotta, G., Iron and Iron-Binding Capacity, In: Pesce, A.J., Kaplan, L.A. eds., Methods in Clinical Chemistry, C.V. Mosby, St. Louis, 1258 - 1261,1987
- Tietz NW (ed). Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1701-1703; 1821t (1999).
- Burtis CA, Brun DE. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 7th ed. St. Louis, p. 306. 2014
- Schreiber WE. „Iron and Porphyrin Metabolism” in Pesce, A.J., Kaplan, L.A. editors. Clinical Chemistry: Theory, Analysis, Correlation. Mosby Inc., an affiliate of Elsevier Inc., St. Louis, 755-770, 2010

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Item name	32	UIBC
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	20	500
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range		
	Lower	OFF
	Higher	OFF
<b>Prozone Range</b>		
OFF		

- Prestige 24i Premium, Biolis 24i Premium

Item No.	32	Item Name	UIBC	Optical
<b>Data information</b>				
Units	µg/dl			
Decimals	0			
<b>Calibration</b>				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Analysis</b>				
Type	END method			
Main Wave Length	570 nm			
Sub Wave Length	800 nm			
Method	Ferrozine			
<b>Correlation</b>				
Slope	Intercept			
Y=	1	X+	0	

Item No.	32	Item Name	UIBC	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample	20		5	µl
Reagent 1	200		10	µl
Reagent 2	50		10	µl
<b>Data Process</b>				
Read	Start	End		
	Main	50	52	
	Sub	30	31	
<b>Abs.Limit</b>				
Low	High			
	-3	~		3
<b>Blank value</b>				
Water Blank				
<b>Correction value</b>				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	32	Item Name	UIBC	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	75	360	80	375
Urine				
Plasma	75	360	80	375
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	32	Item Name	UIBC	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low Value	High Value	
	Re	Dil	Re	Dil
Serum		17	550	
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
Lower				
OFF				
Higher				
OFF				
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:Dil2				

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