

CORMAY VALPROIC ACID CONTROLS



CONTROL SERUM FOR DETERMINATION OF VALPROIC ACID CONCENTRATION

Kit name	Kit size	Cat. No
CORMAY VALPROIC ACID CONTROLS	2 x 5 ml	5-152

CORMAY VALPROIC ACID CONTROLS intended for quality control of valproic acid concentration measurement's.

PACKAGE

	Size
CONTROL 1	1 x 5 ml
CONTROL 2	1 x 5 ml

PROCEDURE

1. Open the vial carefully and add exactly 5 ml of deionised water at 20-25°C.
2. Close the bottle and allow to stand for 30 minutes before use.
3. Swirl gently several times during reconstitution period to ensure that the contents are completely dissolved.
4. Prior to use, mix the contents by inverting the vial several times. Do not shake! Avoid the foam formation! Refrigerate any unused material.

STABILITY AND STORAGE

- Lyophilised serum remains stable when stored at 2-8°C until expiry date given on the product label.
- The constituents in reconstituted control serum are stable for 4 weeks when stored at 2-8°C and kept capped in original container, free from contamination.
- Mix controls thoroughly before use.
- Protect from light and avoid contamination!

WARNINGS AND NOTES

- Products for in vitro diagnostic use only.
- Controls should be used by suitably qualified laboratory personnel only in accordance with intended purpose.
- The controls are based on human serum that has been tested for HBsAg and antibodies to HIV 1, HIV 2 and HCV and found to be non-reactive. However this material and all patient samples should be handled as though capable of transmitting infectious disease.

ASSIGNED VALUES

CONTROL	LOT	MEAN VALUE	RANGE
CONTROL 1	611-10	545 µmol/l 78.6 µg/ml	436 – 654 µmol/l 62.9 – 94.3 µg/ml
CONTROL 2	611-09	918 µmol/l 132 µg/ml	734 – 1102 µmol/l 106 – 158 µg/ml

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 11. 2016.

MANUFACTURER

PZ CORMAY S.A.
22 Wiosenna Street,
05-092 Łomianki, POLAND
tel.: +48 (0) 22 751 79 10
fax: +48 (0) 22 751 79 14
<http://www.cormay.pl>