CORMAY CYSTATIN C CALIBRATORS

Nr kat. 5-185 (PL)

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Cat. No 5-185 (EN)

INTENDED USE

CORMAY CYSTATIN C CALIBRATORS kit is dedicated for calibration of cystatin C determinations at human serum and plasma on automatic analysers.

The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Cystatin C is a low molecular weight protein (13kD), one of the cysteine proteinases inhibitors. Cystatin C is produced in all nucleated cells and secreted into extracellular space at constant rate. Cystatin C molecule stability and dependence of its concentration solely on GFR (Glomerular Filtration Rate) decide about high diagnostic efficiency of cystatin C determination. Cystatin C level is not affected by muscle mass or diet and its increase is observed even at slight GFR reduction. Clinical applications of cystatin C are for monitoring GFR in children, elderly patients, patients with potentially nephrotoxic drug therapy, for assessment of renal transplantation status, for kidney function monitoring in acute and chronic kidney diseases including diabetic nephropathy.

PACKAGE

	Size
CALIBRATOR 1	1 x 1 ml
CALIBRATOR 2	1 x 1 ml
CALIBRATOR 3	1 x 1 ml
CALIBRATOR 4	1 x 1 ml
CALIBRATOR 5	1 x 1 ml
CALIBRATOR 6	1 x 1 ml

STABILITY AND STORAGE

- Calibrators are ready to use.
- Unopened calibrators are stable at 2-8°C until the expiration date printed on the package.
- Calibrators after opening are stable for 6 months, if stored at 2-8°C thoroughly closed.
- Do not allow calibrators to stand for long time periods at room temperature.
- Do not freeze calibrators!

WARNINGS AND NOTES

- Do not use calibrators after expiry date.
- Serum used in the manufacturing of calibrators are tested for HBsAg, anti-HCV, anti-HIV1 and anti-HIV2 and found to be negative. Nevertheless, the materials contain substances of human origin and must be handled with due care.
- Calibrators contain antibiotics as a preservative and must be handled with due care.

ASSIGNED VALUES

The exact concentration of cystatin C is printed on the label of the each vial.

TRACEABILITY

Cystatin C calibrator is calibrated by the use of pure human cystatin C.

WASTE MANAGEMENT

Please refer to local legal requirements.

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