

## Instruction for Malaria Pf/Pv Antigen Test Kit (Immunofluorescence)

### 1. PRODUCT NAME

Generic name: Malaria Pf/Pv Antigen Test Kit (Immunofluorescence)

Trade name: Malaria Pf/Pv

### 2. PACKAGE

Specification 1: 25T/kit REF:52026132

Specification 2: 50T/kit REF:52027123

Quality Control (optional):

Specification: Level 1 : 0.5mL x 1 REF: 52105119

Level 2 : 0.5mL x 1 REF: 52105120

### 3. INTENDED USE & INDICATION

For in vitro qualitative determination of Pf/Pv antigen in human serum, plasma, or whole blood. The test targets the antigen specific to Plasmodium falciparum (P.f.) and a pan-malarial antigen to P. vivax (P.v.). It is intended to aid in the rapid diagnosis of human malaria infections and to aid in the differential diagnosis of Plasmodium falciparum (P.f.) infections from other less virulent malarial infections. Negative results must be confirmed by thin / thick smear microscopy. The test is not intended for use in screening asymptomatic populations.

For professional use only.

### 4. TEST PRINCIPLE

This kit adopts the antibody sandwich method. After the samples containing Pf or Pv are taken and added into the sample well, the Pf in the sample binds to the Pf monoclonal antibody coupled to the fluorescent particles to form a fluorescent particle-antibody-antigen complex. This immune complex is then chromatographed along the nitrocellulose membrane to the determination area (T1), and combined with the pre-coated Pf monoclonal antibody. The Pv in the sample binds to the Pv monoclonal antibody coupled to the fluorescent particles to form a fluorescent particle-antibody-antigen complex. This immune complex is then chromatographed along the nitrocellulose membrane to the determination area (T2), combined with the pre-coated Pv monoclonal antibody, and continues to chromatograph to the quality control area (C). The fluorescent particle-labeled goat anti-rabbit IgG binds to the precoated rabbit IgG and presents the quality control line. The fluorescence intensity of the assay area is directly proportional to the level of Pf or Pv in the sample.

### 5. MAIN COMPONENTS& ADDITIONAL REQUIRED EQUIPMENT

The test kit consists of test card, magcard, sample diluent, quality control (optional) and the instruction.

(1) The test card consists of card shell and test strip. The test strip contains sample pad/marketing pad, nitrocellulose membrane, absorbent paper and PVC plate.

(2) Magcard: load calibration curve information for this batch of reagents.

(3) Sample diluent: The main ingredient is phosphate buffer(PBS).

(4) Quality control (optional): Self-prepared lyophilized powders, mainly consist of Pf/Pv recombinant antigen and PBS. All are free of human-derived substances and have batch specificity. Please find target values in the target value list.

(5) Equipment: Applicable to FA50 and FA120 Quantitative Immunoassay Analyzer manufactured by Genrui Biotech Inc.

Note: Components of kits from different batches should not be used interchangeably.

### 6. ACCESSORIES REQUIRED BUT NOT PROVIDED

(1) Pipettes and pipette tips: 100  $\mu$ L

(2) Timer

### 7. SPECIAL STORAGE &TRANSPORT CONDITIONS

(1) The test kit should be stored at 2-30 $^{\circ}$ C, and the shelf life of test cards and sample diluent is 18 months when sealed. After the test card and sample diluent are opened, the shelf life is 1 hour at 18-30 $^{\circ}$ C and 40%-65% humidity. When the humidity is > 65%, it should be used right after opened.

(2) The unopened QC is stable for 18 months (see the label for specific date) at -25 $^{\circ}$ C to 8 $^{\circ}$ C, the reconstituted QC is stable for 6 days at -20 $^{\circ}$ C or 1 day at 2-8 $^{\circ}$ C in the shade, and can be freeze-thawed once.

(3) Transport: The test kit is at 2-30 $^{\circ}$ C, the QC is at -25-8 $^{\circ}$ C.

### 8. SAMPLE REQUIREMENTS

(1) The optimal sample is fresh non-hemolyzed serum, plasma or whole blood. It is recommended to use sample from venous blood, as results of other body fluids and samples may not be accurate.

(2) Serum/plasma: After sample collection, serum should be separated as soon as possible to avoid hemolysis. Serum and plasma should complete the test within 6 hours at room temperature. The samples that cannot complete the test should be refrigerated at 2-8 $^{\circ}$ C for no more than 7 days; serum and plasma should be frozen below -18 $^{\circ}$ C for no more than 1 month.

(3) Whole blood: It should be used immediately after collection. If it cannot be tested within 4 hours, it should be refrigerated at 2-8 $^{\circ}$ C for no more than 3 days. Samples should not be frozen.

(4) The samples should be brought to room temperature before determination. The frozen samples should be completely thawed, rewarmed and mixed well before use. Do not freeze and thaw repeatedly.

(5) Human serum is preferred for determination, and EDTA-K<sub>2</sub> is recommended as

an anticoagulant for plasma and whole blood testing.

### 9. TEST METHOD

Carefully read the instruction before using the test kit and operate in strict accordance with the instruction to ensure reliable results. Bring all reagents to room temperature (18-30 $^{\circ}$ C) before use.

(1) Startup: Click "STD Mode" in the main menu to enter the measurement interface, click "Item" to select the test item and click "Type" to select the sample type.

(2) Click "Lot No." to enter the card reading interface, place magcard of the corresponding item to the magnetic card reader area, when the magcard is read successfully, check whether the magcard and the test card are of the same batch. (Note: reagents are precalibrated and specific calibration curve parameters for each batch of reagents have been stored in the magcard.)

(3) Quality control procedure: It is recommended to refer to the instrument manual and use the Genrui quality control to verify whether the target value of the test quality control is under control during the measurement procedure after calibration. The quality controls should be used as follows.

a) Bring the quality control to room temperature (18-30 $^{\circ}$ C) before use.

b) Carefully open the bottle cap to avoid spraying of the contents.

c) Add 0.5 mL of purified water.

d) Put on the bottle cap and leave it at room temperature for 15 minutes, gently shake the bottle to fully dissolve the dry powder.

e) After the dry powder is fully dissolved, repeat the operation for the sampling.

If the measured values of quality controls meet the following: the Level 1 is negative, the Level 2 is positive, the determination of clinical samples and data analysis can be continued; otherwise, the causes should be identified before test.

(4) Sampling:

Add 0.1mL of serum, plasma or whole blood into the container with sample diluent, mix thoroughly. Take 0.1mL of diluted sample, and drop it vertically to the sample well on the test card directly and start timing.

(5) Insert it into the analyzer's test card slot (the sample well end towards the inside). Click "Measure", the instrument will automatically detect and print out the results after 15 minutes (If using "Fast Mode", after 15 minutes of external incubation, quickly insert card and click "Measure", then instantly the instrument will detect and print out the results).

Note: For detailed instructions on how to operate the instrument, please refer to the manual of Quantitative Immunoassay Analyzer.

### 10. INTERPRETATION OF TEST RESULT

The results of instrument printing are presented in the form of fluorescent signal value (FSV) as follows.

Fluorescent signal value (FSV)	Interpretation	Note
Pf: $\leq$ 0.9 Pv: $\leq$ 0.9	Negative for Malaria Pf and Pv Ag	No need to additional test
Pf: $>$ 0.9 Pv: $>$ 0.9	Positive for Malaria Pf and Pv Ag	Need for further confirmation test.

### 11. INTERFERING SUBSTANCE

(1) Cross-reactivity

Test 28 pathogenic microorganisms (7 bacteria, 5 protists and 16 viruses) that may be present in sample. All were negative when tested at the concentrations listed below.

Type	Pathogen Tested	Concentration Tested
Bacteria	<i>Borrelia burgdorferi</i> (N40 strain)	2.3 x 10 <sup>6</sup> organisms/mL
	<i>Leptospira interrogans (icterohaemorrhagiae)</i>	1.0 x 10 <sup>7</sup> organisms/mL
	<i>Leptospira biflexa (andamana)</i>	1.0 x 10 <sup>7</sup> organisms/mL
	<i>Treponema pallidum</i>	1.0 x 10 <sup>5</sup> organisms/mL
	<i>Rickettsia conorii</i> (Malish 7)	1.0 x 10 <sup>7</sup> organisms/mL
	<i>Rickettsia typhi</i> (Wilmington)	1.0 x 10 <sup>7</sup> organisms/mL
	<i>Orientia tsutsugamushi-Rickettsia (Karp)</i>	1.0 x 10 <sup>7</sup> organisms/mL
Protists	<i>Babesia microti</i> (RMNS strain)	4.4 x 10 <sup>7</sup> parasites/mL
	<i>Trypanosoma cruzi</i> (Y strain)	1.3 x 10 <sup>6</sup> parasites/mL
	<i>Leishmania donovani</i>	1.0 x 10 <sup>6</sup> parasites/mL
	<i>Leishmania infantum</i>	1.0 x 10 <sup>6</sup> parasites/mL

Type	Pathogen Tested	Concentration Tested
	<i>Leishmania chagasi</i>	1.0 x 10 <sup>6</sup> parasites/mL
Viruses	Cytomegalovirus (CMV) (AD169)	1.2 x 10 <sup>5</sup> PFU/mL
	Epstein-Barr Virus (EBV)	1.1 x 10 <sup>4</sup> copies/mL
	Dengue virus – West Pac 74	1.2 x 10 <sup>5</sup> PFU/mL
	Dengue virus – S16803	3.9 x 10 <sup>4</sup> PFU/mL
	Dengue virus – CH53489	1.3 x 10 <sup>4</sup> PFU/mL
	Dengue virus – TVP360	1.4 x 10 <sup>5</sup> PFU/mL
	Yellow Fever virus	7.9 x 10 <sup>6</sup> PFU/mL
	West Nile virus	1.6 x 10 <sup>5</sup> PFU/mL
	Chikungunya virus	4.0 x 10 <sup>5</sup> PFU/mL
	Ross-River virus	1.0 x 10 <sup>6</sup> PFU/mL
	Influenza A – Bayem/7/95	2.5 x 10 <sup>7</sup> TCID <sub>50</sub> /mL
	Influenza B – Victoria/2/87	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL
	HIV-1 (Subtype B)	1.4 x 10 <sup>5</sup> copies/mL
	Hepatitis B	2.0 x 10 <sup>5</sup> IU/mL
	Hepatitis C	1.9 x 10 <sup>5</sup> IU/mL
Rubella Virus	≥ 2.0 x 10 <sup>2</sup> TCID <sub>50</sub> /mL	

(2) Interference

The following substances that may be artificially introduced into sample were evaluated at the concentrations listed and were found not to affect test performance. Note: The analytical effects of these drugs were studied by taking sample and spiking it with quantities at high therapeutic concentrations and then testing these samples. The effects of the clinical metabolites of these drugs on the test were not studied.

Substance Type	Substance	Concentration
Anti-malarial drugs (prevention)	Mefloquine (Lariam™)	1 mg/mL
	Doxycycline* (Vibramycin™)	1 mg/mL
	Chloroquine	1 mg/mL
	Hydroxychloroquine sulfate	1 mg/mL
	Paludrine™ (Proguanil)	1 mg/mL
	Primaquine	1 mg/mL
	Quinine	1 mg/mL
Antibiotic (treatment)	Sulfadoxine and Pymethamine (Fansidar™)	1 mg/mL
	Amoxicillin (Trimox™)	0.1 mg/mL
	Cephalexin	0.1 mg/mL
	Ciprofloxacin	0.1 mg/mL
Anti- Inflammatory Drugs (treatment)	Erythromycin	0.1 mg/mL
	Aspirin	1 mg/mL
	Acetaminophen	1 mg/mL
	Ibuprofen (NSAID)	1 mg/mL

12. PERFORMANCE CHARACTERISTICS

511 samples which include 173 confirmed as Pf positive and 338 confirmed as Pf negative by contrast reagent, were obtained for testing, and then compared the test results of Genrui Malaria Pf/Pv Ag Test Kit with contrast reagent results. The results are shown below.

Pf	Contrast reagent		Subtotal	
	Positive	Negative		
Malaria Pf/Pv Ag Test Kit	Positive	167	2	169
	Negative	6	336	342
Subtotal		173	338	511

Sensitivity: 96.53% (95%CI: 92.60%~98.72%)

Specificity: 99.41% (95%CI: 97.88%~99.93%)

Overall Percent Agreement: 98.43% (95%CI: 96.94%~99.32%)

525 samples which include 177 confirmed as Pv positive and 348 confirmed as Pv negative by contrast reagent, were obtained for testing, and then compared the test results of Genrui Malaria Pf/Pv Ag Test Kit with contrast reagent results. The results are shown below.

Pv	Contrast reagent		Subtotal	
	Positive	Negative		
Malaria Pf/Pv Ag Test Kit	Positive	172	1	173
	Negative	5	347	352
Subtotal		177	348	525

Sensitivity: 97.18% (95%CI: 93.53%~99.08%)

Specificity: 99.71% (95%CI: 98.41%~99.99%)

Overall Percent Agreement: 98.86% (95%CI: 97.53%~99.58%)

13. PRODUCT PERFORMANCE INDICATORS

(1) Limit of detection:

Titre of cultured virus(pfu/mL)	Fluorescent signal value (FSV)	Result
Pf:2.15x10 <sup>2</sup>	0.93	Pos
Pv:3.21x10 <sup>2</sup>	0.91	Pos

(2) Precision: intra-batch precision: CV ≤ 15%; inter-batch precision of the kit CV ≤ 15%

14. PRECAUTIONS

(1) Once opened, use the test cards as soon as possible, which may be exposed to moisture in the air. Do not reuse the test cards.

(2) Components in test kit of different batches cannot be used interchangeably.

(3) For substances containing sources of infection or suspected of containing sources of infection, there should have proper bio-safety assurance procedures. Pay attention to the following notes:

-- Wear gloves when handling sample or disinfecting the reagent.

-- Disinfect spilled sample or reagent with disinfectant.

-- Disinfect or handle potential contamination sources of all samples or reagents in accordance with local regulations.

15. EXPLANATION OF GRAPHIC SYMBOL

	Consult instructions for use		Temperature limit
	Batch code		Use-by date
	<i>In vitro</i> diagnostic medical device		CE Marking
	Date of manufacture		Biological risks
	Manufacturer		Volume
	Contains sufficient for < n>tests		Keep away from sunlight
	Do not re-use		Keep dry
	Authorized representative in the European community		Catalogue number

16. REFERENCE

(1) Breaman, J.G., M.S. Alilio, and A. Mills. Conquering the intolerable burden of malaria: what's new, what's needed: a summary. American J. of Tropical Medicine and Hygiene, 2004;71 (Suppl 2):1-15.

(2) Centers for Disease Control (CDC). Treatment of Malaria (Guidelines for Clinicians), June 28, 2004.

17. HELP INFORMATION

If you need help, please contact after sales department.

18. MANUFACTURER

Genrui Biotech Inc.

Address: 4-10F, Building 3, Geya Technology Park, Guangming District, 518106, Shenzhen, China.

Web: www.genrui-bio.com

19. INSTRUMENTS & APPLICATIONS

Genrui's Immunofluorescence products are designed to work in automated lab, which are compatible with the FA50/FA120 Quantitative Immunoassay Analyzer. There may or may not be an application developed for your particular instrument, please visit the instrument section of our website.



Lotus NL B.V.  
Koningin Julianaplein 10, 1e Verd  
2595AA, The Hague, Netherlands  
Email:peter@lotusnl.com