## CORMAY ASO LATEX

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### LITERATURA

- Haffjee. Quarterly Journal of Medicine 1992. New series 84; 305: 641-658.
- Ahmed Samir et al. Pediatric Annals 1992; 21: 835-842.
- Spain J et al. Bull Wild Hlth Org 1961; 24: 271-279.
- The association of Clinical Pathologists 1961. Broadsheet 34.
- Picard B et al. La Presse Medicale 1983; 23: 2-6.
- Klein GC. Applied Microbiology 1971; 21: 999-1001.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

Data wydania: 09. 2023.

## CORMAY ASO LATEX

Kit name	Kit size	Cat. No	(EN)
CORMAY ASO Latex 100	100 tests	6-251	

### INTENDED USE

Agglutination test for detection of anti-streptolysin O (ASO) antibodies, intended to use for the qualitative and semi-quantitative manual assays.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Streptolysin O is a toxic immunogenic exoenzyme produced by  $\beta$ -hemolytic Streptococci of groups A, C and G. Measuring the ASO antibodies is useful for the diagnostic of rheumatic fever, acute glomerulonephritis and streptococcal infections. Rheumatic fever is an inflammatory disease affecting connective tissue from several parts of human body as skin, heart, joints, etc and acute glomerulonephritis is a renal infection that affects mainly to renal glomerulus.

### METHOD PRINCIPLE

Latex particles coated with streptolysin O are agglutinated when mixed with samples containing ASO.

### REAGENTS

#### Package

	CORMAY ASO Latex 100
ASO-Latex	1 x 5 ml
ASO- Control (+) (red cap)	1 x 1 ml
ASO-Control (-) (blue cap)	1 x 1 ml
Stirrers	2 x 25 pcs.
Slides (6 circles each)	2 x 9 pcs.

### Reagent preparation and stability

The reagents are ready to use.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

### Concentrations in the test

latex particles coated with streptolysin O	50 mmol/l
human serum solution	150 mmol/l
animal serum solution	150 mmol/l
preservative	

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze.
- Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.
- False positive results may be obtained in conditions such as, rheumatoid arthritis, scarlet fever, tonsillitis.

- Children from 6 months to 2 years may cause false negative results.
- Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- ASO-Control (+) and ASO-Control (-) meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

ASO-Control (+) and ASO-Control (-) contain 2-methyl-2H-isothiazol-3-one.

### Warning

H317 May cause an allergic skin reaction.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P302+P352 IF ON SKIN: Wash with plenty of soap and water.

### ADDITIONAL EQUIPMENT

- mechanical rotator with adjustable speed at 80-100 r.p.m.
- general laboratory equipment.

### SPECIMEN

Serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

It is recommended to perform the assay with freshly collected samples.

### PROCEDURE

#### Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50  $\mu$ l of the sample and one drop of each positive and negative controls into separate circles on the slide test.
- Swirl the ASO-Latex reagent gently before using and add one drop (50  $\mu$ l) next to the sample to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a mechanical rotator at 80-100 r.p.m. **for 2 minutes**. False positive results could appear if the test is read later than two minutes.

### Reading and interpretation

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates an ASO concentration equal or greater than 200 IU/ml.

### Semi-quantitative method

- Make serial two fold dilutions of the sample in 0.9% NaCl solution.
- Proceed for each dilution as in the qualitative method.

### Reading and interpretation

The result (titer), in the semi-quantitative method, is defined as the highest dilution showing a positive result.

The approximate ASO concentration in the patient sample is calculated as follows:

ASO concentration (IU/ml) = 200 x the highest dilution

### REFERENCE VALUES

children < 5 years	< 100 IU/ml
adults	< 200 IU/ml

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

### PERFORMANCE CHARACTERISTICS

- Analytical sensitivity:** 200 IU/ml ( $\pm$ 50).
- Prozone effect:** no prozone effect up to 1500 IU/ml.
- Diagnostic sensitivity:** 98 %.
- Diagnostic specificity:** 97 %.

### Interferences:

Haemoglobin up to 10 g/l, bilirubin up to 20 mg/dl, triglycerides up to 10 g/l and RF up to 300 IU/ml do not interfere with the test.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Haffejee. Quarterly Journal of Medicine 1992. New series 84; 305: 641-658.
- Ahmed Samir et al. Pediatric Annals 1992; 21: 835-842.
- Spain J et al. Bull Wld Hlth Org 1961; 24: 271-279.
- The association of Clinical Pathologists 1961. Broadsheet 34.
- Picard B et al. La Presse Medicale 1983; 23: 2-6.
- Klein GC. Applied Microbiology 1971; 21: 999-1001.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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