CORMAY CRP LATEX

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LITERATURA

- Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10: 196-201.
- M.M. Pepys. The Lancet 1981; March 21: 653 – 656.
- Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139 – 144
- Yoshitsugy Hokama et al. Journal of Clinical Laboratory Status 1987; 1: 15 – 27.
- Yamamoto S et al. Veterinary Immunology and Immunopathology 1993; 36: 257 – 264.
- Charles Wadsworth et al. Clinica Chimica Acta; 1984: 138: 309 – 318.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

Data wydania: 09. 2023.

CORMAY CRP LATEX

Kit name	Kit size	(EN) Cat. No
CORMAY CRP Latex 100	100 tests	6-254

INTENDED USE

Agglutination test for detection of C-reactive protein (CRP) intended to use for the qualitative and semi-quantitative manual assays.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions

INTRODUCTION

C-Reactive protein (CRP) is an acute-phase α 1-globulin synthesized by hepatocytes. The level of CRP can rise up to more than 20 times in association with acute inflammation. High CRP levels can be detected within 24 hours of the onset of the acute phase and the peak serum concentration is achieved within 48 hours. After an inflammatory stimulus, the release of cytokines induces a variety of genes to help control infection. One class of such genes is that of carbohydrate-binding proteins to which pertains CRP. In the present of calcium ions, the protein binds carbohydrates present on the membrane of many bacteria.

The detection of CRP is useful in detecting inflammation and in monitoring its progress or response to treatment.

METHOD PRINCIPLE

Latex particles coated with goat IgG anti-human CRP are agglutinated when mixed with samples containing CRP.

REAGENTS

Package

	CORMAY CRP Latex 100
CRP-Latex	1 x 5 ml
CRP- Control (+) (red cap)	1 x 1 ml
CRP-Control (-) (blue cap)	1 x 1 ml
Stirrers	2 x 25 pcs.
Slides (6 circles each)	2 x 9 pcs.

Reagent preparation and stability

The reagents are ready to use.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the test

latex particles coated with goat IgG anti-human CRP	50 mmol/l
human serum solution	150 mmol/l
animal serum solution preservative	150 mmol/l

WARNINGS AND NOTES

- Do not freeze reagents.
- Protect before direct sunlight and contamination!
- Always keep vials in vertical position. If stored in different position gently mix to dissolve aggregates that may be present.
- Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.
- High CRP concentration samples may give negative results. Re-test the sample again using a drop of 20 μ l.
- Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- CRP-Latex and CRP-Control (+) and CRP-Control (-) meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

CRP-Latex contains boric acid.

Danger

H360FD May damage fertility. May damage the unborn child.

P201 Obtain special instructions before use.

P308 + P313 IF exposed or concerned: Get

medical advice.

P405 Store locked up.

Ingredients:

CRP-Control (+) and CRP-Control (-) contain 2-methyl-2H-isothiazol-3-one.

Warning

H317 May cause an allergic skin reaction..

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P302 + P352 IF ON SKIN: Wash with plenty

of soap and water.

CALIBRATION

The CORMAY CRP Latex sensitivity is calibrated to the Reference Material ERM-DA 474/IFCC.

ADDITIONAL EQUIPMENT

- mechanical rotator with adjustable speed at 80-100 r.p.m.
- general laboratory equipment.

SPECIMEN

Serum. Specimens are stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

It is recommended to perform the assay with freshly collected samples.

PROCEDURE

Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50 μ l of the sample and one drop of each positive and negative controls into separate circles on the slide test.
- Mix the CRP-Latex reagent vigorously before using and add one drop (50 μ l) next to the sample to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a mechanical rotator at 80-100 r.p.m. **for 2 minutes**. False positive results could appear if the test is read later than two minutes.

Reading and interpretation

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates a CRP concentration equal or greater than 6 mg/l.

Semi-quantitative method

- Make serial two fold dilutions of the sample in 0.9% NaCl solution.
- Proceed for each dilution as in the qualitative method.

Reading and interpretation

The result (titer), in the semi-quantitative method, is defined as the highest dilution showing a positive result.

The approximate CRP concentration in the patient sample is calculated as follows:

CRP concentration (mg/l) = 6 x the highest dilution

REFERENCE VALUES

adults	< 6 mg/l
--------	----------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

Any result different from negative control result will be considered as positive.

PERFORMANCE CHARACTERISTICS

- Analytical sensitivity:** 6 (5-10) mg/l.
- Prozone effect:** no prozone effect up to 1600 mg/l.
- Diagnostic sensitivity:** 95.6 %.
- Diagnostic specificity:** 96.2 %.

Interferences:

Haemoglobin up to 10 g/l, bilirubin up to 20 mg/dl, triglycerides up to 10 g/l and RF up to 100 U/ml do not interfere with the test.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10: 196-201.
- M.M. Pepys. The Lancet 1981; March 21: 653 – 656.
- Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139 – 144
- Yoshitsugu Hokama et al. Journal of Clinical Laboratory Status 1987; 1: 15 – 27.
- Yamamoto S et al. Veterinary Immunology and Immunopathology 1993; 36: 257 – 264.
- Charles Wadsworth et al. Clinica Chimica Acta; 1984: 138: 309 – 318.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

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- Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10: 196-201.
- M.M. Pepys. The Lancet 1981; March 21: 653 – 656.
- Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139 – 144
- Yoshitsugu Hokama et al. Journal of Clinical Laboratory Status 1987; 1: 15 – 27.
- Yamamoto S et al. Veterinary Immunology and Immunopathology 1993; 36: 257 – 264.
- Charles Wadsworth et al. Clinica Chimica Acta; 1984: 138: 309 – 318.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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