CORMAY TPHA

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INTENDED USE

CORMAY TPHA 100 test is intended to detect (qualitatively) and measure (semi-quantitatively) anti-*Treponema pallidum* antibodies in serum and plasma. It is used as an aid to diagnosis of syphilis.

CORMAY TPHA 100 test is intended to use manually. It is only for *in vitro* diagnostics, for healthcare professional users.

INTRODUCTION

Syphilis is a venereal disease caused by *Treponema pallidum* infection. A patient suffering this infection produces two kinds of antibodies: non-treponemal antibodies (reagins) and treponemal antibodies. The treponemal antibodies are specific antibodies against surface antigens of these microorganisms. Determination of treponemal antibodies is one of main syphilis tests.

METHOD PRINCIPLE

Stabilized avian erythrocytes sensitized with an antigenic *T. pallidum* solution, agglutinates in the presence of anti-*T. pallidum* antibodies.

REAGENTS

Package	CORMAY TPHA 100
TPHA-Test cells	1 x 7.5 ml
TPHA-Control cells	1 x 7.5 ml
TPHA-Diluent	2 x 10 ml
TPHA-Control (+)	1 x 1 ml
TPHA-Control (-)	1 x 1 ml

Reagent preparation and stability

The reagents and controls are ready to use.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the test

stabilized avian erythrocytes sensitized with <i>T. pallidum</i> antigens	1.5 mmol/l
stabilized avian erythrocytes	1.5 mmol/l
phosphate buffer	1.5 mmol/l
human serum solution	150 mmol/l
animal serum solution	150 mmol/l

Warnings and notes

- Do not freeze reagents.
- Store the vials in vertical position.
- Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

- The reagents contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.
- Shake vigorously the vials of both TPHA-Test cells and TPHA-Control cells immediately before use.
- During the assay keep the microplate away from the vibrations, heat and direct sunlight.
- Do not reuse the microplates.
- The CORMAY TPHA 100 test cannot discriminate antibodies anti-*T. pallidum* from antibodies to other pathogenic treponemas.
- False positive results have been described with samples of patients with mononucleosis, leprosy, borreliosis, autoimmune diseases and drug addiction.
- Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

ADDITIONAL EQUIPMENT

- U-well microtitration plates,
- general laboratory equipment.

SPECIMEN

Serum or plasma. Stable 7 days at 2-8°C, otherwise store frozen at below -20°C. Samples frozen at -20 °C or lower should be thawed and mixed prior to testing.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

It is recommended to perform the assay with freshly collected samples !

PROCEDURE

Reagent bottles should be shaken gently before use.

Before assay dilute serum 1:20 with TPHA-Diluent (50 µl serum + 950 µl TPHA-Diluent).

Controls do not require dilution.

Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Pipette into adjacent wells of a microtitration plate:

Sample 1:20 or Control	25 µl	25 µl
TPHA-Control cells	75 µl	----
TPHA-Test cells	----	75 µl
- Mix thoroughly the microplate till the complete homogenization of the mixing reaction.
- Cover the microplate and incubate at room temperature for 45-60 min. Do not move the microplate during incubation.
- Examine macroscopically the agglutination patterns of the cells.

Reading and interpretation

Read the results by comparing the agglutination patterns of the TPHA-Test cells with the TPHA-Control cells. Readings are scored and reported according to the following criteria:

Degree of hemagglutination	Reading	Result
Smooth mat of cells covering entire well bottom, sometimes with folded edges	4+	Positive
Smooth mat of cells covering part of the well bottom	3+	Positive
Smooth mat of cells surrounded by a red circle	2+	Positive
Smooth mat of cells covering less area and surrounded by a smaller red circle	1+	Positive
Button of cells having a small hole in centre	+/-	Borderline
Definite compact button of cells, sometimes with a very small hole in the centre	-	Negative

The negative control should not show any agglutination pattern with both TPHA-Test cells and TPHA-Control cells.

The positive control should only show agglutination patterns with TPHA-Test cells.

Any agglutination pattern showed by TPHA-Control cells indicates the presence of non-specific antibodies and cannot be interpreted.

Samples with a borderline pattern should be retested and reported as negatives if the same pattern is reproduced.

Reactive samples should be tittered following the semi-quantitative method.

Semi-quantitative method

- Make serial two fold dilutions of the prediluted 1:20 sample in TPHA-Diluent.
- Proceed for each dilution as in the qualitative method.

Reading and interpretation

The result (titer), in the semi-quantitative method, is defined as the highest dilution showing a positive result.

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity:

0.1 IU/ml against the 1st International Standard for human syphilitic plasma IgG and IgM NIBSC 05/132.

Diagnostic sensitivity: 100 %

Diagnostic specificity: 100 %.

Interferences:

Haemoglobin up to 10 g/l, bilirubin up to 20 mg/dl, triglycerides up to 10 g/l and RF up to 300 IU/ml do not interfere with the test.

WASTE MANAGEMENT

Please refer to local legal requirements.

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