FERRITIN

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INTENDED USE

Diagnostic kit for determination of ferritin concentration, intended to use in automatic analysers according to their user manual.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ferritin is an iron-containing protein with a molecular weight of approximately 450 kD. It is found mainly in the human liver and spleen, where its function is to eliminate and store iron in the body, and is also found in small amounts in human serum. This amount varies according to the movement of iron in the body, and hepatitis and malignant tumors, may be seen to increase due to cell destruction or tumor cell production, independent of iron reserves. Consequently, the measurement of ferritin is considered to be useful in the diagnosis, treatment, assessment of disease progression, and postoperative prognosis for such disease conditions.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ferritin in a sample and anti-ferritin antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ferritin in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	CORMAY FERRITIN	OS- FERRITIN	B50- FERRITIN
1-REAGENT	1 x 40.5 ml	1 x 35 ml	1 x 34.5 ml
2-REAGENT	1 x 24 ml	1 x 19 ml	1 x 18.5 ml

Reagent preparation and stability

Reagent ready to use. The reagents are stable up to the kit expiry date printed on the package when stored at 2-10°C. On board stability of the reagents depends on type of analyser used for analysis.

Concentrations in the test

suspension of latex particles sensitized with anti-ferritin (rabbit) antibodies (pH 7.3) ≤ 0.07 w/v%
 glycine buffer solution (pH 8.3)
 preservative

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- EUH210 Safety data sheet available on request.

ADDITIONAL EQUIPMENT

- automated clinical chemistry analyser capable of accommodating two-reagent assays;
- general laboratory equipment;

SPECIMEN

Serum.

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

Applications for analyzers are available on request.

REFERENCE VALUES ⁶

serum	ng/ml
male	20 – 250
female	10 – 120

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY FERRITIN CALBRATORS kit (Cat. No 4-491) is recommended. Calibrators and 0.9% NaCl should be used for calibration.

Calibration stability depends on type of analyser used for analysis. The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analysers Hitachi 912 and Hitachi 917. Results may vary if a different instrument is used.

- Sensitivity:** 9.2 ng/ml.

- Linearity:** up to 1200 ng/ml.

For higher concentration dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 0.98 g/dl, bilirubin up to 62 mg/dl, RF up to 520 IU/ml, triglycerides up to 500 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n=20	Mean [ng/dl]	SD [ng/dl]	CV [%]
level 1	32.35	0.96	2.98
level 2	1028.95	18.35	1.78
Reproducibility (day to day) n=21	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	16.47	0.87	5.31
level 2	105.18	1.60	1.52
level 3	428.71	3.52	0.82

Method comparison

A comparison between ferritin values determined at **Hitachi 912** (y) and at **Advia 1650** (x) using 28 samples gave following results:

$$y = 0.8682x - 0.9027 \text{ ng/ml};$$

$$R = 0.9950 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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