ACCENT-200 HbA1c DIRECT

Nr kat. 7-111 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of haemoglobin A1C concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

The determination of HbA1c is most commonly performed for the evaluation of glycemic control in diabetes mellitus. HbA1c values provide an indication of glucose levels over the preceding 4-8 weeks.

Throughout the circulatory life of the red cell, hemoglobin A_{1c} is formed continuously by the adduction of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al [1] showed hemoglobin A_{1c} in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that hemoglobin A_{1c} serve as an indicator of metabolic control of the diabetic, since hemoglobin A_{1c} levels approach normal values for diabetics in metabolic control [2,3,4].

Hemoglobin A_{1c} has been defined operationally as the "fast fraction" hemoglobins (HbA_{1a}, A_{1b}, A_{1c}) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA₀.

METHOD PRINCIPLE

Method for hemoglobin A_{1c} determination complies with standardized method certified by the National Glycohemoglobin Standardization Program (NGSP).

The present method utilizes the interaction of antigen and antibody to directly determine the HbA1c concentration in whole blood.

Total hemoglobin and HbA1c have the same unspecific absorption rate to latex particles. When mouse antihuman HbA1c monoclonal antibody is added, latex-HbA1c-mouse anti human HbA1c antibody complex is formed. Agglutination

is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA1c absorbed on to the surface of latex particles.

The amount of agglutination is measured as absorbance.

The HbA1c value is obtained from a calibration curve.

REAGENTS

Package

1-REAGENT	1 x 21 ml
2-REAGENT	1 x 7.7 ml
HEMOLYSING REAGENT	1 x 67.5 ml

The reagents (1-REAGENT, 2-REAGENT) stored at 2-8°C and HEMOLYSING REAGENT stored at 2-25°C are stable until expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

Concentrations in the test

latex	0.13%
mouse anti-human HbA1c monoclonal antibody	0.05 mg/ml
goat anti-mouse IgG polyclonal antibody stabilizers buffer	0.08 mg/dl

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin [6, 7, 8, 9].
- It has been reported that elevated levels of HbF may lead to underestimation of HbA1c and, that uremia does not interfere with HbA1c determination by immunoassay [10].
- This assay should not be used for the diagnosis of diabetes mellitus, but for monitor diabetic patients.
- In using Hemoglobin A_{1c} to monitor diabetic patients, results should be interpreted individually.
- Reagent HEMOLYSING REAGENT (Cat. No 4-398) can be ordered separately.
- Any clinical case with shortened erythrocyte survival (e.g. hemolytic anemia, blood loss, pregnancy) might cause decrease in HbA1c values.

SPECIMEN

Venous blood collected with EDTA.

Hemoglobin A_{1c} in whole blood collected with EDTA is stable for 7 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Sample pretreatment:

- Dispense 500 µl HEMOLYSING REAGENT into tubes labeled: Control, Patients, etc.
- Place 10 µl of well mixed whole blood into the appropriately labeled lyse reagent tube. Mix well and allow to stand for minimum 5 minutes, until complete lysis is evident. Next mix sample for 5 minutes.
- The treated sample may be stored up to 10 days at 2-8°C. Mix sample again for 5 minutes before measurement.
- Note:** calibrators and controls should be also hemolysed according to sample pretreatment.

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. 0.9% NaCl is recommended as a reagent blank.

Test result is read automatically and the value is reported in % of haemoglobin unit in accordance with NGSP standardization.

In order to convert the result reported in % haemoglobin (NGSP) to value reported in SI units mmol/mol in accordance with IFCC standardization, the following master equation should be used:

$$\text{HbA1c [mmol/mol IFCC]} = (\text{HbA1c [\% NGSP]} - 2.15) \times 10.929$$

REFERENCE VALUES ¹¹

Non-diabetes	< 6%
Patients with diabetes, control of glycaemia	< 7%

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY HbA1c DIRECT CONTROLS (Cat. No 4-328) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S ACCENT S120, ACCENT MC240, ACCENT M320, the CORMAY HbA1c DIRECT CALIBRATORS (Cat. No 4-308) are recommended. 0.9% NaCl should be used as a calibrator 0.

Controls and calibrators should be treated with HEMOLYSING REAGENT.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument is used.

- Analytical range:** 2 – 16% (up to 151 mmol/mol).

Specificity / Interferences

Bilirubin up to 50 mg/dl, triglycerides up to 2000 mg/dl, ascorbate up to 50 mg/dl, carbamylated Hb up to 7.5 mmol/l and acetylated Hb up to 5.0 mmol/l do not interfere with the test.

Precision

Repeatability (run to run)		Mean [%]	SD [%]	CV [%]
ACCENT-200 n=10	level 1	5.96	0.01	0.20
	level 2	11.20	0.29	2.56
ACCENT MC240 n=20	level 1	6.20	0.04	0.68
	level 2	12.36	0.31	2.52
Reproducibility (day to day)		Mean [%]	SD [%]	CV [%]
ACCENT-200 n=20	level 1	6.03	0.05	0.75
	level 2	12.28	0.16	1.29
ACCENT MC240 n=80	level 1	6.3	0.11	1.8
	level 2	13.2	0.38	2.8

Method comparison

A comparison between HbA1c values determined at **ACCENT 200** (y) and at **ADVIA 1650** (x) using 54 samples gave following results:

$$y = 0.864x + 0.727$$

$$R = 0.984 \quad (R - \text{correlation coefficient})$$

A comparison between HbA1c values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 60 samples gave following results:

$$y = 0.8436x + 1.0831$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Logit-Log 4P
Sensitivity	1
Replicates	2
Interval (day)	35
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	HbA1c	R1	180
No	61	R2	60
Full Name	HbA1c Direct	Sample Volume	3
Standard No	61	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	670 nm	Concentration	
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	0 38	Factor	
Incuba. Time	28	<input type="checkbox"/> Prozone check	
Unit	%	q1	<input type="checkbox"/>
Precision	0.01	q2	<input type="checkbox"/>
		q3	<input type="checkbox"/>
		q4	<input type="checkbox"/>
		PC	<input type="checkbox"/>
		Abs	<input type="checkbox"/>

Calibration Rule

Rule	Logit-Log 4P
Sensitivity	1
Replicates	1
Interval (day)	35
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	HbA1c	No.	061	Sample Type	OTHER
Chemistry	HbA1c DIRECT	Print name	HbA1c	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	670nm	Incubation Time	20	Reaction Time	31 32
Unit	%	Standard	3	Sample Vol	3
Blank Time		Aspirated		Diluent	
Decreased		Reagent Vol	180	R1	180
Increased		R2	60	R2	60
		Sample Blank		Auto Rerun	
Linearity range (Standard)		Linearity Limit		Substrate Depletion	
Linearity Range (Decreased)		Mixed Blank Abs	-40000 40000	On-board Stability	
Linearity Range (Increased)		Blank Response	-40000 40000	Reagent Alarm Limit	
R1 Blank Abs	-40000 40000	Twin Chemistry		Enzyme Linear Extension	
Blank Response	-40000 40000	<input type="checkbox"/> Prozone Check		Q1	
Q1		Q2		V1	
Q2		V1		Q3	
Q3		Q4		V2	
Q4		V2		Q5	
Q5		Q6		V3	
Q6		PC1		PC2	
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
Pretreat Sample Vol		Pretreat Sample Vol		Pretreat Sample Vol	
CALIBRATION SETTINGS			AUTO CALIBRATION		
Math model	Spline	Replicates	2	Bottle Changed	
Factor		Lot Changed		Cal Time	
ACCEPTANCE LIMITS					
Cal Time		Hour		Slope Diff	
Slope Diff		SD		Sensitivity	
Sensitivity		Repeatability	40000	Deter Coeff	
Deter Coeff					

ACCENT-200 HbA1c DIRECT

• ACCENT MC240

Chem HbA1c	No. 061	Sample Type OTHER
Chemistry HbA1c DIRECT	Print name HbA1c	
Reaction Type Endpoint	Reaction Direction positive	
Pri Wave 660nm	Sec Wave	
Unit %	Decimal 0.01	
Blank Time	Incubation Time 21	
	Reaction Time 31 32	
Standard 4 μL	Aspirated μL	Diluent μL
Decreased μL	Reagent Vol R1 180 μL	R2 60 μL
Increased μL		
<input type="checkbox"/> Sample Blank	<input type="checkbox"/> Auto Rerun	

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