## ACCENT-200 UREA

Nr kat. **7-206** (PL)

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Cat. No **7-206** (EN)

### INTENDED USE

Diagnostic kit for determination of urea concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400, ACCENT Neo200 and BS-120.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Urea is a product of amino acids catabolism. It is produced in liver and excreted in urine. Urea in the blood is reported as the blood urea nitrogen (BUN). Increased urea concentration in the serum, called uremia, is observed due to dehydration, renal failure, high-protein diet, increased protein catabolism caused by tissue injury or massive bleeding into the alimentary tract. The reason of reduced urea level could be overhydration, low-protein diet or starvation and severe liver disease.

### METHOD PRINCIPLE

Kinetic, enzymatic method with urease and glutamate dehydrogenase.

The rate of absorbance changing at  $\lambda=340$  nm is proportional to the urea concentration.

### REAGENTS

#### Package

1-REAGENT 2 x 30 ml  
2-REAGENT 1 x 15 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C is 12 weeks (ACCENT 200, ACCENT MC240, ACCENT S120).

### Concentrations in the reagent

#### 1-REAGENT

Tris (pH 7.8) ≤ 144 mmol/l  
ADP ≤ 0.84 mmol/l  
urease ≤ 250 µkat/l  
GLDH ≤ 10.5 µkat/l  
stabilizers, detergents, preservatives

#### 2-REAGENT

2-oxoglutarate ≤ 48.6 mmol/l  
NADH ≤ 1.6 mmol/l  
buffer, preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

- 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Warning

H319 Causes serious eye irritation  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

#### SPECIMEN <sup>9,10,11</sup>

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine.

Do not use heparine ammonium salt and fluoride as anticoagulants.

Specimen can be stored up to 7 days at 2-8°C.

**Urine preparation:** Samples with visible turbidity or the presence of precipitates should be pre-centrifuged.

Urine sample should be mixed well before analysis, diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Bacterial growth in the specimen may cause erroneously elevated results.

24-hours urine samples should be adjusted to pH < 7 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

Deionized water is recommended as a reagent blank.

#### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: UREA – FERRITIN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>8</sup>

serum / plasma	mg/dl	mmol/l
	< 50	< 8.3
24-hours urine	g/24h	mmol/24h
	20 – 35	300 – 550

1 mg of urea corresponds to 0.467 mg of urea nitrogen.

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionized water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT 200, ACCENT MC240) or every 2 weeks (ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-200 and/or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### LoB (Limit of Blank):

1.0 mg/dl (0.17 mmol/l)- ACCENT-200  
0.6 mg/dl (0.1 mmol/l)- ACCENT MC240

#### LoD (Limit of Detection):

1.4 mg/dl (0.23 mmol/l)- ACCENT-200  
0.9 mg/dl (0.15 mmol/l)- ACCENT MC240

#### LoQ (Limit of Quantitation):

4.5 mg/dl (0.75 mmol/l)- ACCENT-200  
2.2 mg/dl (0.37 mmol/l)- ACCENT MC240

#### Linearity

up to 220 mg/dl (36.52 mmol/l) - ACCENT-200  
up to 390 mg/dl (64.74 mmol/l) – ACCENT MC240

#### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	32.8	0.27	0.8
	level 2	100.4	0.50	0.5
ACCENT MC240 n=20	level 1	34.9	0.71	2.0
	level 2	101.7	1.58	1.6
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	33.2	0.72	2.2
	level 2	100.3	1.94	1.9
ACCENT MC240 n=80	level 1	34.1	0.63	1.9
	level 2	99.0	1.55	1.6

#### Method comparison

A comparison between urea values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 115 serum samples gave following results:  
y = 1.0337 x + 0.3901 mg/dl;  
R = 0.999 (R – correlation coefficient)

A comparison between urea values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 34 plasma samples gave following results:  
y = 1.0628 x - 1.2194 mg/dl;  
R = 1.000 (R – correlation coefficient)

A comparison between urea values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 30 urine samples gave following results:  
y = 1.0995 x - 100.05 mg/dl;  
R = 0.993 (R – correlation coefficient)

A comparison between urea values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 60 serum samples gave following results:  
y = 1.066 x - 0.9912 mg/dl;  
R = 1.000 (R – correlation coefficient)

A comparison between urea values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 30 plasma samples gave following results:  
y = 1.0547 x - 1.5399 mg/dl;  
R = 0.999 (R – correlation coefficient)

A comparison between urea values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 30 urine samples gave following results:  
y = 0.9388 x + 33.27 mg/dl;  
R = 0.997 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 UREA

### • ACCENT MC240

Chem	UREA	No.	010	Sample Type	SERUM
Chemistry	UREA	Print name	UREA	Reaction Direction	negative
Reaction Type	Kinetic	Pri Wave	340 nm	Sec Wave	450 nm
Unit	mg/dL	Decimal	0.1	Incubation Time	21
Blank Time		Reaction Time	1		5
Standard	3.6 μL	Aspirated		Diluent	
Decreased	3.6 μL		20		180
Increased					
		Reagent Vol		R1	180 μL
				R2	45 μL
		Sample Blank	V	Auto Rerun	

Linearity range (Standard)	2.2	390	Linearity Limit	0.2			
Linearity Range (Decreased)			Substrate Depletion				
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000			
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)			
Blank Response	-35000	35000	Reagent Alarm Limit				
Twin Chemistry			Enzyme Linear Extension				
		Prozone Check					
Q1		Q2	V1	Q3	Q4	V2	
Q5		Q6	V3	PC1	PC2		
		Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
		Pretreat Sample Vol		μL	Pretreat Sample Vol		μL

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear		Bottle Changed
Factor		Replicates	2
			Lot Changed
			Cal Time

<b>ACCEPTANCE LIMITS</b>			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

### • ACCENT M320

Chem	UREA	No.	010	Sample Type	SERUM
Chemistry	UREA	Print name	UREA	Reaction Direction	negative
Reaction Type	Kinetic	Pri Wave	340 nm	Sec Wave	700 nm
Unit	mg/dL	Decimal	0.1	Incubation Time	18
Blank Time		Reaction Time	3		11
Standard	2.4 μL	Aspirated		Diluent	
Decreased	2.4 μL		20		180
Increased					
		Reagent Vol		R1	160 μL
				R2	40 μL
		Sample Blank	V	Auto Rerun	

Linearity range (Standard)	1.5	350	Linearity Limit	0.2			
Linearity Range (Decreased)			Substrate Depletion				
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000			
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)			
Blank Response	-35000	35000	Reagent Alarm Limit				
Twin Chemistry			Enzyme Linear Extension				
		Prozone Check					
O1		O2	V1	O3	O4	V2	
Q5		Q6	V3	PC1	PC2		
		Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
		Pretreat Sample Vol		μL	Pretreat Sample Vol		μL

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear		Bottle Changed
Factor		Replicates	2
			Lot Changed
			Cal Time

<b>ACCEPTANCE LIMITS</b>			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

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