ACCENT-200 UA

Nr kat. 7-208

(PL)

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Cat. No **7-208** (EN)

INTENDED USE

Diagnostic kit with ascorbate oxidase for determination of uric acid concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, BS-120 / BS-130, ACCENT S120, ACCENT MC240 and ACCENT M320.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

REAGENTS

Package

1-Reagent 2 x 30 ml
2-Reagent 1 x 15 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the reagent

1-Reagent

ascorbate oxidase ≤ 104 µkat/l
peroxidase (POD) ≤ 22.4 µkat/l
4-aminoantipyrine ≤ 1.2 mmol/l
sodium hydroxide ≤ 0.8 %
buffer PIPES (pH 7.0) ≤ 120 mmol/l
stabilizers, preservatives, detergent

2-Reagent

buffer PIPES (pH 7.0) ≤ 60 mmol/l
ADPS ≤ 2 mmol/l
uricase ≤ 9.9 µkat/l
ferricyanide potassium ≤ 22.8 µmol/l
sodium hydroxide ≤ 0.4 %
stabilizers, preservatives, detergent

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H315 Causes skin irritation.
H319 Causes serious eye irritation.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

24- hours urine, serum, heparinized plasma free from hemolysis.
Do not use EDTA, fluoride and oxalate as anticoagulants!
Urine preparation: To prevent precipitation of salts of uric acid, 10 ml of NaOH (500 g/L) should be added to the collection bottle before collection of a 24-hour specimen. Urine should be diluted with distilled water in the ratio of 1 to 4 (multiply the result by 5).
Serum and plasma can be stored 3-5 days at 2-8°C or 6 months at -20°C. 24-hours urine samples can be stored approximately 3 days at room temperature.
Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S there is a probability of **cross-contamination** affecting the tests results: ALBUMIN – UA, CREATININE – UA, HDL DIRECT II GEN – UA, LDL DIRECT II GEN – UA, URINE PROTEINS II GEN – UA. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁵

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458
24-hours urine	mg/24h	mmol/24h
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) - for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) - for determination in urine. For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0. The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.04 mg/dl (2.38 µmol/l) – ACCENT-200
0.01 mg/dl (0.59 µmol/l) – ACCENT MC240

LoD (Limit of Detection):

0.07 mg/dl (4.16 µmol/l) – ACCENT-200
0.03 mg/dl (1.78 µmol/l) – ACCENT MC240

LoQ (Limit of Quantitation):

0.30 mg/dl (17.84 µmol/l) – ACCENT-200 (serum/ plasma)
0.36 mg/dl (21.41 µmol/l) – ACCENT-200 (urine)
0.15 mg/dl (8.92 µmol/l) – ACCENT MC240 (serum/ plasma)
0.28 mg/dl (16.65 µmol/l) – ACCENT MC240 (urine)

Linearity:

up to 33 mg/dl (1962.84 µmol/l) – ACCENT-200 (serum/ plasma)
up to 56 mg/dl (3330.88 µmol/l) – ACCENT-200 (urine)
up to 43 mg/dl (2557.64 µmol/l) – ACCENT MC240 (serum/ plasma)
up to 57 mg/dl (3390.36 µmol/l) – ACCENT MC240 (urine)

For higher concentration of uric acid in serum or plasma, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 30 mg/dl for determinations in serum, ascorbate up to 50 mg/dl for determinations in urine, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	4.95	0.06	1.15
	level 2	9.28	0.06	0.66
ACCENT MC240 n=20	level 1	4.87	0.07	1.49
	level 2	9.04	0.09	1.03
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	4.84	0.11	2.3
	level 2	9.52	0.18	1.9
ACCENT MC240 n=80	level 1	4.94	0.09	1.8
	level 2	9.72	0.11	1.1

Method comparison

A comparison between uric acid values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:
y = 0.9678 x + 0.3409 mg/dl;
R = 0.998 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 35 plasma samples gave following results:
y = 0.9718 x + 0.292 mg/dl;
R = 0.994 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:
y = 1.0354 x - 0.6126 mg/dl;
R = 0.997 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:
y = 0.9657 x + 0.424 mg/dl;
R = 0.998 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 35 plasma samples gave following results:
y = 0.9614 x + 0.352 mg/dl;
R = 0.994 (R – correlation coefficient)

A comparison between uric acid values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:
y = 1.0354 x - 0.7314 mg/dl;
R = 0.996 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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ACCENT-200 UA

• ACCENT MC240

Chem	UA	No.	009	Sample Type	SERUM/URINE
Chemistry	URIC ACID	Print name	UA	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	660 nm	Decimal	0.01
Pri Wave	546 nm	Incubation Time	21	Reaction Time	19 21
Unit	mg/dL	Blank Time	-3 -1	Standard	Sample Vol 3.2 μL Aspirated 20 μL Diluent 180 μL
				Decreased	3.2 μL 20 μL 180 μL
				Increased	μL μL μL
					<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun

Linearity range (Standard)	0.15	43	Linearity Limit					
Linearity Range (Decreased)			Substrate Depletion					
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000				
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)				
Blank Response	-35000	35000	Reagent Alarm Limit					
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension					
<input type="checkbox"/> Prozone Check								
O1		O2	V1	O3	O4	V2		
O5		O6	V3	PC1	PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment				
		Pretreat Sample Vol	μL	Pretreat Sample Vol	μL			

CALIBRATION SETTINGS

Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor	Replicates 2	<input type="checkbox"/> Lot Changed
		<input type="checkbox"/> Cal Time

AUTO CALIBRATION

Bottle Changed
 Lot Changed
 Cal Time

ACCEPTANCE LIMITS

Cal Time	Hour	
Slope Diff		SD
Sensitivity		Repeatability 35000
Deter Coeff		

• ACCENT M320

Chem	UA	No.	009	Sample Type	SERUM/URINE
Chemistry	URIC ACID	Print name	UA	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700 nm	Decimal	0.1
Pri Wave	546 nm	Incubation Time	20	Reaction Time	15 17
Unit	mg/dL	Blank Time	-3 -1	Standard	Sample Vol 2.4 μL Aspirated 20 μL Diluent 180 μL
				Decreased	2.4 μL 20 μL 180 μL
				Increased	μL μL μL
					<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun

Linearity range (Standard)	0.08	46	Linearity Limit					
Linearity Range (Decreased)			Substrate Depletion					
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000				
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)				
Blank Response	-35000	35000	Reagent Alarm Limit					
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension					
<input type="checkbox"/> Prozone Check								
O1		O2	V1	O3	O4	V2		
O5		O6	V3	PC1	PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment				
		Pretreat Sample Vol	μL	Pretreat Sample Vol	μL			

CALIBRATION SETTINGS

Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor	Replicates 2	<input type="checkbox"/> Lot Changed
		<input type="checkbox"/> Cal Time

AUTO CALIBRATION

Bottle Changed
 Lot Changed
 Cal Time

ACCEPTANCE LIMITS

Cal Time	Hour	
Slope Diff		SD
Sensitivity		Repeatability 35000
Deter Coeff		

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).