ACCENT-200 COMPLEMENT C3

Nr kat. 7-211

(PL)

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Cat. No 7-211 (EN)

INTENDED USE

Diagnostic kit for determination of complement C3 concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Complement is a group of 20 immunologically distinct proteins present in blood and tissues. They are able to interact with antigen-antibody complexes, with each other and with cell membranes, in a complex way intended at destroying viruses and bacteria. They are synthesised in liver and are present in serum as functionally inactive molecules. They are activated by antigen-antibody complexes. C3 complement is a α -glycoprotein of 2 subunits. It is an acute phase reactant whose levels are increased during the acute phase. Low levels are found in immune complex diseases and in inherited deficiency which results in recurrent infections.

METHOD PRINCIPLE

The complement C3, present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to complement C3 concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 35 ml
 2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: Imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human complement C3 antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found

to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma.

Serum should be separated from red blood cells as soon as possible after blood collection. If the test cannot be done immediately, the sample should be stored at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ²

adults	0.9 – 1.8 g/l
--------	---------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320, the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzers Hitachi or ACCENT MC240. Results may vary if a different instrument is used.

- LOQ:** 0.015 g/l

- Linearity:** up to 4.31 g/l

- Specificity / Interferences**

Haemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.759	0.005	0.66
level 2	1.615	0.013	0.80
level 3	2.367	0.021	0.87

Reproducibility (day to day) n = 59	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.778	0.019	2.42
level 2	1.634	0.026	1.59
level 3	2.428	0.045	1.87

- Method comparison**

A comparison between C3 complement concentration determined at **ACCENT MC240** (y) and at **BS-400** (x) using 54 serum samples gave following results:

$$y = 1.034 x - 0.204 \text{ g/l;}$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT MC240

Chem	C3	No.	043	Sample Type	SERUM
Chemistry	COMPLEMENT C3	Reaction Type	Endpoint	Reaction Direction	positive
Pri Wave	340 nm	Unit	g/L	Sec Wave	700 nm
Blank Time	-3		-1	Decimal	0.001
Standard	4,5	Sample Vol	4,5	Aspirated	20
Decreased	4,5	Diluent	180	Diluent	180
Increased		Reagent Vol	32	R1	160
		R2	32	R2	32
		Sample Blank	V	Auto Rerun	
Linearity range (Standard)	0,007	2,40	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000	35000
R1 Blank Abs	-35000	35000	On-board Stability		Day(s)
Blank Response	-35000	35000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
	<input type="checkbox"/> Prozone Check				
Q1		Q2		V1	
Q5		Q6		V3	
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		PC1	
				PC2	
		Pretreat Sample Vol		μL	
		Pretreat Sample Vol		μL	
CALIBRATION SETTINGS	Math model	Spline	AUTO CALIBRATION		
Factor		Replicates	2	<input type="checkbox"/> Bottle Changed	
				<input type="checkbox"/> Lot Changed	
				<input type="checkbox"/> Cal Time	
ACCEPTANCE LIMITS	Cal Time		Hour		
Slope Diff		SD			
Sensitivity		Repeatability	35000		
Deter Coeff					

ACCENT-200 COMPLEMENT C3

• ACCENT M320

Chem <input type="text" value="C3"/>	No. <input type="text" value="043"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="COMPLEMENT C3"/>	Print name <input type="text" value="C3"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="700 nm"/>	
Unit <input type="text" value="g/L"/>	Decimal <input type="text" value="0.001"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="25"/>	Reaction Time <input type="text" value="13"/> <input type="text" value="15"/>
Standard <input type="text" value="4"/> <input type="text" value="4"/> <input type="text" value="20"/> <input type="text" value="180"/>	Aspirated <input type="text" value="4"/> <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="4"/> <input type="text" value="20"/> <input type="text" value="180"/>
Reagent Vol <input type="text" value="160"/> <input type="text" value="32"/>	R1 <input type="text" value="160"/> <input type="text" value="32"/>	R2 <input type="text" value="160"/> <input type="text" value="32"/>
<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun		
Linearity range (Standard) <input type="text" value="0,007"/> <input type="text" value="2,45"/>	Linearity Limit <input type="text"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)	
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check		
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/> <input type="text"/>		
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/> <input type="text"/>		
<input type="checkbox"/> Sample Pretreatment <input type="checkbox"/> Control Pretreatment <input type="checkbox"/> Calibrator Pretreatment		
<input type="text"/>	Pretreat Sample Vol <input type="text"/> <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> <input type="text"/> μ L
CALIBRATION SETTINGS	AUTO CALIBRATION	
Math model <input type="text" value="Spline"/>	<input type="checkbox"/> Bottle Changed	
Factor <input type="text"/> Replicates <input type="text" value="2"/>	<input type="checkbox"/> Lot Changed	
	<input type="checkbox"/> Cal Time	
ACCEPTANCE LIMITS		
Cal Time <input type="text"/> Hour		
Slope Diff <input type="text"/> SD <input type="text"/>		
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>		
Deter Coeff <input type="text"/>		

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