ACCENT-200 ALP

Nr kat. 7-212 (PL)

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Cat. No. **7-212** (EN) Zn²⁺ 1.25 mmol/l

- INTENDED USE**


ALP reagent is intended to determine quantitatively alkaline phosphatase activity in serum and plasma. It is used to monitoring and as an aid to diagnosis of clinical conditions associated with abnormal ALP activity. ALP reagent is intended to use on automatic analyzers ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180, ACCENT S120 / BS-230, ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E, BS-120, ACCENT 400 and ACCENT Neo200. It is only for *in vitro* diagnosis, for healthcare professional users.

- SUMMARY¹⁻⁶**

Determinations of alkaline phosphatase activity in serum and plasma are mainly used as an aid in diagnosis, differentiation and monitoring of hepatobiliary disease and bone diseases. Elevated levels of the enzyme in the blood are observed, e.g. in cholestatic liver diseases, osteoporosis or Paget's disease. Decreased ALP levels are associated with Wilson's disease, Hypophosphatasia, and deficiency of magnesium, zinc and vitamin B12. Measurement of ALP activity is also used as a predictive marker in outcomes (death or liver transplantation) of patients with primary sclerosing cholangitis (PSC) and primary biliary cholangitis (PBC). The determination of the level of alkaline phosphatase is also used to monitoring the effectiveness of the PBC and IgG4-related cholangiopathy treatment.

- METHOD PRINCIPLE^{7,8}**

Kinetic method recommended by International Federation of Clinical Chemistry (IFCC).

The rate of 4-nitrophenol formation is directly proportional to the ALP activity.

- REAGENTS**

Package	
1-REAGENT	1 x 30 ml
2-REAGENT	1 x 8,5 ml

- Number of tests:**

ACCENT-200 (II GEN)	130
ACCENT-220S	130
ACCENT S120	160
ACCENT MC240	220
ACCENT M320	220
BS-120	130

- CONCENTRATIONS OF THE ACTIVE INGREDIENTS IN THE REAGENT**

1-REAGENT	
2-amino-2-methyl-1-propanol (AMP)	425 mmol/l
Mg ²⁺	2.5 mmol/l

metal ion chelator	
pH adjuster preservative	
2-REAGENT	
p-nitrophenylphosphate	83.9 mmol/l
2-amino-2-methyl-1-propanol (AMP)	50 mmol/l
pH adjuster preservative	

- REAGENT STABILITY**

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 11 weeks (ACCENT-200 (II GEN)) or 12 weeks (ACCENT MC240, ACCENT S120).

- WARNINGS AND NOTES**

- Protect from direct sunlight and avoid contamination!
- During the reaction p-nitrophenol is produced. Do not swallow or inhale, avoid contact with skin.
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

- Ingredients:**

1-REAGENT and 2-REAGENT contain 2-methyl-2H-isothiazol-3-one.

- Warning**

H317 May cause an allergic skin reaction.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

- SPECIMEN¹**

Serum, heparinized plasma free from hemolysis. Do not use EDTA, citrate and oxalate as anticoagulants because of ALP activity inhibition!

ALP activity remains stable in specimen up to 4 hours at 15-25°C. Freezing of sample causes a decrease of enzyme activity. Frozen specimens should be thawed and kept at room temperature for 18 to 24 hours before measurement to achieve full enzyme reactivation.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

Reagents are ready to use.

For reagent blank deionized water is recommended.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: ALP – MICROALBUMIN, ALP – MG. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

QUALITY CONTROL

Pay attention to preparation of calibrator and control before ALP determination.

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every week (ACCENT-200 (II GEN)), every 2 weeks (ACCENT MC240) or every 3 weeks (ACCENT S120).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ⁹⁻¹¹

Gender	Age	U/l (37°C)	µkat/l (37°C)
female	1 – 30 days	48 – 406	0.80 – 6.77
	31 days – 1 year	124 – 341	2.07 – 5.68
	1 year – 3 years	108 – 317	1.80 – 5.28
	4 – 15 years	54 – 369	0.91 – 6.23
	16 – 18 years	35 – 124	0.58 – 2.07
	19 – 20 years	39 – 118	0.65 – 1.97
	20 – 50 years	42 – 98	0.71 – 1.67
	50 – 60 years	39 – 118	0.65 – 1.97
≥ 60 years	53 – 141	0.90 – 2.40	

Male	1 – 30 days	75 – 316	1.25 – 5.27
	31 days – 1 year	82 – 383	1.37 – 6.38
	1 year – 3 years	104 – 345	1.73 – 5.75
	4 – 15 years	54 – 369	0.91 – 6.23
	16 – 18 years	58 – 331	0.97 – 5.52
	19 – 20 years	41 – 137	0.68 – 2.28
	20 – 50 years	53 – 128	0.90 – 2.18
	50 – 60 years	41 – 137	0.68 – 2.28
≥ 60 years	56 – 119	0.95 – 2.02	

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 (II GEN) and/or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity**
3 U/l (0.05 µkat/l) - ACCENT 200 (II GEN)
- LoQ (Limit of Quantitation)**
6.4 U/l (0.11 µkat/l) – ACCENT MC240
- Linearity**
up to 1000 U/l (16.7 µkat/l) - ACCENT 200 (II GEN)
up to 1050 U/l (17.5 µkat/l) – ACCENT MC240

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

- Measurement range**
3 U/l (0.05 µkat/l) - 1000 U/l (16.7 µkat/l) – ACCENT 200 (II GEN)
6.4 U/l (0.11 µkat/l) – 1050 U/l (17.5 µkat/l) – ACCENT MC240
- Specificity / Interferences**
Haemoglobin up to 0.625 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision (ACCENT MC240)

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	105.1	0.49	0.5
level 2	449.7	1.21	0.3
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	101.6	4.13	4.1
level 2	437.9	16.10	3.7

Method comparison

A comparison between ALP values determined at **ACCENT-200 (II GNE)** (y) and at **COBAS INTEGRA 400** (x) using 33 samples gave following results:

$$y = 0.8703 x + 9.6303 \text{ U/l};$$
$$R = 0.988 \quad (R - \text{correlation coefficient})$$

A comparison between ALP values determined at **ACCENT MC240** (y) and at **BS-800** (x) using 61 serum samples gave following results:

$$y = 0.9072 x + 14.442 \text{ U/l};$$
$$R = 0.996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT ¹²

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹³

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

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- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2290 (2006).
- Soldin S. J., Brugnara C., Wong E. C.: Pediatric reference Ranges, 4th ed, AACCPress, 10-11 (2003).
- Hay, W.W., et al, Current pediatric diagnosis and treatment (15th ed.). New York: Lange MedicalBooks/McGraw Hill. (2000).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 08	Current version: 09
Sections updated: WARNING AND NOTES; WASTE MANAGEMENT	

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
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 417, 1393 (2018).
- European Association for the Study of the Liver, EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis, Journal of Hepatology, 67, 145-172, (2017).
- Hartleb M. i wsp., Choroby cholestatyczne u doroslych - wytyczne postepowania Sekcji Hepatologicznej Polskiego Towarzystwa Gastroenterologii, Gastroenterologia Praktyczna, vol. 2, 5-31 (2018).
- Alvarez L, et al., Usefulness of Biochemical Markers of Bone Turnover in Assessing Response to the Treatment of Paget's Disease, Bone, 29 (5), 441-452 (2001).
- Kravets I., Paget's Disease of Bone: Diagnosis and Treatment, The American Journal of Medicine, 131 (11), 1298-1303, (2018)
- Lammers W.J., et al., Levels of Alkaline Phosphatase and Bilirubin Are Surrogate End Points of Outcomes of Patients With Primary Biliary Cirrhosis: An International Follow-up Study, Gastroenterology, 1:47, 1338-1349, (2014)
- Bessey O.A., Lowry O.H. And Brock M.J., A method for the rapid determination of alkaline phosphates with five cubic millimeters of serum, Biol. Chem., 164, 321-9 (1946).
- Bowers G.N.Jr. and McComb R.B., A continuous spectrophotometric method for measuring the activity of serum alkaline phosphatase, Clin Chem., 12(2), 70-89 (1966)
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2290 (2006)
- Soldin S. J., Brugnara C., Wong E. C.: Pediatric reference Ranges, 4th ed, AACC Press, 10-11 (2003).
- Hay, W.W., et al, Current pediatric diagnosis and treatment (15th ed.). New York: Lange MedicalBooks/McGraw Hill. (2000).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

ACCENT-220S

Parameters

Test	ALP	R1	200
No	18	R2	50
Full Name	Alkaline phosphatase	Sample Volume	5
Standard No	18	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	7.8 1300
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
React. Time	2 13	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

BS-120

Parameters

Test	ALP	R1	200
No	18	R2	50
Full Name	Alkaline phosphatase	Sample Volume	5
Standard No	18	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	3.2 950
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
React. Time	3 11	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem No. Sample Type
 Chemistry Print name
 Reaction Type Reaction Direction
 Pri Wave Sec Wave
 Unit Decimal
 Incubation Time
 Blank Time Reaction Time

Standard	<input type="text" value="6.4"/> μ L	Aspirated	<input type="text"/> μ L	Diluent	<input type="text"/> μ L	Reagent Vol	<input type="text"/> μ L
Decreased	<input type="text" value="6.4"/> μ L	<input type="text" value="20"/> μ L	<input type="text" value="180"/> μ L	R1	<input type="text" value="160"/> μ L	R2	<input type="text" value="40"/> μ L
Increased	<input type="text"/> μ L	<input type="text"/> μ L	<input type="text"/> μ L				

 Sample Blank Auto Rerun

Linearity range (Standard) Linearity Limit
 Linearity Range (Decreased) Substrate Depletion
 Linearity Range (Increased) Mixed Blank Abs
 R1 Blank Abs On-board Stability Day(s)
 Blank Response Reagent Alarm Limit
 Twin Chemistry Enzyme Linear Extension
 Prozone Check
 Q1 Q2 V1 Q3 Q4 V2
 Q5 Q6 V3 PC1 PC2
 Sample Pretreatment Control Pretreatment Calibrator Pretreatment
 Pretreat Sample Vol μ L Pretreat Sample Vol μ L

CALIBRATION SETTINGS

Math model
 Factor Replicates

AUTO CALIBRATION

Bottle Changed
 Lot Changed
 Cal Time

ACCEPTANCE LIMITS

Cal Time Hour
 Slope Diff SD
 Sensitivity Repeatability
 Deter Coeff

• ACCENT MC240

Chem No. Sample Type
 Chemistry Print name
 Reaction Type Reaction Direction
 Pri Wave Sec Wave
 Unit Decimal
 Incubation Time
 Blank Time Reaction Time

Standard	<input type="text" value="5.4"/> μ L	Aspirated	<input type="text"/> μ L	Diluent	<input type="text"/> μ L	Reagent Vol	<input type="text"/> μ L
Decreased	<input type="text" value="5.4"/> μ L	<input type="text" value="20"/> μ L	<input type="text" value="180"/> μ L	R1	<input type="text" value="120"/> μ L	R2	<input type="text" value="30"/> μ L
Increased	<input type="text"/> μ L	<input type="text"/> μ L	<input type="text"/> μ L				

 Sample Blank Auto Rerun

Linearity range (Standard) Linearity Limit
 Linearity Range (Decreased) Substrate Depletion
 Linearity Range (Increased) Mixed Blank Abs
 R1 Blank Abs On-board Stability Day(s)
 Blank Response Reagent Alarm Limit
 Twin Chemistry Enzyme Linear Extension
 Prozone Check
 Q1 Q2 V1 Q3 Q4 V2
 Q5 Q6 V3 PC1 PC2
 Sample Pretreatment Control Pretreatment Calibrator Pretreatment
 Pretreat Sample Vol μ L Pretreat Sample Vol μ L

CALIBRATION SETTINGS

Math model
 Factor Replicates

AUTO CALIBRATION

Bottle Changed
 Lot Changed
 Cal Time

ACCEPTANCE LIMITS

Cal Time Hour
 Slope Diff SD
 Sensitivity Repeatability
 Deter Coeff

• ACCENT M320

Chem	<input type="text" value="ALP"/>	No.	<input type="text" value="18"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="ALP"/>			Print name	<input type="text" value="ALP"/>
Reaction Type	<input type="text" value="Kinetic"/>			Reaction Direction	<input type="text" value="positive"/>
Pri Wave	<input type="text" value="412 nm"/>			Sec Wave	<input type="text" value="660 nm"/>
Unit	<input type="text" value="U/L"/>			Decimal	<input type="text" value="0.1"/>
Blank Time	<input type="text"/>	<input type="text"/>		Incubation Time	<input type="text" value="20"/>
				Reaction Time	<input type="text" value="5"/> <input type="text" value="15"/>
Standard	<input type="text" value="5.4"/> <input type="text" value="μL"/>	Aspirated	<input type="text" value="20"/> <input type="text" value="μL"/>	Diluent	<input type="text" value="180"/> <input type="text" value="μL"/>
Decreased	<input type="text" value="5.4"/> <input type="text" value="μL"/>			Reagent Vol	
Increased	<input type="text" value="μL"/>			R1	<input type="text" value="120"/> <input type="text" value="μL"/>
				R2	<input type="text" value="30"/> <input type="text" value="μL"/>
	<input type="text"/>	Sample Blank	<input checked="" type="checkbox"/>	Auto Rerun	<input type="checkbox"/>

Linearity range (Standard)	<input type="text" value="13"/>	<input type="text" value="1360"/>	Linearity Limit	<input type="text" value="0.2"/>							
Linearity Range (Decreased)	<input type="text"/>	<input type="text"/>	Substrate Depletion	<input type="text" value="35000"/>							
Linearity Range (Increased)	<input type="text"/>	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>							
R1 Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	On-board Stability	<input type="text" value="84"/> Day(s)							
Blank Response	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>							
Twin Chemistry	<input type="text"/>		<input type="checkbox"/> Enzyme Linear Extension								
	<input type="checkbox"/> Prozone Check										
Q1	<input type="text"/>	Q2	<input type="text"/>	V1	<input type="text"/>	Q3	<input type="text"/>	Q4	<input type="text"/>	V2	<input type="text"/>
Q5	<input type="text"/>	Q6	<input type="text"/>	V3	<input type="text"/>	PC1	<input type="text"/>	PC2	<input type="text"/>		<input type="text"/>
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
	<input type="text"/>	Pretreat Sample Vol	<input type="text"/>	μL		Pretreat Sample Vol	<input type="text"/>	μL			

CALIBRATION SETTINGS

Math model	<input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor	<input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates	<input type="text" value="2"/>	<input type="checkbox"/> Cal Time

AUTO CALIBRATION

Bottle Changed

Lot Changed

Cal Time

ACCEPTANCE LIMITS

Cal Time	<input type="text" value="336"/> Hour
Slope Diff	<input type="text"/>
Sensitivity	<input type="text"/>
Deter Coeff	<input type="text"/>
SD	<input type="text"/>
Repeatability	<input type="text" value="35000"/>

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