## ACCENT-200 COMPLEMENT C4

Nr kat. 7-213

(PL)

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Cat. No 7-213 (EN)

### INTENDED USE

Diagnostic kit for determination of complement C4 concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Complement is a group of 20 immunologically distinct proteins present in blood and tissues. They are able to interact with antigen-antibody complexes, with each other and with cell membranes, in a complex way intended at destroying viruses and bacteria. They are synthesised in liver and are present in serum as functionally inactive molecules. They are activated by antigen-antibody complexes. C4 complement is a  $\alpha$ -glycoprotein of 3 subunits.

It is an acute phase reactant whose levels are increased during the acute phase. Low levels are found in immune complex diseases and in inherited angioedema, while C3 complement levels are normal.

### METHOD PRINCIPLE

The complement C4 present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at  $\lambda=340$  nm is proportional to complement C4 concentration in the sample.

### REAGENTS

#### Package

1-Reagent 1 x 35 ml  
2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

#### Reagent components

**1-Reagent:** imidazole buffer, polymer, inorganic salt, buffer, preservative.

**2-Reagent:** anti-human complement C4 antiserum, inorganic salt, buffer, preservative.

#### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma.

Serum should be separated from red blood cells as soon as possible after blood collection. If the test cannot be done immediately, the sample should be stored at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES <sup>2</sup>

adults	0.1 – 0.4 g/l
--------	---------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320 the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzers Hitachi or ACCENT MC240.

Results may vary if a different instrument is used.

- LOQ:** 0.002 g/l
- Linearity:** up to 0.825 g/l.
- Specificity / Interferences**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl triglycerides up to 155.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

### Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.147	0.001	0.77
level 2	0.321	0.003	0.76
level 3	0.462	0.003	0.64
Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.137	0.005	3.48
level 2	0.321	0.008	2.52
level 3	0.472	0.010	2.08

### Method comparison

A comparison between C4 complement concentration determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 66 serum samples gave following results:

$$y = 0.9901x + 0.0093 \text{ g/l};$$

$$R = 0.978 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 COMPLEMENT C4

• ACCENT M320

Chem <input type="text" value="C4"/>	No. <input type="text" value="044"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="COMPLEMENT C4"/>	Print name <input type="text" value="C4"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="700 nm"/>	
Unit <input type="text" value="g/L"/>	Decimal <input type="text" value="0.001"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="25"/>	Reaction Time <input type="text" value="26"/> <input type="text" value="28"/>
Standard <input type="text" value="6"/> <input type="text" value="6"/> <input type="text" value="6"/>	Aspirated <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="6"/> <input type="text" value="20"/>	Reagent Vol R1 <input type="text" value="160"/> <input type="text" value="32"/>	Reagent Vol R2 <input type="text" value="32"/>
Increased <input type="text" value="6"/> <input type="text" value="20"/>	Sample Blank <input type="text" value="V"/>	Auto Rerun <input type="text" value=""/>

Linearity range (Standard) <input type="text" value="0,005"/> <input type="text" value="0,68"/>	Linearity Limit <input type="text" value=""/>
Linearity Range (Decreased) <input type="text" value=""/>	Substrate Depletion <input type="text" value=""/>
Linearity Range (Increased) <input type="text" value=""/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text" value=""/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text" value=""/>
Twin Chemistry <input type="text" value=""/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text" value=""/>	Q2 <input type="text" value=""/>	V1 <input type="text" value=""/>	Q3 <input type="text" value=""/>	Q4 <input type="text" value=""/>	V2 <input type="text" value=""/>
Q5 <input type="text" value=""/>	Q6 <input type="text" value=""/>	V3 <input type="text" value=""/>	PC1 <input type="text" value=""/>	PC2 <input type="text" value=""/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="text" value=""/>	<input type="checkbox"/> Calibrator Pretreatment	<input type="text" value=""/>	
<input type="text" value=""/>	Pretreat Sample Vol <input type="text" value=""/>	<input type="text" value=""/>	Pretreat Sample Vol <input type="text" value=""/>	<input type="text" value=""/>	

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <input type="text" value="Spline"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text" value=""/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <input type="text" value=""/> Hour
Slope Diff <input type="text" value=""/> SD <input type="text" value=""/>
Sensitivity <input type="text" value=""/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text" value=""/>

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