ACCENT-200 ALPHA 1-GLYCOPROTEIN ACID

Nr kat. 7-221

(PL)

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Cat. No **7-221** (EN)

INTENDED USE

Diagnostic kit for determination of a α 1-glycoprotein acid concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT 400 and ACCENY Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

α 1-glycoprotein acid (orosomucoid) is unique amongst plasma proteins because of its low pH and high carbohydrate content. The protein is a serum transporter for steroid hormone and for many drugs. Its physiological role remains unknown but it is an acute phase reactant. The concentration of this protein in serum is used clinically to monitor acute phase responses and tumour recurrence.

METHOD PRINCIPLE

The α 1-glycoprotein acid presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to α -1-glycoprotein acid concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 35 ml
 2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: Glycylglycin buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human α 1-glycoprotein acid antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 6 hours at room temperature, up to 7 days at 2-8°C or up to 3 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

| | |
|-------|-----------------|
| serum | 0.39 – 1.15 g/l |
|-------|-----------------|

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument is used.

- LOQ:** 0.143 g/l
- Linearity:** up to 2.42 g/l.

Interferences:

Haemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 1669.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

| Repeatability (run to run) n = 30 | Mean [g/l] | SD [g/l] | CV [%] |
|--|---------------|-------------|-----------|
| level 1 | 0.34 | 0.005 | 1.36 |
| level 2 | 0.66 | 0.005 | 0.69 |
| level 3 | 0.92 | 0.004 | 0.41 |
| Reproducibility (day to day) n = 60 | Mean [g/l] | SD [g/l] | CV [%] |
| level 1 | 0.42 | 0.012 | 2.48 |
| level 2 | 0.72 | 0.011 | 1.57 |
| level 3 | 0.96 | 0.012 | 1.22 |

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using at least 30 samples gave following results:

$$y = 1.0751 x - 0.0527 \text{ g/l};$$

$$R = 0.990 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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