ACCENT-200 ANTITHROMBIN III

Nr kat. 7-223 (PL)

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Cat. No **7-223** (EN)

INTENDED USE

Diagnostic kit for determination of antithrombin III concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT MC240 and ACCENT M320.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Antithrombin III (AT III) is an α_2 -glycoprotein of MW 58000 and is made in the liver. AT III is one of the most important regulators of the coagulation system. AT III inactivates thrombin and factors Xa, IXa, XIa and XIIa. This anticoagulant activity is enhanced by the presence of heparin, which forms a ternary complex with AT III and these procoagulant factors. Reduced concentration of AT III in blood means a great and well established risk for thrombotic complications. Clinical low value are associated with congenital or acquired deficiencies caused by a decreased biosynthesis (liver disease, medical treatment) or an increased loss (gastrointestinal disease, nephrotic syndrome) or an increased consumption (sepsis, major trauma due to surgery and burns wounds).

METHOD PRINCIPLE

The antithrombin III, presents in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to antithrombin III concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 35 ml
2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human antithrombin III antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN⁷

Citrated plasma.

Samples may be stored up to 4 hours on ice. If the test cannot be done immediately, the sample should be stored frozen.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES⁷

plasma	0.21 – 0.30 g/l
--------	-----------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT MC240, ACCENT M320 the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzers Hitachi or ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- LOQ:** 0.028 g/l
- Linearity:** up to 0.806 g/l.

Interferences:

Hemoglobin up to 0.06 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 155,5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.16	0.002	0.011
level 2	0.32	0.002	0.005
level 3	0.46	0.005	0.012

Reproducibility (day to day) n = 56	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.18	0.005	0.03
level 2	0.34	0.008	0.02
level 3	0.50	0.010	0.02

Method comparison

A comparison between antithrombin III concentration determined at **ACCENT MC240** (y) and at **BS-800** (x) using 60 plasma samples gave following results:

$$y = 1.0198x - 0.0077 \text{ g/l};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT MC240

Chem	AT III	No.	050	Sample Type	PLASMA
Chemistry	ANTITHROMBIN III	Print name	AT III		
Reaction Type	Endpoint	Reaction Direction	positive		
Pri Wave	340 nm	Sec Wave	700 nm		
Unit	g/L	Decimal	0.001		
Blank Time	-3	-1	Incubation Time	21	
			Reaction Time	32	33
Standard	6	μL	Aspirated		μL
Decreased	6	μL	20	μL	180
Increased		μL		μL	
			Diluent		μL
			Reagent Vol	150	μL
			R2	30	μL
			<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	
Linearity range (Standard)	0,031	0,58	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000	35000
R1 Blank Abs	-35000	35000	On-board Stability		Day(s)
Blank Response	-35000	35000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check					
Q1		Q2		V1	
Q5		Q6		V3	
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		PC1	
				PC2	
				Calibrator Pretreatment	
				Pretreat Sample Vol	μL
				Pretreat Sample Vol	μL
CALIBRATION SETTINGS			AUTO CALIBRATION		
Math model	Spline		<input type="checkbox"/> Bottle Changed		
Factor		Replicates	2	<input type="checkbox"/> Lot Changed	
ACCEPTANCE LIMITS			<input type="checkbox"/> Cal Time		
Cal Time		Hour			
Slope Diff		SD			
Sensitivity		Repeatability	35000		
Deter Coeff					

ACCENT-200 ANTITHROMBIN III

• **ACCENT M320**

Chem <input type="text" value="AT III"/>	No. <input type="text" value="050"/>	Sample Type <input type="text" value="Plasma"/>
Chemistry <input type="text" value="ANTITHROMBIN III"/>	Print name <input type="text" value="AT III"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="g/L"/>	Decimal <input type="text" value="0.001"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="20"/>	Reaction Time <input type="text" value="30"/> <input type="text" value="32"/>
Standard <input type="text" value="8"/> <input type="text" value="8"/> <input type="text" value="8"/> <input type="text" value="8"/>	Aspirated <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="8"/> <input type="text" value="20"/>	Reagent Vol R1 <input type="text" value="200"/>	Reagent Vol R2 <input type="text" value="40"/>
Increased <input type="text" value="8"/> <input type="text" value="20"/>		
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="0,03"/> <input type="text" value="0,56"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Dav(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/>	Q4 <input type="text"/>	V2 <input type="text"/>
Q5 <input type="text"/>	O6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
	Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L			

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Logit-log 5P"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

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