ACCENT-200 CRP ULTRA

Nr kat. **7-225** (PL)

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Cat. No **7-225** (EN)

INTENDED USE

Diagnostic kit for determination of C-reactive protein concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

CRP (C-reactive protein) is an acute phase protein whose concentration is seen to increase as a result of the inflammatory process, most notably in response to pneumococcal (bacterial) infectious, histolytic disease and a variety of disease states. CRP to be used as a marker or general diagnostic indicator of infections and inflammation, in addition to serving as a monitor of patient response to therapy and surgery. Furthermore, regular measurements of CRP in infants can be a useful aid in the early diagnosis of infectious disease.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between CRP in a sample and anti-CRP antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of CRP in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-REAGENT 1 x 25 ml
 2-REAGENT 1 x 25 ml

The reagents, stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 4 weeks.

Concentrations in the test

suspension of latex particles sensitized with anti-CRP antibodies (rabbit) (pH 7.3) 0.20 w/v%
 glycine buffer solution (pH 7.0)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- EUH210 Safety data sheet available on request.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citrate).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. 0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES^{3,4}

serum, plasma	
adults	< 0.5 mg/dl (< 5 mg/l)
children (2 months – 15 years)	0.01 – 0.28 mg/dl (0.1 – 2.8 mg/l)
newborns (0 – 3 weeks)	0.01 – 0.41 mg/dl (0.1 – 4.1 mg/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, the CORMAY CRP ULTRA CALIBRATORS kit (Cat. No 4-276) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 4 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-220S and ACCENT MC240. Results may vary if a different instrument is used.

Sensitivity

0.17 mg/dl – ACCENT-220S
 0.09 mg/dl – ACCENT MC240

Linearity

up to 32 mg/dl – ACCENT-220S
 up to 31.5 mg/dl – ACCENT MC240

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 500 mg/dl and RF up to 560 IU/ml do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-220S n=10	level 1	0.93	0.02	2.35
	level 2	2.39	0.03	1.22
ACCENT MC240 n=20	level 1	0.87	0.02	2.70
	level 2	2.37	0.03	1.40
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-220S n=20	level 1	0.90	0.03	3.15
	level 2	2.30	0.03	1.19
ACCENT MC240 n=80	level 1	0.9	0.03	3.7
	level 2	2.5	0.07	2.9

Method comparison

A comparison between CRP ULTRA values determined at ACCENT-220S (y) and at ADVIA 1650 (x) using 30 samples gave following results:

$$y = 1.1745x - 0.5587 \text{ mg/dl};$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

A comparison between CRP ULTRA values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU 680 (x) using 51 serum samples gave following results:

$$y = 0.9601x + 0.0386 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule	
Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	28
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters			
Test	CRPu	R1	200
No	27	R2	200
Full Name	CRP ULTRA	Sample Volume	4
Standard No	27	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	0.17 32
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	4 30	Factor	
Incuba. Time	20	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule	
Rule	Logit-Log 5P
Sensitivity	1
Replicates	2
Interval (day)	28
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	CRP U	No.	027	Sample Type	SERUM
Chemistry	CRP ULTRA			Print name	CRP U
Reaction Type	Endpoint			Reaction Direction	positive
Pri Wave	670 nm			Sec Wave	
Unit	mg/dl			Decimal	0.01
Blank Time	-3		-1	Incubation Time	20
				Reaction Time	29 30
Standard	3	Sample Vol	Aspirated	Diluent	Reagent Vol
Decreased	3	μL	20	μL	180
Increased		μL		μL	160
		μL		μL	160
		μL		μL	
		Sample Blank	V	Auto Rerun	

Linearity range (Standard)	0.15	31.7	Linearity Limit			
Linearity Range (Decreased)			Substrate Depletion			
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000		
R1 Blank Abs	-40000	40000	On-board Stability			
Blank Response	-40000	40000	Reagent Alarm Limit			
Twin Chemistry			Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check						
O1		O2	V1	O3	O4	V2
Q5		Q6	V3	PC1	PC2	
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment		
		Pretreat Sample Vol	μL	Pretreat Sample Vol	μL	

CALIBRATION SETTINGS	
Math model	Spline
Factor	
Replicates	2

AUTO CALIBRATION	
<input type="checkbox"/> Bottle Changed	
<input type="checkbox"/> Lot Changed	
<input type="checkbox"/> Cal Time	

ACCEPTANCE LIMITS	
Cal Time	Hour
Slope Diff	
Sensitivity	Repeatability
Deter Coeff	