## ACCENT-200 CK-MB

Nr kat. 7-227

(PL)

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Cat. No 7-227 (EN)

### INTENDED USE


Diagnostic kit for determination of CK-MB fraction activity used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased CK-MB serum level is a strong marker of myocardial infarction.

### METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC) with use of antibodies against CK-M fraction. Specific antibodies against CK-M inhibit the complete CK-MM activity (which is the main part of total CK activity) and the CK-M subunit of CK-MB. Only CK-B activity is measured.

The rate of absorbance changes at  $\lambda=340$  nm is directly proportional to half of CK-MB activity (B subunit activity).

### REAGENTS

<b>Package</b>	
1-REAGENT	1 x 30 ml
2-REAGENT	1 x 7 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

### Concentrations in the test

<b>1-REAGENT</b>	
imidazole buffer pH 6.7	100 mmol/l
glucose	20 mmol/l
N-acetylcysteine	20 mmol/l
magnesium acetate	10 mmol/l
EDTA	2 mmol/l
NADP	2 mmol/l
ADP	2 mmol/l
AMP	5 mmol/l
HK	> 2.5 U/ml
polyclonal antibodies against CK-M; inhibiting capacity	8000 U/l

### 2-REAGENT

diadenosinepentaphosphate	10 μmol/l
glucose-6-phosphate-dehydrogenase (G6P-DH)	> 1.5 U/ml
creatine phosphate preservatives	30 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- Do not use reagents past the expiry date.
- Do not interchange caps among reagents.
- Results CK-MB can be falsely high in case of prostate, kidney, ovary, breast and bladder cancer when isoenzyme CK-BB appears in the blood.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT contains imidazole.

### Danger

H360D: May damage the unborn child.  
P201: Obtain special instructions before use.  
P308+P313: IF exposed or concerned: Get medical advice.  
P405: Store locked up.

### SPECIMEN

Serum, free from hemolysis. CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light. Specimens can be stored up to 4-8 hours at 15-25°C or 1-2 days at 2-8°C or 1 month at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionised water is recommended as a reagent blank.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, there is a probability of **cross-contamination** affecting the tests results: CK-MB - CK, CK-MB - CREATININE, GLUCOSE HEX - CK-MB, ASAT - CK-MB, ALAT - CK-MB. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER

### REFERENCE VALUES<sup>9</sup>

serum	37°C	
adults	up to 24 U/l	up to 0.401 μkat/l

The probability that cardiac infarction has occurred is high when CK-MB and total CK activities are above normal values and CK-MB activity is between 6 and 25% of the total CK activity. It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY CK-MB CONTROL N (Cat. No 5-183) and CORMAY CK-MB CONTROL P (Cat. No 5-184) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY CK-MB CALIBRATOR (Cat. No 5-182) is recommended. The calibration curve should be prepared every 6 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 21 U/l (0.35 μkat/l) - ACCENT-200
- LOQ:** 2.6 U/l (0.043 μkat/l) - ACCENT MC240
- Linearity**  
up to 2300 U/l (38.3 μkat/l) - ACCENT-200  
up to 2870 U/l (47.83 μkat/l) - ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin interfere even in small amounts, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=10	level 1	42.9	0.48	1.13
	level 2	97.9	1.04	1.07
ACCENT MC240 n=20	level 1	37.95	0.49	1.29
	level 2	132.42	0.98	0.74
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=80	level 1	22.13	2.47	11.15
	level 2	125.79	5.58	4.44
ACCENT MC240 n=80	level 1	36.4	0.75	2.1
	level 2	133.1	2.03	1.5

### Method comparison

A comparison between CK-MB values determined at ACCENT-200 (y) and at COBAS INTEGRA 400 (x) using 30 samples gave following results:  
y = 0.9436 x + 0.9208 U/l;  
R = 0.960 (R – correlation coefficient)

A comparison between CK-MB values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU 680 (x) using 67 serum samples gave following results:  
y = 1.0495 x – 0.8043 U/l;  
R = 0.999 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 CK-MB

### • ACCENT MC240

Chem	CK-MB	No.	020	Sample Type	SERUM
Chemistry	CK-MB	Print name	CK-MB	Reaction Direction	positive
Reaction Type	Kinetic	Sec Wave	412 nm	Decimal	0.1
Pri Wave	340 nm	Incubation Time	21	Reaction Time	10 20
Unit	U/l	Sample Vol	9	Aspirated	20
Blank Time		Diluent	180	Reagent Vol	30
Standard	9	R1	150	R2	30
Decreased	9	Increased			
Increased		Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>

Linearity range (Standard)	2.6	2870	Linearity Limit	0.2							
Linearity Range (Decreased)			Substrate Depletion	35000							
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/>	Enzyme Linear Extension							
<input type="checkbox"/>	Prozone Check										
O1		O2		V1		O3		O4		V2	
O5		O6		V3		PC1		PC2			
<input type="checkbox"/>	Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment						
		Pretreat Sample Vol		µL	Pretreat Sample Vol		µL				

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Two-point linear	<input type="checkbox"/>	Bottle Changed
Factor		<input type="checkbox"/>	Lot Changed
Replicates	2	<input type="checkbox"/>	Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

### • ACCENT M320

Chem	CK-MB	No.	020	Sample Type	SERUM
Chemistry	CK-MB	Print name	CK-MB	Reaction Direction	positive
Reaction Type	Kinetic	Sec Wave	412 nm	Decimal	0.1
Pri Wave	340 nm	Incubation Time	20	Reaction Time	25 35
Unit	U/l	Sample Vol	9	Aspirated	20
Blank Time		Diluent	180	Reagent Vol	30
Standard	9	R1	150	R2	30
Decreased	9	Increased			
Increased		Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>

Linearity range (Standard)	2.6	2870	Linearity Limit	0.2							
Linearity Range (Decreased)			Substrate Depletion	35000							
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/>	Enzyme Linear Extension							
<input type="checkbox"/>	Prozone Check										
O1		O2		V1		O3		O4		V2	
O5		O6		V3		PC1		PC2			
<input type="checkbox"/>	Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment						
		Pretreat Sample Vol		µL	Pretreat Sample Vol		µL				

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Two-point linear	<input type="checkbox"/>	Bottle Changed
Factor		<input type="checkbox"/>	Lot Changed
Replicates	2	<input type="checkbox"/>	Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

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