## ACCENT-200 MG

Nr kat. 7-229 (PL)

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(EN)

**Danger**

H314 Causes severe skin burns and eye damage.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P310 Immediately call a POISON CENTER or doctor.

### SPECIMEN <sup>6,7</sup>

Serum, heparinized plasma free from hemolysis, 24-hours urine. Recommended anticoagulants: heparine lithium, sodium or ammonium salt.  
Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain approximately 3 times the magnesium concentration found in normal serum.  
Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.  
**Urine preparation:** Acidify urine with some drops of concentrated hydrochloride acid to pH 1.0. Then dilute 1 part of acidified urine with 4 parts of distilled water. Multiply the result by 5. Mix well samples before analysis. 24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT is ready to use.  
Avoid foaming.  
Deionised water is recommended as a reagent blank.  
**Actions required:**  
On analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, BS-120, the reagent should be used in combination with the WASHING PROGRAM using ACCENT WASHING SOLUTION (Cat. No. 3-119) or ACCENT-200 ACID WASHING SOLUTION (Cat. No. 3-109) according to the guidelines contained in the instruction: 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>6</sup>

serum / plasma	mg/dl	mmol/l
newborn 2 – 4 d	1.5 – 2.2	0.62 – 0.91
children 5 mo – 6 y	1.7 – 2.3	0.70 – 0.95
6 – 12 y	1.7 – 2.1	0.70 – 0.86
12 – 20 y	1.7 – 2.2	0.70 – 0.91
adults	1.6 – 2.6	0.66 – 1.07
24-hours urine:	mg/24h	mmol/24h
	72.9 – 145.8	3 – 5

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) - for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) - for determination in urine.  
For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.  
The calibration curve should be prepared every week (ACCENT-200) or every 2 weeks (ACCENT M320 ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- **Sensitivity**  
0.12 mg/dl (0.05 mmol/l) - ACCENT 200  
0.06 mg/dl (0.025 mmol/l) – ACCENT MC240
- **Linearity**  
up to 5.0 mg/dl (2.05 mmol/l) - ACCENT 200  
up to 5.0 mg/dl (2.05 mmol/l) – ACCENT MC240
- **Specificity / Interferences**  
Haemoglobin up to 0.313 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl, triglycerides up to 1000 mg/dl and calcium up to 20 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	1.97	0.02	1.22
	level 2	4.26	0.03	0.63
ACCENT MC240 n=20	level 1	1.97	0.02	1.09
	level 2	4.31	0.05	1.10
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	2.06	0.05	2.31
	level 2	4.21	0.09	2.04
ACCENT MC240 n=80	level 1	1.9	0.05	2.8
	level 2	4.2	0.13	3.2

### Method comparison

A comparison between magnesium values determined at **ACCENT-200** (y) and at **ADVIA 1650** (x) using 25 samples gave following results:  
 $y = 0.9218x + 0.0901$  mg/dl;  
 $R = 0.964$  (R – correlation coefficient)

A comparison between magnesium values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 60 serum samples gave following results:  
 $y = 0.9871x + 0.0268$  mg/dl;  
 $R = 0.967$  (R – correlation coefficient)

A comparison between magnesium values determined at **ACCENT MC240** (y) and at **BS-400** (x) using 20 urine samples gave following results:  
 $y = 0.9847x + 2.1093$  mg/dl;  
 $R = 0.996$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

1. C.K., Mann, J.H., Yoe, Anal. Chem. 28 (2), 202-205 (1956).
2. C. Bohuon, Clin. Chim. Acta 7, 811-817 (1962).
3. C.A. Burtis, E.R. Ashwood, eds. Tietz Fundamentals of Chemical Clinical Chemistry 4th ed. W.B. Saunders, 690-691, (1996).
4. A.J. Pesce, L.A. Kaplan, ed., Methods in Clinical Chemistry, Mosby, 1021-1026 (1987).
5. W.R. Faulkner, S. Meites, ed. Selected Methods for the Small Clinical Laboratory, AAC Press, 277-281 (1982).
6. C.A. Burtis, E.R. Ashwood, ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: Moss D. W., Henderson A. R. 1825 (1999).
7. Regulation of the Minister of Health of 15 December 2017 amending the regulation on quality standards for medical diagnostic and microbiological laboratories.

**Date of issue:** 08. 2023.

### INTENDED USE

Diagnostic kit for determination of magnesium concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Magnesium in human organism occurs mainly in bone (about 50%) but is present also intracellularly in other tissues. Magnesium serves as a cofactor for multiple enzymatic reactions involved in nucleic acids synthesis, transport and production of energy. Magnesium is important in neuromuscular conduction and activation. Reduced magnesium level generates: concentration disturbances, fatigue, muscle tremor, anxiety state.

### METHOD PRINCIPLE

Magnesium forms a purple coloured complex in alkaline solution. In the presence of EGTA, the reaction is specific. The intensity of the purple colour is proportional to the magnesium concentration.

### REAGENTS

**Package**  
1-REAGENT 1 x 31 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. After first opening the reagent's stability on board of the analyser at 2-10°C: 3 weeks (ACCENT-200), 10 weeks (ACCENT M320, ACCENT S120).

### Concentrations in the test

Xylidyl blue ≤ 0.18 mmol/l  
EGTA ≤ 0.12 mmol/l  
CAPS ≤ 0.06 mol/l  
potassium hydroxide ≤ 0.07 mol/l  
buffer (pH 11.5) ≤ 0.14 mol/l  
detergent  
preservative

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- The reagent is air sensitive, to extend reagents stability it is recommended to keep reagent's bottles recapped on the board of analyser
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT contains potassium hydroxide.

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## ACCENT-200 MG

### • ACCENT MC240

Chem	MG	No.	012	Sample Type	SERUM
Chemistry	MAGNESIUM	Print name	MG	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700 nm	Decimals	0.01
Pri Wave	505 nm	Incubation Time	0	Reaction Time	10 11
Unit	mg/dl	Blank Time	-3 -1	Reagent Vol	
Standard	2 μL	Aspirated		Diluent	
Decreased	2 μL	20 μL		180 μL	
Increased					
		<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard)	0.06	5	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	70 Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol		μL		Pretreat Sample Vol					

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
	Replicates 2	<input type="checkbox"/> Cal Time	

<b>ACCEPTANCE LIMITS</b>	
Cal Time	336 Hour
Slope Diff	
Sensitivity	Repeatability 35000
Deter Coeff	

### • ACCENT M320

Chem	MG	No.	012	Sample Type	SERUM
Chemistry	MAGNESIUM	Print name	MG	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700 nm	Decimals	0.01
Pri Wave	505 nm	Incubation Time	0	Reaction Time	16 17
Unit	mg/dl	Blank Time	-3 -1	Reagent Vol	
Standard	2 μL	Aspirated		Diluent	
Decreased	2 μL	20 μL		180 μL	
Increased					
		<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard)	0.13	5	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	70 Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol		μL		Pretreat Sample Vol					

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
	Replicates 2	<input type="checkbox"/> Cal Time	

<b>ACCEPTANCE LIMITS</b>	
Cal Time	336 Hour
Slope Diff	
Sensitivity	Repeatability 35000
Deter Coeff	

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