## ACCENT-200 CREATININE

Nr kat. 7-233 (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of creatinine concentration intended to use in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

### METHOD PRINCIPLE

Modified Jaffe's method, without deproteinization. In alkaline solution picrate reacts with creatinine to form a yellow-red 2,4,6-trinitrocyclohexadienate. The colour intensity is proportional to the creatinine concentration.

### REAGENTS

#### Package

1-REAGENT 4 x 15 ml  
2-REAGENT 2 x 7.5 ml

The reagents, stored at 15-25°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 4 weeks (ACCENT-200) or 6 weeks (ACCENT MC240, ACCENT S120).

### Concentrations in the test

sodium hydroxide ≤ 450 mmol/l  
carbonate buffer ≤ 150 mmol/l  
picric acid ≤ 38.8 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT contains sodium hydroxide.  
2-REAGENT contains picric acid.

### Danger

H314 Causes severe skin burns and eye damage.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

### SPECIMEN

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine collected without preservatives.<sup>9</sup>

Urine preparation: Before analysis urine sample should be diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Mix well probes before measurement.

Specimen can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionised water is recommended as a reagent blank.

#### Actions required:

When performing assays in analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results CREATININE - BIL TOTAL II GEN, CREATININE - BIL DIRECT III GEN, CREATININE - CHOL, CREATININE - MICROALBUMIN, CREATININE - UA / UA PLUS, URINE PROTEINS II GEN - CREATININE, CK - CREATININE, CK-MB - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>7</sup>

serum / plasma	mg/dl	μmol/l
female	0.6 – 1.1	53 – 97
male	0.7 – 1.3	62 – 115
24-hours urine	mg/kg/24h	μmol/kg/24h
female	11 – 20	97 – 177
male	14 – 26	124 – 230

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every week (ACCENT-200) or 6 weeks (ACCENT MC240, ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

0.32 mg/dl (28.3 μmol/l) – ACCENT-200  
0.8 mg/dl (70.7 μmol/l) – ACCENT MC240

#### Linearity

up to 20 mg/dl (1768 μmol/l) – ACCENT-200  
up to 37 mg/dl (3271 μmol/l) – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, triglycerides up to 500 mg/dl, ascorbate up to 62 mg/l and bilirubin up to 20 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	0.88	0.06	6.65
	level 2	3.48	0.08	2.33
ACCENT MC240 n=20	level 1	1.75	0.02	1.28
	level 2	4.58	0.03	0.63
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	1.44	0.07	4.86
	level 2	3.92	0.16	4.00
ACCENT MC240 n=80	level 1	1.5	0.04	2.7
	level 2	4.3	0.10	2.3

#### Method comparison

A comparison between creatinine values determined at ACCENT-200 (y) and at COBAS INTEGRA 400 (x) using 63 serum samples gave following results:

$y = 1.1032x - 0.0409$  mg/dl;  
 $R = 0.999$  (R – correlation coefficient)

A comparison between creatinine values determined at ACCENT MC240 (y) and at BS-400 (x) using 59 serum samples gave following results:

$y = 0.9856x - 0.0035$  mg/dl;  
 $R = 0.999$  (R – correlation coefficient)

A comparison between creatinine values determined at ACCENT MC240 (y) and at BS-400 (x) using 32 urine samples gave following results:

$y = 0.9768x + 1.0125$  mg/dl;  
 $R = 1.000$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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**Calibration Rule**

Rule	Two-Point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • ACCENT-200 II GEN

**Parameters**

Test Name	CREA	R1	200
Test No	8	R2	50
Full Name	Creatinine	Sample Volume	25
Reference No	8	R1 Blank	
Analy. Type	Fixed Time	Mixed Reag. Blank	
Pri. Wave	510 nm	Concentration	0.32   20
Secou. Wave	670 nm	Linearity Limit	0.2
Trend	Ascending	Substrate Limit	
Reac. Time	1   4	Factor	
Incuba. Time	4	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Two-Point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • ACCENT-220S

**Parameters**

Test	CREA	R1	200
No	8	R2	50
Full Name	Creatinine	Sample Volume	25
Standard No	8	R1 Blank	
Reac. Type	Fixed-time	Mixed Rtg. Blank	
Pri. Wave	510 nm	Linearity Range	0.72   30
Sec. Wave	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	2   5	Factor	
Incuba. Time	5	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • BS-120

**Parameters**

Test	CREA	R1	200
No	8	R2	50
Full Name	Creatinine	Sample Volume	25
Standard No	8	R1 Blank	
Reac. Type	Fixed time	Mixed Rtg. Blank	
Pri. Wave	510 nm	Linearity Range	
Sec. Wave	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	1   4	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • ACCENT S120

Chem  No.  Sample Type

Chemistry  Print name

Reaction Type  Reaction Direction

Pri Wave  Sec Wave

Unit  Decimal

Blank Time

Incubation Time  Reaction Time

Standard	<input type="text" value="30"/> μL	Aspirated	<input type="text" value="20"/> μL	Diluent	<input type="text" value="180"/> μL	Reagent Vol	R1 <input type="text" value="200"/> μL
Decreased	<input type="text" value="30"/> μL	<input type="text" value="20"/> μL	<input type="text" value="180"/> μL	R2	<input type="text" value="50"/> μL		
Increased	<input type="text" value=""/> μL	<input type="text" value=""/> μL	<input type="text" value=""/> μL				

Sample Blank  Auto Rerun

Linearity range (Standard)   Linearity Limit

Linearity Range (Decreased)  Substrate Depletion

Linearity Range (Increased)  Mixed Blank Abs

R1 Blank Abs   On-board Stability  Day(s)

Blank Response   Reagent Alarm Limit

Twin Chemistry   Enzyme Linear Extension

Prozone Check

Q1  Q2  V1  Q3  Q4  V2

O5  O6  V3  PC1  PC2

Sample Pretreatment  Control Pretreatment  Calibrator Pretreatment

Pretreat Sample Vol  μL Pretreat Sample Vol  μL

**CALIBRATION SETTINGS**

Math model   Bottle Changed

Factor  Replicates   Lot Changed

Cal Time

**ACCEPTANCE LIMITS**

Cal Time  Hour

Slope Diff  SD

Sensitivity  Repeatability

Deter Coeff

## ACCENT-200 CREATININE

### • ACCENT MC240

Chem <input type="text" value="CREA"/>	No. <input type="text" value="008"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CREATININE"/>	Print name <input type="text" value="CREA"/>	
Reaction Type <input type="text" value="Fixed-time"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="505nm"/>	Sec Wave <input type="text" value="700nm"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	
	Reaction Time <input type="text" value="2"/> <input type="text" value="5"/>	
Standard <input type="text" value="25"/> <input type="text" value="25"/> $\mu$ L	Aspirated <input type="text" value="20"/> $\mu$ L	Diluent <input type="text" value="180"/> $\mu$ L
Decreased <input type="text" value="25"/> $\mu$ L	<input type="text" value="20"/> $\mu$ L	<input type="text" value="180"/> $\mu$ L
Increased <input type="text" value="25"/> $\mu$ L	<input type="text" value="20"/> $\mu$ L	<input type="text" value="180"/> $\mu$ L
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="0.8"/> <input type="text" value="37"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

O1 <input type="text"/>	O2 <input type="text"/>	V1 <input type="text"/>	O3 <input type="text"/>	O4 <input type="text"/>	V2 <input type="text"/>
O5 <input type="text"/>	O6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	<input type="text"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
<input type="text"/>	Pretreat Sample Vol <input type="text"/> $\mu$ L	Pretreat Sample Vol <input type="text"/> $\mu$ L			

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

### • ACCENT M320

Chem <input type="text" value="CREA"/>	No. <input type="text" value="008"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CREATININE"/>	Print name <input type="text" value="CREA"/>	
Reaction Type <input type="text" value="Fixed-time"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="505nm"/>	Sec Wave <input type="text" value="700nm"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="7"/>	
	Reaction Time <input type="text" value="1"/> <input type="text" value="5"/>	
Standard <input type="text" value="25"/> $\mu$ L	Aspirated <input type="text" value="20"/> $\mu$ L	Diluent <input type="text" value="180"/> $\mu$ L
Decreased <input type="text" value="25"/> $\mu$ L	<input type="text" value="20"/> $\mu$ L	<input type="text" value="180"/> $\mu$ L
Increased <input type="text" value="25"/> $\mu$ L	<input type="text" value="20"/> $\mu$ L	<input type="text" value="180"/> $\mu$ L
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="0.8"/> <input type="text" value="37"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/>	Q4 <input type="text"/>	V2 <input type="text"/>
Q5 <input type="text"/>	Q6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	<input type="text"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
<input type="text"/>	Pretreat Sample Vol <input type="text"/> $\mu$ L	Pretreat Sample Vol <input type="text"/> $\mu$ L			

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

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