ACCENT-200 PHOSPHORUS

Nr kat. 7-243 (PL)

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Cat. No 7-243 (EN)

INTENDED USE

PHOSPHORUS reagent is intended to determine quantitatively phosphorus levels in plasma, serum and urine. It is intended to monitor and as an aid to the diagnosis of clinical conditions associated with abnormal phosphorus levels. PHOSPHORUS reagent is intended to use on automatic analysers ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180, ACCENT S120 / BS-230, ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E, BS-120, ACCENT 400 and ACCENT Neo200. It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY^{1,2,3,4}

Determinations of phosphorus levels in plasma, serum and urine are used to monitor and as an aid to the diagnosis of parathyroid and kidney diseases or vitamin D imbalance. Increased phosphorus levels may be associated with increased phosphate intake, cell lysis (leukemia), transcellular phosphate shift (acidosis). Decreased phosphorus levels may occur in refeeding syndrome, decreased net intestinal phosphate absorption, lowered phosphate threshold in the kidneys.

METHOD PRINCIPLE^{5,6}

Direct phosphomolybdate reaction without deproteinization. Phosphate ions form with molybdate ions in acid solution proportional amounts of unreduced phosphomolybdate complex. The concentration of the complex formed is determined by measuring its absorbance.

REAGENTS

Package
1-REAGENT 1 x 32 ml

Number of tests:
ACCENT-200 (II GEN) 90
ACCENT-220S 90
ACCENT S120 140
ACCENT MC240 140
ACCENT M320 140
BS-120 90

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT
ammonium molybdate 0.4 mmol/l
sulphuric acid
hydrochloric acid
detergents

REAGENT STABILITY

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyser at 2-10°C is stable for 12 weeks (ACCENT-200 (II GEN), ACCENT S120, ACCENT MC240).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains sulfuric acid (VI) and hydrochloric acid.

Danger

H314 Causes severe skin burns and eye damage.

H318 Causes serious eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

SPECIMEN^{7,8,9}

Serum, heparinized plasma free from hemolysis, 24-hours urine.

Serum is the preferred specimen. Level of inorganic phosphate in heparinized plasma is about 0.2 - 0.3 mg/dl (0.06 - 0.10 mmol/l) lower than in serum.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain several times higher phosphate concentration than normal serum.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Urine preparation: For collection of 24-hour urine specimens, 20-30 ml of 6M HCl should be added into the container to avoid phosphate precipitations. Dilute 1 part of acidified urine with 9 parts of distilled water before determination. Multiply the result by the dilution factor.

24-hours urine samples can be stored up to 2 days at 20-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

1-REAGENT is ready to use.

Deionized water is recommended as a reagent blank.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO - PHOSPHORUS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionized water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200 (II GEN), ACCENT S120, ACCENT MC240).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES⁷

serum / plasma	mg/dl	mmol/l
age:		
0 - 9 d	4.5 - 9.0	1.45 - 2.91
10 d - 2 y	4.0 - 6.5	1.29 - 2.10
3 - 9 y	3.2 - 5.8	1.03 - 1.87
10 - 15 y	3.3 - 5.4	1.07 - 1.74
16 - 59 y	2.4 - 4.4	0.78 - 1.42
60 - 89 y males	2.3 - 3.7	0.74 - 1.20
60 - 89 y females	2.8 - 4.0	0.90 - 1.26
>89 y females	2.5 - 4.2	0.81 - 1.36
>89 y males	2.2 - 3.9	0.71 - 1.26
24-hours urine	g/24h	mmol/24h
	0.4 - 1.3	12.9 - 42.0

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

Phosphorus concentration in 24-hours urine - calculation

phosphorus concentration in 24-hours urine [g/24h]	=	phosphorus concentration in sample of 24-hours urine [mg/dl]	x	urine volume of 24-hours urine [dl/24h]	÷	1000
--	---	--	---	---	---	------

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 (II GEN), ACCENT MC240 and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.03 mg/dl (0.01 mmol/l) - BS-400

LoD (Limit of Detection):

0.04 mg/dl (0.013 mmol/l) - BS-400

LoQ (Limit of Quantitation):

0.25 mg/dl (0.081 mmol/l) - ACCENT-200 (II GEN)
0.11 mg/dl (0.036 mmol/l) - ACCENT MC240

Linearity

up to 20.0 mg/dl (6.46 mmol/l) - ACCENT-200 (II GEN)
up to 19.0 mg/dl (6,137 mmol/l) - ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Measurement range

0.25 mg/dl (0.081 mmol/l) - 20.0 mg/dl (6.46 mmol/l) ACCENT-200 (II GEN)
0.11 mg/dl (0.036 mmol/l) - 19.0 mg/dl (6,137 mmol/l) ACCENT MC240

Specificity / Interferences

Haemoglobin up to 0,12 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

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LIST OF CHANGES

Previous version: 06	Current version: 07
Sections updated: <i>INTENDED USE</i>	

Date of issue: 05. 2022.

Method comparison

A comparison between phosphorus values determined at ACCENT-200 (II GEN) (y) and at BECKMAN COULTER AU680 (x) using 80 serum samples gave following results:

$$y = 1.0137 x + 0.0364 \text{ mg/dl};$$
$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at ACCENT-200 (II GEN) (y) and at BECKMAN COULTER AU680 (x) using 30 plasma samples gave following results:

$$y = 0.9878 x + 0.1429 \text{ mg/dl};$$
$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at ACCENT-200 (II GEN) (y) and at BECKMAN COULTER AU680 (x) using 28 urine samples gave following results:

$$y = 0.9538 x + 0.9787 \text{ mg/dl};$$
$$R = 0.993 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 71 serum samples gave following results:

$$y = 0.9858 x + 0.0843 \text{ mg/dl};$$
$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 30 plasma samples gave following results:

$$y = 0.9504 x + 0.2636 \text{ mg/dl};$$
$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 48 urine samples gave following results:

$$y = 0.9191 x + 1.5334 \text{ mg/dl};$$
$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents:
 - 1-REAAGENT – 18 01 06*
- Empty packages:
 - 1-REAAGENT – 15 01 10*
- Wastewater from the analyzer:
 - 1-REAAGENT – 18 01 03*

INCIDENTS¹¹

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

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- Dalay J.A., Ertinghausen G. Direct method for determining inorganic phosphate in serum with the „CentrifChem“ Clin. Chem. 18, 263-265 (1972).
- M.A Muñoz, M Balón, C Fernandez, Direct determination of inorganic phosphorus in serum with a single reagent., Clin. Chem. 29, 2, 372-374 (1983).
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 852-855 (2006).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 430 (2018).
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 883 (1989).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters		R1		R2	
Test	PHOS		300		
No	13				
Full Name	Phosphorus	Sample Volume	3		
Standard No	13	R1 Blank			
Reac. Type	Endpoint	Mixed Rtg. Blank			
Pri. Wave.	340 nm	Linearity Range	0.31	18	
Sec. Wave.	670 nm	Linearity Limit			
Direction	Increase	Substrate Limit			
Reac. Time	0 22	Factor			
Incuba. Time		<input type="checkbox"/> Prozone check			
Unit	mg/dl	q1	<input type="checkbox"/>	q2	<input type="checkbox"/>
Precision	0.01	q3	<input type="checkbox"/>	q4	<input type="checkbox"/>
		PC	<input type="checkbox"/>	Abs	<input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters		R1		R2	
Test	PHOS		300		
No	13				
Full Name	Phosphorus	Sample Volume	3		
Standard No	13	R1 Blank			
Reac. Type	Endpoint	Mixed Rtg. Blank			
Pri. Wave.	340 nm	Linearity Range	0.31	18	
Sec. Wave.	670 nm	Linearity Limit			
Direction	Increase	Substrate Limit			
Reac. Time	1 22	Factor			
Incuba. Time	16	<input type="checkbox"/> Prozone check			
Unit	mg/dl	q1	<input type="checkbox"/>	q2	<input type="checkbox"/>
Precision	0.01	q3	<input type="checkbox"/>	q4	<input type="checkbox"/>
		PC	<input type="checkbox"/>	Abs	<input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	<input type="text" value="PHOS"/>	No.	<input type="text" value="013"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="PHOSPHORUS"/>	Print name	<input type="text" value="PHOS"/>	Reaction Type	<input type="text" value="Endpoint"/>
Reaction Type	<input type="text" value="Endpoint"/>	Reaction Direction	<input type="text" value="positive"/>	Pri Wave	<input type="text" value="340 nm"/>
Pri Wave	<input type="text" value="340 nm"/>	Sec Wave	<input type="text" value="670 nm"/>	Unit	<input type="text" value="mg/dL"/>
Unit	<input type="text" value="mg/dL"/>	Decimal	<input type="text" value="0.01"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>
Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time	<input type="text" value="0"/>	Reaction Time	<input type="text" value="28"/> <input type="text" value="30"/>
Standard	Sample Vol <input type="text" value="2"/> uL	Aspirated	<input type="text" value="20"/> uL	Diluent	<input type="text" value="180"/> uL
Decreased	<input type="text" value="2"/> uL	<input type="text" value="20"/> uL	<input type="text" value="180"/> uL	Reagent Vol	R1 <input type="text" value="200"/> uL
Increased	<input type="text" value=""/> uL	<input type="text" value=""/> uL	<input type="text" value=""/> uL	R2	<input type="text" value=""/> uL
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			

Linearity range (Standard)	<input type="text" value="0.23"/> <input type="text" value="18"/>	Linearity Limit	<input type="text" value=""/>
Linearity Range (Decreased)	<input type="text" value=""/>	Substrate Depletion	<input type="text" value=""/>
Linearity Range (Increased)	<input type="text" value=""/>	Mixed Blank Abs	<input type="text" value="-40000"/> <input type="text" value="40000"/>
R1 Blank Abs	<input type="text" value="-40000"/> <input type="text" value="40000"/>	On-board Stability	<input type="text" value=""/> Dav(s)
Blank Response	<input type="text" value="-40000"/> <input type="text" value="40000"/>	Reagent Alarm Limit	<input type="text" value=""/>
Twin Chemistry	<input type="text" value=""/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check			
O1 <input type="text" value=""/>	O2 <input type="text" value=""/>	V1 <input type="text" value=""/>	O3 <input type="text" value=""/>
O5 <input type="text" value=""/>	O6 <input type="text" value=""/>	V3 <input type="text" value=""/>	PC1 <input type="text" value=""/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	
<input type="text" value=""/>	Pretreat Sample Vol <input type="text" value=""/> uL	Pretreat Sample Vol <input type="text" value=""/> uL	

CALIBRATION SETTINGS		AUTO CALIBRATION
Math model	<input type="text" value="Multi-point linear"/>	<input type="text" value=""/> Bottle Changed
Factor	<input type="text" value=""/> Replicates <input type="text" value="2"/>	<input type="text" value=""/> Lot Changed
		<input type="text" value=""/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value=""/> Hour
Slope Diff <input type="text" value=""/> SD <input type="text" value=""/>
Sensitivity <input type="text" value=""/> Repeatability <input type="text" value="40000"/>
Deter Coeff <input type="text" value=""/>

• ACCENT MC240

Chem	<input type="text" value="PHOS"/>	No.	<input type="text" value="013"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="PHOSPHORUS"/>	Print name	<input type="text" value="PHOS"/>	Reaction Type	<input type="text" value="Endpoint"/>
Reaction Type	<input type="text" value="Endpoint"/>	Reaction Direction	<input type="text" value="positive"/>	Pri Wave	<input type="text" value="340 nm"/>
Pri Wave	<input type="text" value="340 nm"/>	Sec Wave	<input type="text" value="700 nm"/>	Unit	<input type="text" value="mg/dL"/>
Unit	<input type="text" value="mg/dL"/>	Decimal	<input type="text" value="0.01"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>
Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time	<input type="text" value="0"/>	Reaction Time	<input type="text" value="28"/> <input type="text" value="30"/>
Standard	Sample Vol <input type="text" value="2"/> uL	Aspirated	<input type="text" value="20"/> uL	Diluent	<input type="text" value="180"/> uL
Decreased	<input type="text" value="2"/> uL	<input type="text" value="20"/> uL	<input type="text" value="180"/> uL	Reagent Vol	R1 <input type="text" value="200"/> uL
Increased	<input type="text" value=""/> uL	<input type="text" value=""/> uL	<input type="text" value=""/> uL	R2	<input type="text" value=""/> uL
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			

Linearity range (Standard)	<input type="text" value="0.11"/> <input type="text" value="19"/>	Linearity Limit	<input type="text" value=""/>
Linearity Range (Decreased)	<input type="text" value=""/>	Substrate Depletion	<input type="text" value=""/>
Linearity Range (Increased)	<input type="text" value=""/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text" value=""/> Dav(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text" value=""/>
Twin Chemistry	<input type="text" value=""/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check			
O1 <input type="text" value=""/>	O2 <input type="text" value=""/>	V1 <input type="text" value=""/>	O3 <input type="text" value=""/>
O5 <input type="text" value=""/>	O6 <input type="text" value=""/>	V3 <input type="text" value=""/>	PC1 <input type="text" value=""/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	
<input type="text" value=""/>	Pretreat Sample Vol <input type="text" value=""/> uL	Pretreat Sample Vol <input type="text" value=""/> uL	

CALIBRATION SETTINGS		AUTO CALIBRATION
Math model	<input type="text" value="Multi-point linear"/>	<input type="text" value=""/> Bottle Changed
Factor	<input type="text" value=""/> Replicates <input type="text" value="2"/>	<input type="text" value=""/> Lot Changed
		<input type="text" value=""/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text" value=""/> SD <input type="text" value=""/>
Sensitivity <input type="text" value=""/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text" value=""/>

• ACCENT M320

Chem	PHOS		No.	013		Sample Type	SERUM	
Chemistry	PHOSPHORUS		Print name	PHOS		Reaction Direction	positive	
Reaction Type	Endpoint		Pri Wave	340 nm		Sec Wave	700 nm	
Unit	mg/dL		Decimal	0.01		Incubation Time	0	
Blank Time	-3		Reaction Time	42		44		
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol	R1	R2		
Decreased	2	uL	20	uL	200	uL		
Increased	uL	uL	180	uL	uL	uL		
	uL	uL	uL	uL				
	<input type="checkbox"/>	Sample Blank	<input checked="" type="checkbox"/>	Auto Rerun				
Linearity range (Standard)	0.23	18	Linearity Limit					
Linearity Range (Decreased)			Substrate Depletion					
Linearity Range (Increased)			Mixed Blank Abs	-35000	35000			
R1 Blank Abs	-35000	35000	On-board Stability	Dav(s)				
Blank Response	-35000	35000	Reagent Alarm Limit					
Twin Chemistry			<input type="checkbox"/>	Enzyme Linear Extension				
	<input type="checkbox"/>	Prozone Check						
O1	O2	V1	O3	O4	V2			
O5	O6	V3	PC1	PC2				
<input type="checkbox"/>	Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment			
		Pretreat Sample Vol	uL	Pretreat Sample Vol	uL			
CALIBRATION SETTINGS			AUTO CALIBRATION					
Math model	Multi-point linear		<input type="checkbox"/>	Bottle Changed				
Factor		Replicates	2	<input type="checkbox"/>	Lot Changed			
				<input type="checkbox"/>	Cal Time			
ACCEPTANCE LIMITS								
Cal Time		Hour						
Slope Diff		SD						
Sensitivity		Repeatability	35000					
Deter Coeff								

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