## ACCENT-200 PHOSPHORUS

Nr kat. 7-243

(PL)

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Cat. No 7-243 (EN)

### INTENDED USE

Diagnostic kit for determination of inorganic phosphorus concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 and BS-120 / BS-130. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Phosphorus is present in all body cells as a component of nucleic acids, phospholipids and phosphoproteins. Phosphorus is essential for intracellular storage and conversion of energy (ATP, creatine phosphate) and participates in carbohydrates metabolism. In the blood phosphorus is present as a mixture of inorganic phosphates  $\text{HPO}_4^{2-}$  and  $\text{H}_2\text{PO}_4^-$ . Besides phosphorus and calcium constitute mineral portion of bone. Continuous flux of phosphorus in organism is controlled by parathyroid hormone (PTH), vitamin D and calcitonin. Phosphorus serum level abnormalities are caused usually by disorders of vitamin D metabolism or parathyroid and kidney diseases.

### METHOD PRINCIPLE

Direct phosphomolybdate reaction without deproteinization. Phosphate ions form with molybdate ions in acid solution proportional amounts of unreduced phosphomolybdate complex. The concentration of the complex formed is determined by measuring its absorbance.

### REAGENTS

**Package**  
1-Reagent 1 x 32 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyser at 2-10°C is stable for 12 weeks (ACCENT-200, ACCENT MC240).

### Concentrations in the test

ammonium molybdate  $\leq 0.48$  mmol/l  
sulphuric acid  $\leq 180$  mmol/l  
hydrochloric acid  $\leq 120$  mmol/l  
detergents

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. Disposable plastic ware is recommended for the test.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1- REAGENT contains sulfuric acid (VI) and hydrochloric acid.

### Danger

H314 Causes severe skin burns and eye damage.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.

P305 +P351 +P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 - Immediately call a POISON CENTER or doctor.

### SPECIMEN

Serum, heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis, 24-hours urine. Serum is the preferred specimen. Level of inorganic phosphate in heparinized plasma is about 0.2 to 0.3 mg/dl (0.06-0.10 mmol/l) lower than in serum.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain several times higher phosphate concentration than normal serum.

**Urine preparation:** To prevent phosphate precipitation in urine, specimens should be collected in HCl, 20-30 ml of 6 mol/L for 24-h specimen. Then dilute 1 part of acidified urine with 9 parts of distilled water. Multiply the result by the dilution factor.

Serum and plasma can be stored up to 7 days at 2 – 8°C. For longer storage samples should be frozen at -20°C. 24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.  
Deionized water is recommended as a reagent blank.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 / BS-130 there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO – PHOSPHORUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>7</sup>

serum / plasma	mg/dl	mmol/l
age:		
0 – 10 d	4.5 – 9.0	1.45 – 2.91
10 d – 24 mo	4.5 – 6.7	1.45 – 2.16
24 mo – 12 y	4.5 – 5.5	1.45 – 1.78
12 – 60 y	2.7 – 4.5	0.87 – 1.45
> 60 y males	2.3 – 3.7	0.74 – 1.20
> 60 y females	2.8 – 4.1	0.90 – 1.32
24-hours urine	g/24h	mmol/24h
	0.4 – 1.3	12.9 – 42.0

It is recommended for each laboratory to establish its own reference ranges for local population.

### Phosphorus concentration in 24-hours urine – calculation

phosphorus concentration in 24-hours urine [g/24h]	phosphorus concentration in sample of 24-hours urine [mg/dl]	urine volume of 24-hours urine [dl/24h]	=	x	÷	1000
--	--	---	---	---	---	------

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP

(Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT MC240) with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### LoB (Limit of Blank):

0.0 mg/dl (0.0 mmol/l)- ACCENT-200  
0.02 mg/dl (0.006 mmol/l)- ACCENT MC240

#### LoD (Limit of Detection):

0.03 mg/dl (0.0097 mmol/l)- ACCENT-200  
0.06 mg/dl (0.019 mmol/l)- ACCENT MC240

#### LoQ (Limit of Quantitation):

0.31 mg/dl (0.1 mmol/l)- ACCENT-200  
0.23 mg/dl (0.074 mmol/l)- ACCENT MC240

#### Linearity

up to 18.0 mg/dl (5.8 mmol/l) - ACCENT-200, ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.12 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	3.37	0.06	1.8
	level 2	7.63	0.08	1.1
ACCENT MC240 n=20	level 1	3.33	0.04	1.23
	level 2	7.16	0.07	1.03
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	3.38	0.09	2.7
	level 2	7.29	0.08	1.0
ACCENT MC240 n=80	level 1	3.32	0.11	3.4
	level 2	7.20	0.07	1.0

#### Method comparison

A comparison between phosphorus values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 72 serum samples gave following results:

$y = 1.0288x - 0.0517$  mg/dl;  
 $R = 0.997$  (R – correlation coefficient)

A comparison between phosphorus values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 30 urine samples gave following results:

$y = 0.92x + 3.0939$  mg/dl;  
 $R = 0.998$  (R – correlation coefficient)

A comparison between phosphorus values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 71 serum samples gave following results:

$y = 0.9858x + 0.0843$  mg/dl;  
 $R = 0.997$  (R – correlation coefficient)

A comparison between phosphorus values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 50 urine samples gave following results:

$y = 0.9016x + 2.4503$  mg/dl;  
 $R = 0.999$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 PHOSPHORUS

### • ACCENT MC240

Chem	PHOS	No.	013	Sample Type	SERUM
Chemistry	PHOSPHORUS	Print name	PHOS	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700 nm	Incubation Time	0
Pri Wave	340 nm	Decimal	0.01	Reaction Time	28 30
Unit	mg/dL	Sample Vol	-3 -1	Aspirated	-1
Standard	2 μL	Diluent	μL	Reagent Vol	200 μL
Decreased	2 μL	20 μL	180 μL	R1	200 μL
Increased	μL	μL	μL	R2	μL
Sample Blank		<input checked="" type="checkbox"/>		Auto Rerun	

Linearity range (Standard)	0.23	18	Linearity Limit			
Linearity Range (Decreased)			Substrate Depletion			
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000		
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)		
Blank Response	-35000	35000	Reagent Alarm Limit			
Twin Chemistry			Enzyme Linear Extension			
Prozone Check						
Q1	Q2	V1	Q3	Q4	V2	
Q5	Q6	V3	PC1	PC2		
Sample Pretreatment	Control Pretreatment	Calibrator Pretreatment	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear		Bottle Changed
Factor	Replicates 2		Lot Changed
			Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

### • ACCENT M320

Chem	PHOS	No.	013	Sample Type	SERUM
Chemistry	PHOSPHORUS	Print name	PHOS	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700 nm	Incubation Time	0
Pri Wave	340 nm	Decimal	0.01	Reaction Time	42 44
Unit	mg/dL	Sample Vol	-3 -1	Aspirated	-1
Standard	2 μL	Diluent	μL	Reagent Vol	200 μL
Decreased	2 μL	20 μL	180 μL	R1	200 μL
Increased	μL	μL	μL	R2	μL
Sample Blank		<input checked="" type="checkbox"/>		Auto Rerun	

Linearity range (Standard)	0.23	18	Linearity Limit			
Linearity Range (Decreased)			Substrate Depletion			
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000		
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)		
Blank Response	-35000	35000	Reagent Alarm Limit			
Twin Chemistry			Enzyme Linear Extension			
Prozone Check						
Q1	Q2	V1	Q3	Q4	V2	
Q5	Q6	V3	PC1	PC2		
Sample Pretreatment	Control Pretreatment	Calibrator Pretreatment	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL

<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear		Bottle Changed
Factor	Replicates 2		Lot Changed
			Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

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