ACCENT-200 BIL TOTAL

Nr kat. **7-245** (PL)

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LITERATURA

1. Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2); 116-122.
2. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
3. Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

Data wydania: 10. 2023.

ACCENT-200 BIL TOTAL

Cat. No **7-245** (EN)

INTENDED USE

Diagnostic kit for determination of total bilirubin concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Bilirubin is a yellow pigment – product of heme degradation. For clinical purposes, bilirubin is expressed as two fractions: conjugated and unconjugated. In hepatocytes bilirubin is enzymatically conjugated with glucuronic acid residues. This form is called direct or conjugated. Bilirubin without glucuronic acid modification is bound to albumin and is termed unconjugated or indirect. Indirect bilirubin is calculated as the difference between total and direct bilirubin. Serum bilirubin measurement is widely used as a screening test for liver functions. Hiperbilirubinemia is usually the result of jaundice (mechanical, hemolytic), Dubin-Jonson syndrome, Gilbert's syndrome, Crigler-Najjar syndrome, bile ducts disease.

METHOD PRINCIPLE

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in an acidic solution, total bilirubin (both conjugated – direct, and unconjugated bilirubin) is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the total bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

REAGENTS

Package

1-REAGENT 3 x 29 ml
2-REAGENT 2 x 11.5 ml

The reagents, stored at 10-25°C are stable up to the kit expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

Concentrations in the test

1-REAGENT
citrate buffer (pH 2.8) 90 mmol/l
detergent
2-REAGENT
phosphate buffer (pH 7.0) 4.6 mmol/l
sodium metavanadate 3.0 mmol/l

Warnings and notes

- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
H410 Very toxic to aquatic life with long lasting effects.

P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P273 Avoid release to the environment.

P391 Collect spillage.

SPECIMEN

Serum free from hemolysis.

Serum should be separated from red blood cells as soon as possible after blood collection. Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended. It is recommended to follow CLSI procedures regarding specimen collecting and handling.

Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark at 2-8°C, at the most 3 days.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: BIL TOTAL II GENERACJA – MICROALBUMIN, CREATININE - BIL TOTAL II GENERACJA. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES⁴

Age	mg/dl	µmol/l
0 - 1 days	< 8	< 137
1 - 2 days	< 12	< 205
3 - 5 days	< 16	< 274
5 days - 60 years	0.3 - 1.2	5 - 21
60 - 90 years	0.2 - 1.1	3 - 19
> 90 years	0.2 - 0.9	3 - 15

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 and/or ACCENT-200 II GEN and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity**
0.13 mg/dl (2.2 µmol/l) – ACCENT-200
0.14 mg/dl (2.395 µmol/l) – ACCENT MC240
- Linearity**
do 49 mg/dl (838 µmol/l) – ACCENT-200
do 50 mg/dl (855,2 µmol/l) – ACCENT MC240
- Specificity / Interferences**
Haemoglobin up to 0.25 g/dl, ascorbic acid up to 500 mg/l and intralipid up to 250 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 II GEN n=10	level 1	0.93	0.01	0.72
	level 2	4.08	0.04	1.08
ACCENT MC240 n=20	level 1	1.14	0.01	0.66
	level 2	4.48	0.02	0.50
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 II GEN n=20	level 1	1.00	0.05	4.49
	level 2	4.27	0.19	4.35
ACCENT MC240 n=80	level 1	1.1	0.02	1.8
	level 2	4.5	0.07	1.5

Method comparison

A comparison between total bilirubin values determined at ACCENT-200 (y) and at OLYMPUS AU 400 (x) using 53 samples gave following results:

y = 1.1355 x - 0.1294 mg/dl;

R = 0.998 (R – correlation coefficient)

A comparison between total bilirubin values determined at ACCENT MC240 (y) and at BS-800 (x) using 58 serum samples gave following results:

y = 0.9905 x + 0.0112 mg/dl;

R = 1.000 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2); 116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters			
Test Name	BIL T	R1	280
Test No	6	R2	70
Full Name	Bil Total	Sample Volume	20
Reference No	6	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	450 nm	Concentration	0.13 51
Secon. Wave.	546 nm	Linearity Limit	
Trend	Descending	Substrate Limit	
Reac. Time	-1 20	Factor	
Incuba. Time	10	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters			
Test Name	BIL T	R1	280
Test No	6	R2	70
Full Name	Bil Total	Sample Volume	20
Standard No	6	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	450 nm	Linearity Range	0.12 50
Sec. Wave.	546 nm	Linearity Limit	
Direction	Decrease	Substrate Limit	
Reac. Time	-1 22	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters			
Test Name	BIL T	R1	280
Test No	6	R2	70
Full Name	Bil Total	Sample Volume	10
Standard No	6	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	450 nm	Linearity Range	0.13 47
Sec. Wave.	546 nm	Linearity Limit	
Direction	Descending	Substrate Limit	
Reac. Time	0 18	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	BIL TOTAL	No.	006	Sample Type	SERUM
Chemistry	BILIRUBIN TOTAL	Print name	BIL T	Reaction Direction	negative
Reaction Type	Endpoint	Sec Wave	546 nm	Decimal	0.01
Pri Wave	450 nm	Incubation Time	12	Reaction Time	28 30
Unit	mg/dL	Blank Time	-3 -1	Standard	Sample Vol 8 μL Aspirated μL Diluent μL
Decreased	8 μL	Standard	8 μL	Decreased	8 μL
Increased	μL	Standard	20 μL	Increased	μL
		Standard	180 μL		
		Standard	μL		
		Sample Blank	<input type="checkbox"/>	Auto Rerun	<input type="checkbox"/>
Linearity range (Standard)	0.14	50	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000	
R1 Blank Abs	-40000	40000	On-board Stability		
Blank Response	-40000	40000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension	<input type="checkbox"/>	
<input type="checkbox"/> Prozone Check			Q1 <input type="checkbox"/> Q2 <input type="checkbox"/> V1 <input type="checkbox"/> Q3 <input type="checkbox"/> Q4 <input type="checkbox"/> V2 <input type="checkbox"/>		
Q5 <input type="checkbox"/> O6 <input type="checkbox"/> V3 <input type="checkbox"/> PC1 <input type="checkbox"/> PC2 <input type="checkbox"/>			<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	
<input type="checkbox"/> Calibrator Pretreatment			Pretreat Sample Vol <input type="checkbox"/> μL	Pretreat Sample Vol <input type="checkbox"/> μL	
CALIBRATION SETTINGS	Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed		
Factor <input type="checkbox"/>	Replicates	2	<input type="checkbox"/> Lot Changed		
ACCEPTANCE LIMITS	Cal Time	<input type="checkbox"/> Hour	<input type="checkbox"/> Cal Time		
Slope Diff <input type="checkbox"/>	SD	<input type="checkbox"/>			
Sensitivity <input type="checkbox"/>	Repeatability	40000			
Deter Coeff <input type="checkbox"/>					

ACCENT-200 BIL TOTAL

ACCENT MC240

Chem	BIL T	No.	006	Sample Type	SERUM		
Chemistry	BILIRUBIN TOTAL	Print name	BIL T	Reaction Direction	negative		
Reaction Type	Endpoint	Pri Wave	450 nm	Sec Wave	546 nm		
Unit	mg/dL	Decimal	0.01	Incubation Time	21		
Blank Time	-3	-1	Reaction Time	28	30		
Standard	9	Aspirated	20	Diluent	180	Reagent Vol	200
Decreased	9	20	180	R1	200	R2	50
Increased							
		Sample Blank	V	Auto Rerun			

Linearity range (Standard)	0,14	50,00	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability		Day(s)						
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry		Enzyme Linear Extension									
		Prozone Check									
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
		Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment					
		Pretreat Sample Vol		μL	Pretreat Sample Vol		μL				

CALIBRATION SETTINGS	AUTO CALIBRATION				
Math model	Multi-point linear		Bottle Changed		
Factor		Replicates	2		Lot Changed
					Cal Time

ACCEPTANCE LIMITS			
Cal Time		Hour	
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

ACCENT M320

Chem	BIL T	No.	006	Sample Type	SERUM		
Chemistry	BILIRUBIN TOTAL	Print name	BIL T	Reaction Direction	negative		
Reaction Type	Endpoint	Pri Wave	450 nm	Sec Wave	546 nm		
Unit	mg/dL	Decimal	0.01	Incubation Time	15		
Blank Time	-3	-1	Reaction Time	28	30		
Standard	8	Aspirated	20	Diluent	180	Reagent Vol	200
Decreased	8	20	180	R1	200	R2	50
Increased							
		Sample Blank	V	Auto Rerun			

Linearity range (Standard)	0,12	50	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability		Day(s)						
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry		Enzyme Linear Extension									
		Prozone Check									
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
		Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment					
		Pretreat Sample Vol		μL	Pretreat Sample Vol		μL				

CALIBRATION SETTINGS	AUTO CALIBRATION				
Math model	Multi-point linear		Bottle Changed		
Factor		Replicates	2		Lot Changed
					Cal Time

ACCEPTANCE LIMITS			
Cal Time		Hour	
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

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