ACCENT-200 CALCIUM ARSENAZO

Nr kat. 7-247

(PL)

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Cat. No 7-247 (EN)

INTENDED USE

Diagnostic kit for determination of calcium concentration, used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, BS-120, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

METHOD PRINCIPLE

Photometric test using arsenazo III.

Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

REAGENTS

PACKAGE

1-Reagent 2 x 36 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyzer at 2-10°C is stable for 12 weeks (ACCENT-200, ACCENT MC240, ACCENT S120).

Concentrations in the test

phosphate buffer (pH 7.5) ≤ 80 mmol/l
8-hydroxyquinoline-5-sulfonic acid ≤ 6 mmol/l
arsenazo III ≤ 120 μmol/l
detergents, preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended.

- Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

SPECIMEN

Serum. Random or 24-hours urine.

Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months.

Discard contaminated specimens.

24-hours urine preparation: To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

Deionized water is recommended as a reagent blank.

For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO – PHOSPHORUS, RF - CALCIUM ARSENAZO, CALCIUM ARSENAZO - LIPASE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁴

serum	mg/dl	mmol/l
adult	8.6 – 10.3	2.15 – 2.57
random urine	mg/dl	mmol/l
male	0.9 – 37.9	0.225 – 9.47
female	0.5 – 35.7	0.125 – 8.92
24-hours urine	mg/24h	mmol/24h
adult	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

Calcium concentration 24-hours urine – calculation

calcium concentration in 24-hours urine [mg/24h]	=	calcium concentration in sample of 24-hours urine [mg/dl]	×	urine volume of 24-hours urine [dl/24h]
--	---	---	---	---

QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT MC240, ACCENT S120), with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.00 mg/dl (0.00 mmol/l) - ACCENT-200
0.00 mg/dl (0.00 mmol/l) – ACCENT MC240

LoD (Limit of Detection):

0.08 mg/dl (0.02 mmol/l) - ACCENT-200
0.03 mg/dl (0.0075 mmol/l) – ACCENT MC240

LoQ (Limit of Quantitation):

0.34 mg/dl (0.085 mmol/l) - ACCENT-200
0.33 mg/dl (0.0825 mmol/l) – ACCENT MC240

Linearity:

up to 22 mg/dl (5.5 mmol/l) – ACCENT-200
up to 20.1 mg/dl (5.025 mmol/l) – ACCENT MC240

For higher concentrations dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	9.73	0.08	0.77
	level 2	12.4	0.13	1.06
ACCENT MC240 n=20	level 1	10.1	0.17	1.65
	level 2	12.5	0.06	0.45
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	9.8	0.25	2.6
	level 2	12.6	0.24	1.9
ACCENT MC240 n=80	level 1	10.3	0.20	1.9
	level 2	12.7	0.18	1.4

Method comparison

A comparison between calcium values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 60 serum samples gave following results:
y = 0.9846 x + 0.2834 mg/dl;
R = 0.975 (R – correlation coefficient)

A comparison between calcium values determined at ACCENT-200 (y) and at BECKMAN COULTER AU680 (x) using 67 urine samples gave following results:
y = 0.956 x + 0.2535 mg/dl;
R = 0.998 (R – correlation coefficient)

A comparison between calcium values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 60 serum samples gave following results:
y = 0.939 x + 0.9582 mg/dl;
R = 0.981 (R – correlation coefficient)

A comparison between calcium values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 70 urine samples gave following results:
y = 0.9233 x + 0.7828 mg/dl;
R = 0.997 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters	
Test Name	Ca ARS
Test No	59
Full Name	Calcium Arsenazo
Reference No	59
Analy. Type	Endpoint
Pri. Wave.	670 nm
Secon. Wave.	
Trend	Increase
Reac. Time	0 19
Incuba. Time	
Unit	mg/dl
Precision	0.01
R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Reag. Blank	
Concentration	0.61 22
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters	
Test Name	Ca ARS
Test No	59
Full Name	Calcium Arsenazo
Standard No	59
Reac. Type	Endpoint
Pri. Wave.	670 nm
Sec. Wave.	
Direction	Increase
Reac. Time	0 30
Incuba. Time	
Unit	mg/dl
Precision	0.01
R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Linearity Range	0.33 20.4
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters	
Test Name	Ca ARS
Test No	67
Full Name	Calcium Arsenazo
Standard No	67
Reac. Type	Endpoint
Pri. Wave.	670 nm
Sec. Wave.	
Direction	Increase
Reac. Time	0 19
Incuba. Time	16
Unit	mg/dl
Precision	0.01
R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Linearity Range	0.31 20.9
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	Ca Ars	No.	059	Sample Type	SERUM / URINE
Chemistry	Calcium Arsenazo	Print name	Ca Ars	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	670 nm	Incubation Time	0	Reaction Time	24 26
Unit	mg/dL	Blank Time	-3 -1	Reagent Vol	
Standard	Sample Vol 2 μL Aspirated μL Diluent μL	R1	200 μL	R2	
Decreased	2 μL 20 μL 180 μL	Factor			
Increased	μL μL μL	PC	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>
Linearity range (Standard)	0.31 19.6	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-40000 40000	On-board Stability	84 Day(s)
R1 Blank Abs	-40000 40000	Reagent Alarm Limit			
Blank Response	-40000 40000	Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check					
Q1		Q2		V1	
Q3		Q4		V2	
Q5		Q6		PC1	
PC2					
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
Pretreat Sample Vol		μL		Pretreat Sample Vol	
μL					
CALIBRATION SETTINGS			AUTO CALIBRATION		
Math model	Multi-point linear				<input type="checkbox"/> Bottle Changed
Factor		Replicates	2		<input type="checkbox"/> Lot Changed
ACCEPTANCE LIMITS					
Cal Time	2016	Hour			
Slope Diff		SD			
Sensitivity		Repeatability	40000		
Deter Coeff					

ACCENT-200 CALCIUM ARSENAZO

• ACCENT MC240

Chem	Ca Ars	No.	059	Sample Type	SERUM / URINE
Chemistry	Calcium Arsenazo	Print name	Ca Ars	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	660 nm	Incubation Time	0	Reaction Time	28 30
Unit	mg/dL	Blank Time	-3 -1	Standard	Sample Vol 2 μL Aspirated μL Diluent μL
		Decreased	2 μL 20 μL 180 μL	Reagent Vol	R1 200 μL R2 μL
		Increased	μL μL μL	Reagent Alarm Limit	
			<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun		
Linearity range (Standard)	0.33 20.1	Linearity Limit		Linearity Range (Decreased)	
Linearity Range (Increased)		Substrate Depletion		Linearity Range (Increased)	
R1 Blank Abs	-35000 35000	Mixed Blank Abs	-35000 35000	R1 Blank Abs	-35000 35000
Blank Response	-35000 35000	On-board Stability	84 Day(s)	Blank Response	-35000 35000
Twin Chemistry		Reagent Alarm Limit		Twin Chemistry	
<input type="checkbox"/> Prozone Check		<input type="checkbox"/> Enzyme Linear Extension		<input type="checkbox"/> Prozone Check	
Q1 Q2 V1 Q3 Q4 V2		Q5 Q6 V3 PC1 PC2		Q1 Q2 V1 Q3 Q4 V2	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment		Q5 Q6 V3 PC1 PC2	
	Pretreat Sample Vol μL	Pretreat Sample Vol μL		<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
				<input type="checkbox"/> Calibrator Pretreatment	Pretreat Sample Vol μL Pretreat Sample Vol μL
CALIBRATION SETTINGS		AUTO CALIBRATION		CALIBRATION SETTINGS	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed		Math model	Multi-point linear
Factor	Replicates 2	<input type="checkbox"/> Lot Changed		Factor	Replicates 2
		<input type="checkbox"/> Cal Time			
ACCEPTANCE LIMITS				ACCEPTANCE LIMITS	
Cal Time	2016 Hour			Cal Time	2016 Hour
Slope Diff	SD			Slope Diff	SD
Sensitivity	Repeatability 35000			Sensitivity	Repeatability 35000
Deter Coeff				Deter Coeff	

• ACCENT M320

Chem	Ca Ars	No.	059	Sample Type	SERUM / URINE
Chemistry	Calcium Arsenazo	Print name	Ca Ars	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	700 nm	Decimal	0.01
Pri Wave	660 nm	Incubation Time	0	Reaction Time	20 22
Unit	mg/dL	Blank Time	-3 -1	Standard	Sample Vol 2 μL Aspirated μL Diluent μL
		Decreased	2 μL 20 μL 180 μL	Reagent Vol	R1 200 μL R2 μL
		Increased	μL μL μL	Reagent Alarm Limit	
			<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun		
Linearity range (Standard)	0.31 20.2	Linearity Limit		Linearity Range (Decreased)	
Linearity Range (Increased)		Substrate Depletion		Linearity Range (Increased)	
R1 Blank Abs	-35000 35000	Mixed Blank Abs	-35000 35000	R1 Blank Abs	-35000 35000
Blank Response	-35000 35000	On-board Stability	84 Day(s)	Blank Response	-35000 35000
Twin Chemistry		Reagent Alarm Limit		Twin Chemistry	
<input type="checkbox"/> Prozone Check		<input type="checkbox"/> Enzyme Linear Extension		<input type="checkbox"/> Prozone Check	
Q1 Q2 V1 Q3 Q4 V2		Q5 Q6 V3 PC1 PC2		Q1 Q2 V1 Q3 Q4 V2	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment		Q5 Q6 V3 PC1 PC2	
	Pretreat Sample Vol μL	Pretreat Sample Vol μL		<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
				<input type="checkbox"/> Calibrator Pretreatment	Pretreat Sample Vol μL Pretreat Sample Vol μL
CALIBRATION SETTINGS		AUTO CALIBRATION		CALIBRATION SETTINGS	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed		Math model	Multi-point linear
Factor	Replicates 2	<input type="checkbox"/> Lot Changed		Factor	Replicates 2
		<input type="checkbox"/> Cal Time			
ACCEPTANCE LIMITS				ACCEPTANCE LIMITS	
Cal Time	2016 Hour			Cal Time	2016 Hour
Slope Diff	SD			Slope Diff	SD
Sensitivity	Repeatability 35000			Sensitivity	Repeatability 35000
Deter Coeff				Deter Coeff	

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