ACCENT-200 CALCIUM

Nr kat. 7-251

(PL)

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Cat. No 7-251 (EN)

INTENDED USE

Diagnostic kit for determination of calcium concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium and phosphorus as a hydroxyapatite constitute mineral portion of bone. Calcium occurs also as divalent cations (free or bound with negatively charged proteins) which participate in blood coagulation, neuromuscular excitability, skeletal and cardiac muscle contractility and in multiple cellular functions. Calcium flux in organism is controlled by action of parathyroid hormone (PTH), vitamin D and calcitonin. Calcium serum level abnormalities are caused usually by parathyroid or thyroid disease, disorders of vitamin D metabolism or acute pancreatitis.

METHOD PRINCIPLE

Calcium ions form a violet complex with o-cresolphthalein complexone in alkaline solution. The intensity of violet colour of this complex measured at 570-580 nm is proportional to the calcium concentration in the sample.

REAGENTS

Package

1-REAGENT 1 x 30 ml
2-REAGENT 1 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 10 weeks.

Concentrations in the test

o-cresolphthalein complexone 0.06 mmol/l
8-quinolinol 8.6 mmol/l
hydrochloric acid 30 mmol/l
ethanolamine 377 mmol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains hydrochloric acid.

Danger

H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

Ingredients:

2-REAGENT contains 2-aminoethanol.

Danger

H314 Causes severe skin burns and eye damage.
H335 May cause respiratory irritation.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

SPECIMEN

Serum, heparinized plasma free from hemolysis, 24-hours urine. Recommended anticoagulants: heparine lithium, sodium or ammonium salt.

Urine preparation: To prevent calcium salt precipitation in urine, specimens should be collected in HCl, 20-30 ml of 6M for 24-h specimen. Any specimens collected without acid should be acidified using 20-30 ml of 6M HCl, well mixed and allowed to stand for 1 h before aliquotting. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor.

Serum and plasma can be stored up to 8 hours at 15-25°C or up to 1 day at 2-8°C. Samples frozen at -20°C can be stored up to 1 year.

24-hours urine samples should be kept at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

Deionised water is recommended as a reagent blank.

For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: RF - CALCIUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES⁸

serum, plasma	mg/dl	mmol/l
premature	6.2 – 11.0	1.55 – 2.75
adults 18 – 60 yr	8.6 – 10.0	2.15 – 2.50
60 – 90 yr	8.8 – 10.2	2.20 – 2.55
> 90 yr	8.2 – 9.6	2.05 – 2.40
24-hours urine	mg/24h	mmol/24h
	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples:

CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analyser ACCENT-220S, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

For the calibration of automatic analyser BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 4 days, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or ACCENT-200 II GEN, and/or ACCENT-220S. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

2.00 mg/dl (0.50 mmol/l) - ACCENT 200
0.19 mg/dl (0.05 mmol/l) - ACCENT 200 II GEN
0.20 mg/dl (0.05 mmol/l) - ACCENT 220S

Linearity

up to 15 mg/dl (3.8 mmol/l) - ACCENT 200
up to 15 mg/dl (3.8 mmol/l) - ACCENT 200 II GEN
up to 21 mg/dl (5.3 mmol/l) - ACCENT 220S

For higher concentrations dilute the sample with 0.9% NaCl in the ratio of 1 to 1 and reassay. Multiply the result by 2.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n=20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	3.56	0.12	3.38
level 2	7.58	0.08	1.07

Reproducibility (day to day) n=20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.96	0.48	4.80
level 2	12.74	0.56	4.38

Method comparison

A comparison between calcium values determined at ACCENT-200 (y) and at COBAS INTEGRA 400 (x) using 66 samples gave following results:

y = 0.9589 x + 0.357 mg/dl;

R = 0.952 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	4
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	CALC	R1	200
No	11	R2	50
Full Name	Calcium	Sample Volume	6
Standard No	11	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	0.2 21
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-2 3	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	4
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters

Test Name	CALC	R1	200
Test No	11	R2	50
Full Name	Calcium	Sample Volume	6
Reference No	11	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	578 nm	Concentration	0.19 15
Secon. Wave.	670 nm	Linearity Limit	
Trend	Ascending	Substrate Limit	
Reac. Time	-2 4	Factor	
Incuba. Time	10	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	4
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters

Test	CALC	R1	200
No	11	R2	50
Full Name	Calcium	Sample Volume	3
Standard No	11	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	1.55 22.5
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1 8	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

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