## ACCENT-200 AMYLASE

Nr kat. 7-255 (PL)

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Cat. No **7-255** (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -amylase activity intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$ -Amylase is a digestive enzyme secreted by salivary glands and pancreas. Low level of amylase is also found in skeletal muscle, adipose tissue and fallopian tubes.  $\alpha$ -Amylase is measured generally in pancreas diseases. Elevation of amylase activity is observed also due to inflammation of abdominal cavity or salivary glands.

### METHOD PRINCIPLE

2-Chloro-4-nitrophenyl- $\alpha$ -maltotrioidide (CNP-G3) is a direct substrate for determination of  $\alpha$ -amylase activity, which does not require the presence of ancillary enzymes.

10 CNP-G3  $\alpha$ -amylase  $\rightarrow$  9 CNP + CNP-G2 + 9 maltotriose + glucose

The rate of 2-chloro-4-nitrophenol formation can be monitored at 405 nm and is proportional to the  $\alpha$ -amylase activity.

### REAGENTS

**Package**  
1-REAGENT 2 x 30 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagent stored on board of the analyzer at 2-10°C is stable for 9 weeks (ACCENT-200) or 12 weeks (ACCENT M320, ACCENT S120)..

### Concentrations in the test

MES buffer < 120 mmol/l  
calcium acetate < 7 mmol/l  
potassium hydroxide < 40 mmol/l  
potassium thiocyanate < 1100 mmol/l  
2-chloro-4-nitrophenyl- $\alpha$ -maltotrioidide < 2 mmol/l  
preservative, stabilizer

### Warnings and notes

- Protect from direct sunlight!
- Prevent the reagent from microbiological contamination and from saliva and sweat  $\alpha$ -amylase! Saliva and sweat contain  $\alpha$ -amylase. Do not pipette by the mouth, avoid skin contact with reagent, specimens, tips, cuvettes. Ensure to use automatic pipettes and laboratory gloves.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.<sup>7</sup>

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.<sup>9</sup> Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT is ready to use.  
Deionised water is recommended as a reagent blank.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: AMYLASE-MICROALBUMIN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>8</sup>

serum / plasma	U/l	$\mu$ kat/l
	20 – 104	0.34 – 1.77
urine	U/l	$\mu$ kat/l
	32 – 641	0.54 – 10.90

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 9 weeks (ACCENT-200) or 12 weeks (ACCENT M320, ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

8.4 U/l (0.14  $\mu$ kat/l) - ACCENT-200  
1.7 U/l (0.03  $\mu$ kat/l) - ACCENT MC240

#### Linearity

up to 1750 U/l (29.2  $\mu$ kat/l) - ACCENT-200  
up to 1760 U/l (29.33  $\mu$ kat/l) - ACCENT MC240

#### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=20	level 1	61.30	1.14	1.87
	level 2	358.75	5.88	1.64
ACCENT MC240 n=20	level 1	72.04	0.29	0.40
	level 2	396.67	2.06	0.52
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=80	level 1	130.74	2.86	2.18
	level 2	280.50	5.53	1.97
ACCENT MC240 n=80	level 1	72.9	2.82	3.9
	level 2	401.5	20.06	5.0

#### Method comparison

A comparison between amylase values determined at **ACCENT-200** (y) and at **BS-400** (x) using 30 samples gave following results:

$$y = 0.9977x + 0.6639 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between amylase values determined at **ACCENT MC240** (y) and at **BS-800** (x) using 60 serum samples gave following results:

$$y = 0.974x + 0.3367 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between amylase values determined at **ACCENT MC240** (y) and at **BS-800** (x) using 60 urine samples gave following results:

$$y = 0.9714x + 2.5395 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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#### Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	63
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • ACCENT-200 II GEN

#### Parameters

Test Name	AMYL	R1	250
Test No	17	R2	
Full Name	Amylase	Sample Volume	5
Reference No	17	R1 Blank	
Analy. Type	Kinetic	Mixed Reag. Blank	
Pri. Wave	405 nm	Concentration	12.5 1200
Secon. Wave	670 nm	Linearity Limit	0.2
Trend	Ascending	Substrate Limit	
Reac. Time	4 23	Factor	
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

#### Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	63
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • ACCENT-220S

#### Parameters

Test	AMYL	R1	250
No	17	R2	
Full Name	Amylase	Sample Volume	5
Standard No	17	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave	405 nm	Linearity Range	7.5 1900
Sec. Wave	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	4 23	Factor	
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

#### Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	63
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • BS-120

#### Parameters

Test	AMYL	R1	250
No	17	R2	
Full Name	Amylase	Sample Volume	5
Standard No	17	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave	405 nm	Linearity Range	8.5 2100
Sec. Wave	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	4 15	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

#### Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

### • ACCENT S120

Chem	AMY	No.	017	Sample Type	SERUM
Chemistry	AMYLASE	Print name	AMYLASE	Reaction Direction	positive
Reaction Type	Kinetic	Sec Wave	670 nm	Decimal	0.1
Pri Wave	405 nm	Incubation Time	0	Reaction Time	6 25
Unit	U/L	Standard	Sample Vol 3 Aspirated Diluent	Reagent Vol	R1 200 R2
Blank Time		Decreased	3 20 180	Increased	
		Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>
Linearity range (Standard)	3.5 1700	Linearity Limit	0.2	Substrate Depletion	40000
Linearity Range (Decreased)		Linearity Range (Increased)		Mixed Blank Abs	-40000 40000
R1 Blank Abs	-40000 40000	Blank Response	-40000 40000	On-board Stability	
Twin Chemistry		Enzyme Linear Extension	<input type="checkbox"/>	Reagent Alarm Limit	
<input type="checkbox"/> Prozone Check		O1 O2 V1 O3 O4 V2		PC1 PC2	
O5 O6 V3		Sample Pretreatment		Control Pretreatment	
Calibrator Pretreatment		Pretreat Sample Vol		Pretreat Sample Vol	
<b>CALIBRATION SETTINGS</b>	Math model Multi-point linear	<b>AUTO CALIBRATION</b>		Bottle Changed	
Factor		Lot Changed		Cal Time	
Replicates	2				
<b>ACCEPTANCE LIMITS</b>	Cal Time				
Slope Diff		SD			
Sensitivity		Repeatability	40000		
Deter Coeff					

## ACCENT-200 AMYLASE

### ACCENT MC240

Chem	AMY	No.	017	Sample Type	SERUM						
Chemistry	AMYLASE	Print name	AMYLASE	Reaction Direction	positive						
Reaction Type	Kinetic	Sec Wave	700 nm	Decimal	0.1						
Pri Wave	412 nm	Incubation Time	0	Reaction Time	10 30						
Unit	U/L	Sample Vol	3	Aspirated	20	Diluent	180	Reagent Vol	200		
Blank Time		Standard	3	Decreased	3	Increased		R1	200	R2	
		Decreased	3	Increased		Sample Blank	V	Auto Rerun			

Linearity range (Standard)	1.7	1760	Linearity Limit	0.2							
Linearity Range (Decreased)			Substrate Depletion	35000							
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			Enzyme Linear Extension								
Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment							
Pretreat Sample Vol		uL	Pretreat Sample Vol		uL						

Math model	Multi-point linear	Bottle Changed	
Factor		Lot Changed	
Replicates	2	Cal Time	

Cal Time	Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

### ACCENT M320

Chem	AMY	No.	017	Sample Type	SERUM						
Chemistry	AMYLASE	Print name	AMYLASE	Reaction Direction	positive						
Reaction Type	Kinetic	Sec Wave	700 nm	Decimal	0.1						
Pri Wave	412 nm	Incubation Time	0	Reaction Time	10 30						
Unit	U/L	Sample Vol	3	Aspirated	20	Diluent	180	Reagent Vol	200		
Blank Time		Standard	3	Decreased	3	Increased		R1	200	R2	
		Decreased	3	Increased		Sample Blank	V	Auto Rerun			

Linearity range (Standard)	3	1800	Linearity Limit	0.2							
Linearity Range (Decreased)			Substrate Depletion	35000							
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			Enzyme Linear Extension								
Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment							
Pretreat Sample Vol		uL	Pretreat Sample Vol		uL						

Math model	Multi-point linear	Bottle Changed	
Factor		Lot Changed	
Replicates	2	Cal Time	

Cal Time	Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

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