ACCENT-200 CHOLINESTERASE

Nr kat. **7-256** (PL)

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Cat. No **7-256** (EN)

INTENDED USE

Diagnostic kit for determination of cholinesterase activity intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

There are two cholinesterases (CHE and ACHE) differing in substrate specificity, tissue of origin and biological role. The term cholinesterase (ACHE), also known as acetylcholine acetylhydrolase, is found in erythrocytes, in the lungs, spleen and in grey matter of the brain. The pseudocholinesterase (CHE), also referred to as acylcholine acylhydrolase, is found in serum, the liver, pancreas, heart and in the white matter of brain. The assay of serum cholinesterase (CHE) is useful to diagnose: liver disorders, hepatitis, cirrhosis, carcinoma with metastasis, sensitivity to succinylcholine administration and pesticide poisoning. Levels decrease in all of the diseases above.

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The method uses butyrylthiocholine as the specific substrate for cholinesterase (CHE). Cholinesterase catalyses the hydrolysis of butyrylthiocholine substrate forming butyrate and thiocholine, in presence of potassium hexacyanoferrate (III). Thiocholine reduces potassium hexacyanoferrate (III) (yellow colour) to potassium hexacyanoferrate (II) (colourless). The decrease in absorbance is directly proportional to the CHE activity in the sample.

REAGENTS

Package
1-REAGENT 1 x 25.2 ml
2-REAGENT 1 x 7.1 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 11 weeks (ACCENT-200) or 12 weeks (ACCENT MC240, ACCENT S120).

Concentrations in the test

1-REAGENT
pyrophosphate buffer, pH 7.7 65 mmol/l
hexacyanoferrate (III) 2 mmol/l
2-REAGENT
Good's Buffer pH 4.0 20 mmol/l
butyrylthiocholine iodide 65 mmol/l

Warnings and notes

- Do not freeze reagents.
- Protect from direct sunlight, evaporation and avoid contamination!
- Immediately after use, recap the bottles and store at 2-8°C.
- Do not use after expiry date.
- Do not interchange caps.
- Reagents should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of reagent instability.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains tetrasodium pyrophosphate decahydrate.

Danger

H318 Causes serious eye damage.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Fresh serum free from haemolysis, plasma (EDTA, heparin) not hemolyzed.

Serum / plasma should be separated from red blood cells as soon as possible after blood collection.

Do not use citrate, borate, oxalate and fluoride as an anticoagulant because it inhibits cholinesterase activity. It is recommended to follow CLSI procedures regarding specimen collecting and handling.

Sample may be stored for up to 15 days at 2-8°C or 12 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionized water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: CHOLINESTERASE - MICROALBUMIN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ^{4,5}

serum / plasma	37°C	
female	4000 – 12600 U/l	67 – 210 µkat/l
male	5100 – 11700 U/l	85 – 195 µkat/l

In infants up to 6 months of age, cholinesterase activity is 40% to 50% higher than in adults. In young adult (< 35 years) women, the enzyme activity is approximately 64% to 74% of that in adult males. The activity decreases during pregnancy. It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, BS-120 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number.

For the calibration of automatic analyser ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 11 weeks (ACCENT-200) or 12 weeks (ACCENT MC240, ACCENT S120) with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity**
55 U/l (0.92 µkat/l) - ACCENT-200
38 U/l (0.63 µkat/l) - ACCENT MC240
- Linearity**
up to 21000 U/l (350 µkat/l) - ACCENT-200
up to 16400 U/l (273 µkat/l) - ACCENT MC240

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**
Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=10	level 1	4884	35.10	0.72
	level 2	6830	43.06	0.63
ACCENT MC240 n=20	level 1	4843	35.38	0.7
	level 2	5966	43.65	0.7
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=20	level 1	6903	245.11	3.55
	level 2	4931	156.35	3.17
ACCENT MC240 n=80	level 1	6118	210.47	3.4
	level 2	4981	175.7	3.5

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 76 samples gave following results:

$$y = 1.0545 x - 207.49 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between cholinesterase values determined at ACCENT MC240 (y) and at BECKMAN COULTER AU680 (x) using 60 samples gave following results:

$$y = 0.9315 x + 504.68 \text{ U/l};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters

Test Name	CHE	R1	200
Test No	62	R2	50
Full Name	Cholinesterase	Sample Volume	4
Reference No	62	R1 Blank	
Analy. Type	Kinetic	Mixed Reag. Blank	
Pri. Wave.	405 nm	Concentration	97 20000
Secon. Wave.	670 nm	Linearity Limit	0.2
Trend	Descending	Substrate Limit	
Reac. Time	2 17	Factor	
Incuba. Time	4	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test Name	CHE	R1	200
Test No	62	R2	50
Full Name	Cholinesterase	Sample Volume	4
Standard No	62	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	78 19000
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Decrease	Substrate Limit	
Reac. Time	6 16	Factor	
Incuba. Time	4	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters

Test Name	CHE	R1	200
Test No	62	R2	50
Full Name	Cholinesterase	Sample Volume	5
Standard No	62	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	110 15400
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Decrease	Substrate Limit	
Reac. Time	2 17	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	CHE	No.	062	Sample Type	SERUM						
Chemistry	CHOLINESTERASE			Print name	CHE						
Reaction Type	Kinetic			Reaction Direction	negative						
Pri Wave	405nm			Sec Wave	670nm						
Unit	U/L			Decimal	0						
Blank Time				Incubation Time	10						
				Reaction Time	5 15						
Standard	4	Aspirated		Diluent		Reagent Vol					
Decreased	4	20	180			R1	160				
Increased						R2	40				
Linearity range (Standard)	75	16600		Linearity Limit	0,2						
Linearity Range (Decreased)				Substrate Depletion							
Linearity Range (Increased)				Mixed Blank Abs	-40000 40000						
R1 Blank Abs	-40000	40000		On-board Stability	Day(s)						
Blank Response	-40000	40000		Reagent Alarm Limit							
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension							
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment				<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment					
				Pretreat Sample Vol		µL		Pretreat Sample Vol		µL	
CALIBRATION SETTINGS				AUTO CALIBRATION							
Math model	Two-point Linear			<input type="checkbox"/> Bottle Changed							
Factor		Replicates	2	<input type="checkbox"/> Lot Changed							
				<input type="checkbox"/> Cal Time							
ACCEPTANCE LIMITS											
Cal Time		Hour									
Slope Diff				SD							
Sensitivity		Repeatability	40000								
Deter Coeff											

ACCENT-200 CHOLINESTERASE

• ACCENT MC240

Chem	<input type="text" value="CHE"/>	No.	<input type="text" value="062"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="CHOLINESTERASE"/>	Print name	<input type="text" value="CHE"/>		
Reaction Type	<input type="text" value="Kinetic"/>	Reaction Direction	<input type="text" value="negative"/>		
Pri Wave	<input type="text" value="412nm"/>	Sec Wave	<input type="text" value="700nm"/>		
Unit	<input type="text" value="U/L"/>	Decimal	<input type="text" value="0"/>		
Blank Time	<input type="text"/>	Incubation Time	<input type="text" value="21"/>	Reaction Time	<input type="text" value="5"/> <input type="text" value="15"/>
Standard	Sample Vol <input type="text" value="5"/> μL	Aspirated	<input type="text"/>	Diluent	<input type="text"/>
Decreased	<input type="text" value="5"/> μL	<input type="text" value="20"/> μL	<input type="text" value="180"/> μL	Reagent Vol	<input type="text"/>
Increased	<input type="text"/>	<input type="text"/>	<input type="text"/>	R1	<input type="text" value="160"/> μL
				R2	<input type="text" value="40"/> μL
	<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun		

Linearity range (Standard)	<input type="text" value="38"/> <input type="text" value="16400"/>	Linearity Limit	<input type="text" value="0,2"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text"/> Day(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check			
O1 <input type="text"/>	O2 <input type="text"/>	V1 <input type="text"/>	O3 <input type="text"/>
O5 <input type="text"/>	O6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	
	Pretreat Sample Vol <input type="text"/> μL	Pretreat Sample Vol <input type="text"/> μL	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/>
Sensitivity <input type="text"/>
Deter Coeff <input type="text"/>
SD <input type="text"/>
Repeatability <input type="text" value="35000"/>

• ACCENT M320

Chem	<input type="text" value="CHE"/>	No.	<input type="text" value="062"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="CHOLINESTERASE"/>	Print name	<input type="text" value="CHE"/>		
Reaction Type	<input type="text" value="Kinetic"/>	Reaction Direction	<input type="text" value="negative"/>		
Pri Wave	<input type="text" value="412nm"/>	Sec Wave	<input type="text" value="700nm"/>		
Unit	<input type="text" value="U/L"/>	Decimal	<input type="text" value="0"/>		
Blank Time	<input type="text"/>	Incubation Time	<input type="text" value="12"/>	Reaction Time	<input type="text" value="10"/> <input type="text" value="30"/>
Standard	Sample Vol <input type="text" value="2,4"/> μL	Aspirated	<input type="text"/>	Diluent	<input type="text"/>
Decreased	<input type="text" value="2,4"/> μL	<input type="text" value="20"/> μL	<input type="text" value="180"/> μL	Reagent Vol	<input type="text"/>
Increased	<input type="text"/>	<input type="text"/>	<input type="text"/>	R1	<input type="text" value="120"/> μL
				R2	<input type="text" value="30"/> μL
	<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun		

Linearity range (Standard)	<input type="text" value="151"/> <input type="text" value="16300"/>	Linearity Limit	<input type="text" value="0,2"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text"/> Day(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check			
O1 <input type="text"/>	O2 <input type="text"/>	V1 <input type="text"/>	O3 <input type="text"/>
O5 <input type="text"/>	O6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	
	Pretreat Sample Vol <input type="text"/> μL	Pretreat Sample Vol <input type="text"/> μL	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/>
Sensitivity <input type="text"/>
Deter Coeff <input type="text"/>
SD <input type="text"/>
Repeatability <input type="text" value="35000"/>

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