## ACCENT-200 UA PLUS

Nr kat. 7-263 (PL)

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(EN)

### INTENDED USE


Diagnostic kit with ascorbate oxidase for determination of uric acid concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN and ACCENT-220S, ACCENT S120, ACCENT MC240 and ACCENT M320. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculosis. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

### METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

### REAGENTS

#### Package

1-Reagent 2 x 31 ml  
2-Reagent 1 x 16 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

### Concentrations in the test

buffer PIPES (pH 7.0) 100 mmol/l  
4-aminoantipyrine 0.78 mmol/l  
ADPS 0.67 mmol/l  
ferricyanide potassium 3.8 μmol/l  
peroxidase (POD) > 38.34 μkat/l  
uricase > 1.65 μkat/l  
ascorbate oxidase > 66.7 μkat/l  
sodium hydroxide < 1 %  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

24- hours urine, serum, heparinized plasma free from hemolysis.

Do not use EDTA, fluoride and oxalate as anticoagulants! Urine preparation: To prevent precipitation of salts of uric acid, 10 ml of NaOH (500 g/L) should be added to the collection bottle before collection of a 24-hour specimen. Urine should be diluted with distilled water in the ratio of 1 to 4 (multiply the result by 5).

Serum and plasma can be stored 3-5 days at 2-8°C or 6 months at -20°C. 24-hours urine samples can be stored approximately 3 days at room temperature.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

#### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S there is a probability of **cross-contamination** affecting the tests results: ALBUMIN – UA PLUS, CREATININE – UA PLUS, HDL DIRECT II GEN – UA PLUS, LDL DIRECT II GEN – UA PLUS, URINE PROTEINS II GEN – UA PLUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>5</sup>

serum / plasma	mg/dl	μmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458
24-hours urine	mg/24h	mmol/24h
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples:

CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) - for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) - for determination in urine. For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240 and ACCENT M320, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or ACCENT-200 II GEN and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

0.31 mg/dl (18.4 μmol/l) - ACCENT-200

#### LOQ

0.16 mg/dl (9.5 μmol/l) – ACCENT MC240

#### Linearity

up to 26.5 mg/dl (1576 μmol/l) - ACCENT-200

up to 42 mg/dl (2498.2 μmol/l) – ACCENT MC240

For higher concentration of uric acid in serum or plasma, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 II GEN n=20	level 1	4.74	0.03	0.72
	level 2	8.70	0.05	0.53
ACCENT MC240 n=20	level 1	4.97	0.04	0.77
	level 2	9.62	0.05	0.55
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 II GEN n=80	level 1	4.65	0.13	2.76
	level 2	8.56	0.28	3.30
ACCENT MC240 n=80	level 1	4.9	0.07	1.4
	level 2	9.5	0.13	1.4

#### Method comparison

A comparison between uric acid values determined at ACCENT-200 (y) and at ADVIA 1650 (x) using 30 serum samples gave following results:

y = 1.0119 x + 0.041 mg/dl;

R = 0.998 (R – correlation coefficient)

A comparison between uric acid values determined at ACCENT-200 (y) and at ADVIA 1650 (x) using 22 urine samples gave following results:

y = 0.8454 x + 1.4363 mg/dl;

R = 0.991 (R – correlation coefficient)

A comparison between uric acid values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 58 serum samples gave following results:

y = 0.9768 x + 0.2281 mg/dl;

R = 0.998 (R – correlation coefficient)

A comparison between uric acid values determined at ACCENT MC240 (y) and at BS-800 (x) using 65 urine samples gave following results:

y = 1.0165 x – 0.3591 mg/dl;

R = 0.998 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 UA PLUS

### • ACCENT MC240

Chem	UA PLUS		No.	058		Sample Type	SERUM/URINE	
Chemistry	UA PLUS		Print name	UA PLUS		Reaction Direction	positive	
Reaction Type	Endpoint		Pri Wave	546 nm		Sec Wave	660 nm	
Unit	mg/dL		Decimal	0.1		Incubation Time	21	
Blank Time	-3	-1	Reaction Time	19	21	Reagent Vol	R1 160 µL R2 40 µL	
Standard	Sample Vol	Aspirated	Diluent			Reagent Vol	R1 160 µL R2 40 µL	
Decreased	3.2	20	180			Reagent Vol	R1 160 µL R2 40 µL	
Increased						Reagent Vol	R1 160 µL R2 40 µL	
	<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun					

Linearity range (Standard)	0.16	42	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000	
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)	
Blank Response	-35000	35000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check					
Q1	Q2	V1	Q3	Q4	V2
O5	O6	V3	PC1	PC2	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
Pretreat Sample Vol	µL	Pretreat Sample Vol	µL		

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor	Replicates 2	<input type="checkbox"/> Lot Changed
		<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

### • ACCENT M320

Chem	UA PLUS		No.	058		Sample Type	SERUM/URINE	
Chemistry	UA PLUS		Print name	UA PLUS		Reaction Direction	positive	
Reaction Type	Endpoint		Pri Wave	546 nm		Sec Wave	700 nm	
Unit	mg/dL		Decimal	0.1		Incubation Time	20	
Blank Time	-3	-1	Reaction Time	15	17	Reagent Vol	R1 120 µL R2 30 µL	
Standard	Sample Vol	Aspirated	Diluent			Reagent Vol	R1 120 µL R2 30 µL	
Decreased	2.4	20	180			Reagent Vol	R1 120 µL R2 30 µL	
Increased						Reagent Vol	R1 120 µL R2 30 µL	
	<input type="checkbox"/> Sample Blank		<input checked="" type="checkbox"/> Auto Rerun					

Linearity range (Standard)	0.16	59	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000	
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)	
Blank Response	-35000	35000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check					
Q1	Q2	V1	Q3	Q4	V2
O5	O6	V3	PC1	PC2	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
Pretreat Sample Vol	µL	Pretreat Sample Vol	µL		

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor	Replicates 2	<input type="checkbox"/> Lot Changed
		<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).