## ACCENT-200 dTIBC

Nr kat. **7-265** (PL)

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Cat. No **7-265** (EN)

### INTENDED USE

Diagnostic kit for determination of total iron binding capacity intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240 and ACCENT M320. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION


Iron is transported as Fe (III) bound to the apotransferrin. The apotransferrin-Fe (III) complex is called transferrin. Normally only about one third of the iron binding sites of transferrin are occupied by Fe (III). TIBC (Total iron-binding capacity) is a measurement for the maximum iron concentration that transferrin can bind. Serum TIBC levels vary in disorders of iron metabolism, e.g. elevated TIBC level when iron levels are low is useful in the diagnosis and monitoring of iron deficiency anaemia and also anaemia in late pregnancy. Decreased level of TIBC may indicate disorders such as anaemia associated with chronic inflammatory disorders, haemochromatosis, other anaemias, thalassemia, cirrhosis and malignancies.

### METHOD PRINCIPLE


Colorimetric method, direct, with chromazurol B.

The assay consists of 2 reaction steps:

1. The sample is added to 1-Reagent acidic buffer which contains excess iron and iron-binding dye (Chromazurol B). The low pH causes transferrin to release all iron atoms which then are bound by the iron-binding dye along with the excess iron, producing a coloured dye-iron complex.

2. Neutral buffer 2-Reagent addition results in pH change and causes a large increase in affinity of transferrin for iron. Then serum transferrin rapidly binding iron ions separating it from the dye-iron complex. The observed decrease in the absorbance of a coloured dye-iron complex is directly proportional to the total iron-binding capacity TIBC in the serum sample.

### REAGENTS

**Package**  
1-Reagent 1 x 23 ml  
2-Reagent 1 x 6 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 6 weeks on board the analyzer.

### Concentrations in the test

<b>1-Reagent</b>	
sodium acetate	200 mmol/l
iron chloride	0.02 mmol/l
chromazurol B	< 0.25 mmol/l
hexadecyltrimethylammonium bromide (CTAB)	< 3.0 mmol/l
preservatives, conservants	
<b>2-Reagent</b>	
sodium bicarbonate	325 mmol/l
buffer, preservatives, conservants	

### Warnings and notes

- Do not freeze the reagents.
- Protect from direct sunlight, air and avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended.
- Because of the great quantities of iron in the environment, scrupulous care is necessary to ensure that glassware, water and reagents do not become contaminated with iron.
- For internal quality control non frozen control serum is recommended.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum free from haemolysis. Presence of fibrin in the serum sample caused by incomplete clot formation causes non-specific interference with TIBC determination.

Do not use plasma samples as anticoagulants affect TIBC test results.

Serum can be stored up to 3 days at 4°C or up to 6 months at -20°C or also at -70°C indefinitely. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S there is a probability of **cross-contamination** affecting the tests results. When it is possible, the tests for determination of total iron binding capacity, using the **ACCENT-200 dTIBC**, should be performed **in the separate order** (follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER).

### REFERENCE VALUES <sup>2</sup>

serum	µg/dl	µmol/l
adults	250 – 425	44.8 – 71.6

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY MULTICALIBRATOR LEVEL 1 (Nr kat. 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is **recommended depending on the calibrator lot number**. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 6 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument is used.

- Sensitivity**  
101 µg/dl (18.08 µmol/l) - ACCENT-200
- LOQ**  
110 µg/dl (19.7 µmol/l) – ACCENT MC240
- Linearity**  
up to 600 µg/dl (107.4 µmol/l) - ACCENT-200  
up to 400 µg/dl (71.6 µmol/l) – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 31 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl, copper up to 3 mg/dl, zinc up to 250 µg/dl, iron up to 357.5 µg/dl, Desferal up to 11.5 µg/ml, Cuprimine up to 250 µg/dl and Iron Dextran (Imferon) up to 1430 µg/ml do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [µg/dl]	SD [µg/dl]	CV [%]
ACCENT-200 n=10	level 1	226	2.24	0.99
	level 2	341	3.30	0.97
ACCENT MC240 n=20	level 1	224.87	1.99	0.88
	level 2	322.59	1.91	0.59
Reproducibility (day to day)		Mean [µg/dl]	SD [µg/dl]	CV [%]
ACCENT-200 n=10	level 1	362	4.31	1.19
	level 2	248	4.73	1.91
ACCENT MC240 n=80	level 1	229.6	3.19	1.4
	level 2	325.1	6.82	2.1

### Method comparison

A comparison between TIBC values determined at **ACCENT-200** (y) and at **ADVIA 1650** (x) using 49 samples gave following results:

$$y = 0.9454 x + 9.7886 \mu\text{g/dl};$$

$$R = 0.952 \quad (R - \text{correlation coefficient})$$

A comparison between TIBC values determined at **ACCENT MC240** (y) and at **BS-400** (x) using 60 serum samples gave following results:

$$y = 1.0578 x - 11.677 \mu\text{g/dl};$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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## ACCENT-200 dTIBC

### • ACCENT MC240

Chem <input type="text" value="dTIBC"/>	No. <input type="text" value="065"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="dTIBC"/>	Print name <input type="text" value="dTIBC"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="negative"/>	
Pri Wave <input type="text" value="660 nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="µg/dL"/>	Decimal <input type="text" value="0"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="28"/> <input type="text" value="30"/>
Standard <input type="text" value="9"/> µL <input type="text" value="Aspirated"/> µL <input type="text" value="Diluent"/> µL	Reagent Vol <input type="text" value="120"/> µL	R1 <input type="text" value="120"/> µL
Decreased <input type="text" value="9"/> µL <input type="text" value="20"/> µL <input type="text" value="180"/> µL	R2 <input type="text" value="30"/> µL	
Increased <input type="text"/> µL <input type="text"/> µL <input type="text"/> µL		
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="110"/> <input type="text" value="400"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
O5 <input type="text"/> O6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Calibrator Pretreatment	
Pretreat Sample Vol <input type="text"/> µL	Pretreat Sample Vol <input type="text"/> µL

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/> Replicates <input type="text" value="2"/>	<input type="checkbox"/> Lot Changed
	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

### • ACCENT M320

Chem <input type="text" value="dTIBC"/>	No. <input type="text" value="065"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="dTIBC"/>	Print name <input type="text" value="dTIBC"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="negative"/>	
Pri Wave <input type="text" value="660 nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="µg/dL"/>	Decimal <input type="text" value="0"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="25"/>	Reaction Time <input type="text" value="23"/> <input type="text" value="25"/>
Standard <input type="text" value="9"/> µL <input type="text" value="Aspirated"/> µL <input type="text" value="Diluent"/> µL	Reagent Vol <input type="text" value="120"/> µL	R1 <input type="text" value="120"/> µL
Decreased <input type="text" value="9"/> µL <input type="text" value="20"/> µL <input type="text" value="180"/> µL	R2 <input type="text" value="30"/> µL	
Increased <input type="text"/> µL <input type="text"/> µL <input type="text"/> µL		
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="117"/> <input type="text" value="590"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
O5 <input type="text"/> O6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Calibrator Pretreatment	
Pretreat Sample Vol <input type="text"/> µL	Pretreat Sample Vol <input type="text"/> µL

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/> Replicates <input type="text" value="2"/>	<input type="checkbox"/> Lot Changed
	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

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