## ACCENT-200 HDL DIRECT

Nr kat. **7-279** (PL)

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### LITERATURA

- Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl: 1 4-13.
- Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85: 1234-41.
- Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
- Camps, J, Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol. Clinical Chemistry, 1999; 45:685-688

Data wydania: 06. 2023.

## ACCENT-200 HDL DIRECT

Cat. No **7-279** (EN)

### INTENDED USE

Diagnostic kit for determination of HDL-cholesterol concentration (direct method) used in automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCEN MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.


### INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL). The principle role of HDL in lipid metabolism is the uptake and transport of cholesterol from peripheral tissues to the liver. Low HDL cholesterol (HDL-C) levels are strongly associated with an increased risk of coronary artery disease. The HDL-C determination is used to identify high-risk patients.

### METHOD PRINCIPLE

The assay is a homogeneous method for directly measuring HDL-cholesterol concentration in serum or plasma, without any off-line pretreatment or centrifugation steps. Accelerator selective detergent methodology. During the first phase, LDL, VLDL particles and Chylomicrons generate free non-HDL cholesterol, which through an enzymatic reaction, produce hydrogen peroxide. The generated peroxide is consumed by a peroxidase reaction with DSBmT yielding a colourless product.

During the second phase, specific detergent solubilises HDL-Cholesterol. In conjunction with cholesterol oxidase (CO) and cholesterol esterase (CE) action, peroxidase and 4-AAP develop a coloured reaction which is proportional to HDL-Cholesterol concentration.

### REAGENTS

#### Package

1-REAGENT 2 x 25 ml  
2-REAGENT 2 x 9 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks (ACCENT-200, ACCENT M320).

### Concentrations in the test

#### 1-REAGENT

Buffer  
Cholesterol oxidase (*E.coli*) < 1000 U/l  
Peroxidase (horseradish) < 1300 ppg U/l  
N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT) < 1 mM  
Accelerator < 1 mM  
Preservative < 0.06 %  
Ascorbic acid oxidase (*Curcubita* spp.) < 3000 U/l

#### 2-REAGENT

Buffer  
Cholesterol esterase (*Pseudomonas* spp.)  
4-aminoantipyrine (4-AAP)  
Detergent  
Preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

#### Warning

 H317 May cause an allergic skin reaction.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

### SPECIMEN

Serum, heparinized or EDTA plasma.  
Anticoagulants containing citrate should not be used.  
Blood should be collected only if the patient has been fasting for 12-14 hours.

Serum: Collect whole blood by venepuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).

Plasma: Specimens may be collected in EDTA or lithium or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).

Serum and plasma should not remain at 15-30°C longer than 14 hours. If assays are not completed within 14 hours, serum or plasma should be stored at 2-8°C for up to 1 week. If specimens need to be stored for more than 1 week, they may be preserved at less than -20°C for up to 3 months. Samples may be frozen once.

Nevertheless it is recommended to perform the assay with freshly collected samples.

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionized water is recommended as a reagent blank. For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT II GEN - LIPASE, HDL DIRECT II GEN - TG, HDL DIRECT II GEN - TG mono, HDL DIRECT II GEN - UA, HDL DIRECT II GEN - UA PLUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

serum / plasma	40 – 60 mg/dl 1.04 – 1.55 mmol/l
----------------	-------------------------------------

As HDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT, the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended. Deionised water should be used as a calibrator 0.

For the calibration of automatic analyser BS-120 the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended.

The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT M320) with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

0.5 mg/dl (0.01 mmol/l) - ACCENT-200

#### LOQ

2.0 mg/dl (0.05 mmol/l) – ACCENT MC240

#### Linearity

up to 200 mg/dl (5.2 mmol/l) - ACCENT-200  
up to 200 mg/dl (5.2 mmol/l) – ACCENT MC240

ACCENT-200 HDL DIRECT (II GENERACJA/ II GENERATION/ II ПОКОЛЕНИЕ)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Bilirubin conjugated up to 60 mg/dl, bilirubin total up to 60 mg/dl, haemoglobin up to 1 g/dl, ascorbic acid up to 100 mg/dl, Intralipid up to 1800 mg/dl, triglycerides up to 2000 mg/dl and gamma-globulins up to 5000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	42.07	0.71	1.69
	level 2	55.03	0.78	1.41
ACCENT MC240 n=20	level 1	69.86	0.74	1.06
	level 2	28.96	0.12	0.42
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	42.12	0.59	1.40
	level 2	60.35	1.35	2.24
ACCENT MC240 n=80	level 1	28.8	0.56	1.9
	level 2	73.1	1.71	2.3

### Method comparison

A comparison between HDL cholesterol values determined at **ACCENT-200** (y) and at **ADVIA 1650** (x) using 54 samples gave following results:

$y = 0.8381x + 4.7706$  mg/dl;

$R = 0.978$  (R – correlation coefficient)

A comparison between HDL cholesterol values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 62 serum samples gave following results:

$y = 0.8628x + 4.4527$  mg/dl;

$R = 0.986$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl:1 4-13.
- Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85:1234-41.
- Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
- Camps, J. Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol, Clinical Chemistry, 1999; 45:685-688.

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#### Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • ACCENT-220S

<b>Parameters</b>	
Test	HDL D
No	25
Full Name	HDL Direct
Standard No	25
Reac. Type	Endpoint
Pri. Wave.	578 nm
Sec. Wave.	670 nm
Direction	Increase
Reac. Time	0   17
Incuba. Time	11
Unit	mg/dl
Precision	0.1
R1	180
R2	60
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Linearity Range	0.6   200
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

#### Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • ACCENT-200 II GEN

<b>Parameters</b>	
Test Name	HDL D
Test No	25
Full Name	HDL Direct
Reference No	25
Analy. Type	Endpoint
Pri. Wave.	578 nm
Secon. Wave.	670 nm
Trend	Increase
Reac. Time	0   15
Incuba. Time	10
Unit	mg/dl
Precision	0.1
R1	180
R2	60
Sample Volume	3
R1 Blank	
Mixed Reag. Blank	
Concentration	0.8   149
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

#### Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • BS-120

<b>Parameters</b>	
Test	HDL D
No	25
Full Name	HDL Direct
Standard No	25
Reac. Type	Endpoint
Pri. Wave.	578 nm
Sec. Wave.	670 nm
Direction	Increase
Reac. Time	4   18
Incuba. Time	16
Unit	mg/dl
Precision	0.1
R1	210
R2	70
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Linearity Range	0.9   340
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

#### Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • ACCENT S120

Chem	HDL D	No.	025	Sample Type	SERUM		
Chemistry	HDL DIRECT	Print name	HDL D	Reaction Direction	positive		
Reaction Type	Endpoint	Sec Wave	670 nm	Decimal	0.1		
Pri Wave	578 nm	Incubation Time	11	Reaction Time	16   18		
Unit	mg/dL	Blank Time	-3   -1	Standard	Sample Vol 2.5 μL Aspirated μL Diluent μL		
Decreased	2.5 μL	20 μL	180 μL	Reagent Vol	R1 150 μL R2 50 μL		
Increased	μL	μL	μL	Sample Blank	<input type="checkbox"/> V Auto Rerun		
Linearity range (Standard)	1.0	230	Linearity Limit				
Linearity Range (Decreased)			Substrate Depletion				
Linearity Range (Increased)			Mixed Blank Abs	-40000   40000			
R1 Blank Abs	-40000	40000	On-board Stability	Day(s)			
Blank Response	-40000	40000	Reagent Alarm Limit				
Twin Chemistry			Enzyme Linear Extension	<input type="checkbox"/>			
<input type="checkbox"/> Prozone Check		Q1	Q2	V1	Q3	Q4	V2
O5	O6	V3	PC1	PC2	Calibrator Pretreatment	Pretreat Sample Vol	μL
Sample Pretreatment	Control Pretreatment	Calibrator Pretreatment	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL	
<b>CALIBRATION SETTINGS</b>				<b>AUTO CALIBRATION</b>			
Math model	Two-point linear					Bottle Changed	
Factor		Replicates	2			Lot Changed	
						Cal Time	
<b>ACCEPTANCE LIMITS</b>							
Cal Time	Hour						
Slope Diff		SD					
Sensitivity		Repeatability	40000				
Deter Coeff							

## ACCENT-200 HDL DIRECT

### • ACCENT MC240

Chem	<input type="text" value="HDL D"/>	No.	<input type="text" value="025"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="HDL DIRECT"/>	Print name	<input type="text" value="HDL D"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Sec Wave	<input type="text" value="800 nm"/>	Decimal	<input type="text" value="0.1"/>
Pri Wave	<input type="text" value="605 nm"/>	Incubation Time	<input type="text" value="21"/>	Reaction Time	<input type="text" value="18"/> <input type="text" value="20"/>
Unit	<input type="text" value="mg/dL"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Standard	<input type="text" value="3.3"/> <input type="text" value="3.3"/> <input type="text" value="3.3"/>
Sample Vol	<input type="text" value="3.3"/> <input type="text" value="3.3"/>	Aspirated	<input type="text" value="20"/>	Diluent	<input type="text" value="180"/>
Decreased	<input type="text" value="3.3"/> <input type="text" value="3.3"/>	Reagent Vol	<input type="text" value="150"/>	R1	<input type="text" value="50"/> <input type="text" value="50"/>
Increased	<input type="text" value="3.3"/> <input type="text" value="3.3"/>	R2	<input type="text" value="50"/>	Sample Blank	<input type="checkbox"/>
		Auto Rerun	<input checked="" type="checkbox"/>		

Linearity range (Standard)	<input type="text" value="2.0"/> <input type="text" value="200"/>	Linearity Limit	<input type="text"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text"/> Day(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	Enzyme Linear Extension	<input type="checkbox"/>
Prozone Check	<input type="checkbox"/>		
O1	<input type="text"/>	O2	<input type="text"/>
O3	<input type="text"/>	O4	<input type="text"/>
O5	<input type="text"/>	O6	<input type="text"/>
Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>
Calibrator Pretreatment	<input type="checkbox"/>	Pretreat Sample Vol	<input type="text"/> $\mu$ L

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>		
Math model	<input type="text" value="Two-point linear"/>	Bottle Changed	<input type="checkbox"/>
Factor	<input type="text"/>	Lot Changed	<input type="checkbox"/>
Replicates	<input type="text" value="2"/>	Cal Time	<input type="text"/>

<b>ACCEPTANCE LIMITS</b>	
Cal Time	<input type="text"/> Hour
Slope Diff	<input type="text"/>
Sensitivity	<input type="text"/>
Repeatability	<input type="text" value="35000"/>
Deter Coeff	<input type="text"/>

### • ACCENT M320

Chem	<input type="text" value="HDL D"/>	No.	<input type="text" value="025"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="HDL DIRECT"/>	Print name	<input type="text" value="HDL D"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Sec Wave	<input type="text" value="800 nm"/>	Decimal	<input type="text" value="0.1"/>
Pri Wave	<input type="text" value="605 nm"/>	Incubation Time	<input type="text" value="20"/>	Reaction Time	<input type="text" value="28"/> <input type="text" value="30"/>
Unit	<input type="text" value="mg/dL"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Standard	<input type="text" value="3.3"/> <input type="text" value="3.3"/> <input type="text" value="3.3"/>
Sample Vol	<input type="text" value="3.3"/> <input type="text" value="3.3"/>	Aspirated	<input type="text" value="20"/>	Diluent	<input type="text" value="180"/>
Decreased	<input type="text" value="3.3"/> <input type="text" value="3.3"/>	Reagent Vol	<input type="text" value="150"/>	R1	<input type="text" value="50"/> <input type="text" value="50"/>
Increased	<input type="text" value="3.3"/> <input type="text" value="3.3"/>	R2	<input type="text" value="50"/>	Sample Blank	<input type="checkbox"/>
		Auto Rerun	<input checked="" type="checkbox"/>		

Linearity range (Standard)	<input type="text" value="1.9"/> <input type="text" value="179"/>	Linearity Limit	<input type="text"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text"/> Day(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	Enzyme Linear Extension	<input type="checkbox"/>
Prozone Check	<input type="checkbox"/>		
Q1	<input type="text"/>	Q2	<input type="text"/>
Q3	<input type="text"/>	Q4	<input type="text"/>
Q5	<input type="text"/>	Q6	<input type="text"/>
Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>
Calibrator Pretreatment	<input type="checkbox"/>	Pretreat Sample Vol	<input type="text"/> $\mu$ L

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>		
Math model	<input type="text" value="Two-point linear"/>	Bottle Changed	<input type="checkbox"/>
Factor	<input type="text"/>	Lot Changed	<input type="checkbox"/>
Replicates	<input type="text" value="2"/>	Cal Time	<input type="text"/>

<b>ACCEPTANCE LIMITS</b>	
Cal Time	<input type="text"/> Hour
Slope Diff	<input type="text"/>
Sensitivity	<input type="text"/>
Repeatability	<input type="text" value="35000"/>
Deter Coeff	<input type="text"/>

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