ACCENT-300 GLUCOSE

Nr kat. 7-301 (PL)

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LITERATURA

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2017, Diabetologia Kliniczna, tom 3, suplement A, 2017.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

Data wydania: 10.2020.

ACCENT-300 GLUCOSE

Cat. No **7-301** (EN)

INTENDED USE


Diagnostic kit for determination of glucose concentration, used in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Glucose is a simple six-carbon sugar. Oxidative metabolism of glucose provides the energy for most cellular processes. Glucose level in the blood is tightly controlled by several hormones. Elevated glucose level is the classic sign of diabetes mellitus. Glucose level abnormalities (hyper- or hypoglycemia) might be caused also by pancreas tumors and diseases of liver, thyroid gland or adrenal glands.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glucose oxidase.

The colour intensity is proportional to the glucose concentration.

REAGENTS

Package
 1-Reagent 4 x 40 ml

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

Concentrations in the test

phosphate buffer (pH 7.0) < 240 mmol/l
 phenol < 6 mmol/l
 glucose oxidase (GOD) < 480 μkat/l
 peroxidase (POD) < 44 μkat/l
 4-aminoantipyrine (4-AA) < 0,9 mmol/l
 stabilizers and preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination!

SPECIMEN

EDTA or heparinized plasma in tubes containing sodium fluoride or sodium iodoacetate additive/ serum, free from hemolysis, cerebrospinal fluid.

Plasma / Serum. Serum and plasma specimens should be separated from cells within 30 minutes after collection. Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level.

ACCENT-300 GLUCOSE

Serum and plasma can be stored up to 2 days at 4°C.³ Plasma is the specimen recommended for the glucose determination in the blood.⁵

Cerebrospinal fluid. Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.⁴

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.
 For reagent blank deionized water is recommended.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE – CREA ENZYMATIC. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT_300_CARRYOVER.

REFERENCE VALUES

	mg/dl	mmol/l
plasma, serum ^{5,6,7}	70 – 99	3.9 – 5.5
cerebrospinal fluid ⁸	40 – 70	2.2 – 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and/or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

ATTENTION: When the biological material is plasma collected in EDTA or heparin, for the calibration of automatic analysers systems only CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No. 5-175; 5-175) is recommended.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 13.7 mg/dl (0.754 mmol/l).
- Linearity:** up to 500 mg/dl (27.5 mmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

str. / page / crp. 3/7

Specificity / Interferences

Haemoglobin up to 2.50 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	107.20	1.16	1.08
level 2	247.10	1.77	0.71
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	117.17	4.53	3.87
level 2	285.44	6.99	2.45

Method comparison

A comparison between glucose values determined at ACCENT-300 (y) and at COBAS INTEGRA 400 (x) using 100 serum samples gave following results:

$$y = 1.0458x + 0.8661 \text{ mg/dl};$$

$$R = 0.9901 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 23 plasma samples gave following results:

$$y = 1.0615x - 1.2848 \text{ mg/dl};$$

$$R = 0.9965 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecena kliniczne dotyczące postępowania u chorych na cukrzycę 2017, Diabetologia Kliniczna, tom 3, suplement A, 2017.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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PROGRAM NA ANALIZATOR/ APPLICATION/ АДАПТАЦИЯ:

Parameters

No.	1	Prim.Wave.	510
Test	GLUC	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	300
Unit	mg/dl	R2 Vol.	
Decimals	1	Line. Limit	
Incubation	0	Antigen Check	
Reaction	0 40	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	13.7
Upper	2.5	Upper	500
Sample Vol.	45	Full Name	Glucose
Dilution	5	Print No.	1

Calibration

Rule	Two Point Linear / One Point Linear	
K Factor	0	
Replicates	3	
Interval	77	
Sensitivity	0	
Correlation	0	
Difference	2.5	
Blank Response	0 2.5	
Coefficient	0	
Difference		
Non-linear SD	0	

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