ACCENT-300 UREA

Nr kat. 7-306

(PL)

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(EN)

INTENDED USE



Diagnostic kit for determination of urea concentration, intended to use in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Urea is a product of amino acids catabolism. It is produced in liver and excreted in urine. Urea in the blood is reported as the blood urea nitrogen (BUN). Increased urea concentration in the serum, called uremia, is observed due to dehydration, renal failure, high-protein diet, increased protein catabolism caused by tissue injury or massive bleeding into the alimentary tract. The reason of reduced urea level could be overhydration, low-protein diet or starvation and severe liver disease.

METHOD PRINCIPLE

Kinetic, enzymatic method with urease and glutamate dehydrogenase.

The rate of absorbance changing at $\lambda=340$ nm is proportional to the urea concentration.

REAGENTS

Package

1-Reagent 4 x 40 ml
2-Reagent 2 x 20,5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the reagent

1-Reagent

Tris (pH 7.8) ≤ 144 mmol/l
ADP ≤ 0.84 mmol/l
urease ≤ 250 $\mu\text{kat/l}$
GLDH ≤ 10.5 $\mu\text{kat/l}$
stabilizers, detergents, preservatives

2-Reagent

2-oxoglutarate ≤ 48.6 mmol/l
NADH ≤ 1.6 mmol/l
buffer, preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ^{9,10,11}

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine.

Do not use heparine ammonium salt and fluoride as anticoagulants.

Specimen can be stored up to 7 days at 2-8°C.

Urine preparation: Samples with visible turbidity or the presence of precipitates should be pre-centrifuged. Urine sample should be mixed well before analysis, diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Bacterial growth in the specimen may cause erroneously elevated results.

24-hours urine samples should be adjusted to pH < 7 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionized water is recommended as a reagent blank.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: UREA - FERRITIN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT-300_CARRYOVER.

REFERENCE VALUES ⁸

serum / plasma	mg/dl	mmol/l
	< 50	< 8,3
24-hours urine	g/24h	mmol/24h
	20 – 35	300 – 550

1 mg of urea corresponds to 0.467 mg of urea nitrogen (BUN). It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.6 mg/dl (0.1 mmol/l)

LoD (Limit of Detection):

1.0 mg/dl (0.17 mmol/l)

LoQ (Limit of Quantitation):

2.5 mg/dl (0.42 mmol/l)

Linearity:

up to 300 mg/dl (49.8 mmol/l)

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	33.4	0.38	1.1
level 2	102.1	0.93	0.9

Method comparison

A comparison between urea values determined at ACCENT-300 (y) and at BECKMAN COULTER AU680 (x) using 60 serum samples gave following results:

$$y = 1.0309 x - 0.2445 \text{ mg/dl};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between urea values determined at ACCENT-300 (y) and at BECKMAN COULTER AU680 (x) using 30 plasma samples gave following results:

$$y = 1.0435 x - 1.3412 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between urea values determined at ACCENT-300 (y) and at BECKMAN COULTER AU680 (x) using 30 urine samples gave following results:

$$y = 0.9803 x + 56.289 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION for / АДАПТАЦИЯ для:

Parameters

No.	10	Prim.Wave.	340
Test	UREA	Sec.Wave.	450
Method	Kinetic	Sample Vol.	3
Direction	Descend	R1 Vol.	200
Unit	mg/dl	R2 Vol.	50
Decimals	1	Line. Limit	20
Incubation	10	Antigen Check	
Reaction	2 9	Substrat	
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	2.5
Upper	2.5	Upper	300
Sample Vol.	45	Full Name	UREA
Dilution	5	Print No.	10

Calibration

Rule	Multi-point Linear
K Factor	0
Replicates	2
Interval	
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient Difference	0
Non-linear SD	0

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