## ACCENT-300 APOLIPOPROTEIN B

Nr kat. 7-307 (PL)

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### LITERATURA

- Marcovina, S.M., Albers, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A1 and B. III.
- Tietz, N.W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 146, (2006).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 1802, (1999).

Data wydania: 10. 2020.

## ACCENT-300 APOLIPOPROTEIN B

Cat. No **7-307** (EN)

### INTENDED USE

Diagnostic kit for determination of apolipoprotein B concentration used in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lipids are transported in serum under the form of micelle known as lipoproteins. Lipoproteins are macromolecular complexes containing proteins (**apolipoproteins**), cholesterol and phospholipids in their outer layer, triglycerides and cholesterol esters in their inner structure. Lipoproteins are classified according to their increasing density respectively as chylomicrons, very low density lipoproteins (VLDL), low density lipoproteins (LDL) and high density lipoproteins (HDL). Apolipoprotein B is the major protein moiety of LDL. Apolipoprotein measurements are more discriminating than HDL and LDL cholesterol measurements in allowing to identify patients with coronary heart diseases and in assessing atherosclerotic risk.

### METHOD PRINCIPLE

The apolipoprotein B presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at  $\lambda=340$  nm is proportional to apolipoprotein B concentration in the sample.

### REAGENTS

**Package**  
1-REAGENT 1 x 40 ml  
2-REAGENT 1 x 9 ml

Buffer (1-REAGENT) stored at 2-8°C and antiserum (2-REAGENT) stored at 2-8°C are stable until expiry date printed on the package.

### Reagent components

**1-Reagent:** TRIS buffer, polymer, inorganic salt, preservative.

**2-Reagent:** anti-human apolipoprotein B antiserum, buffer, inorganic salt, preservative.

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN<sup>3</sup>

Serum or plasma.  
Serum and plasma can be stored 4-7 days at 4°C or at least 6 months at -70°C.

It is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.  
For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES<sup>3</sup>

children (4-11 years)	0.56 – 1.13 g/l
children (12-19 years)	0.55 – 1.19 g/l
adults	0.59 – 1.73 g/l

It is recommended for each laboratory to establish its own reference ranges for local population. These ranges are sex and age dependent.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY APOLIPOPROTEIN CONTROL (Cat. No 4-293).

For the calibration of automatic analysers systems the CORMAY APOLIPOPROTEIN CALIBRATORS (Cat. No 4-289) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument is used.

- Measurement range:** 0.019 g/l – 4 g/l.

- Specificity / Interferences**

Hemoglobin up to 0.32 g/dl, bilirubin up to 29,5 mg/dl, triglycerides up to 264 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean	SD	CV
	[g/l]	[g/l]	[%]
level 1	0.42	0.006	1.53
level 2	0.88	0.011	1.26
level 3	1.35	0.014	1.04
Reproducibility (day to day) n = 66	Mean	SD	CV
	[g/l]	[g/l]	[%]
level 1	0.33	0.001	3.44
level 2	0.95	0.003	3.04
level 3	1.51	0.004	2.71

- Method comparison**

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$$y = 0.8998x + 0.0867 \text{ g/l;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Marcovina, S.M., Albers, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A1 and B. III.
- Tietz, N.W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 146, (2006).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 1802, (1999).

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### PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

#### Parameters

No.	53	Prim.Wave.	340
Test	APO B	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	250
Unit	g/l	R2 Vol.	50
Decimals	3	Line. Limit	
Incubation	25	Antigen Check	v
Reaction	-1   25	Substrat	0
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower	0	Lower	0
Upper	0	Upper	0
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample Vol.	5	Full Name	APO B
Dilution	45	Print No.	53

#### Calibration

Rule	Logistic 5P
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0

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