ACCENT-300 HAPTOGLOBIN

Nr kat. 7-315 (PL)

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Cat. No 7-315 (EN)

INTENDED USE

Diagnostic kit for determination of haptoglobin concentration, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Haptoglobin is an acute phase protein whose primary function consists in binding free haemoglobin in serum. The complex is removed within minutes by the reticulo-endothelial system where its components are metabolised to free aminoacids and iron. Haptoglobin consequently plays a major role in preventing the loss of haemoglobin in urine and the consequent iron loss from the iron pool. The haptoglobin levels are increased during the acute phase and in such conditions as burns or nephrotic syndrome. Haptoglobin levels are abnormally high in intravascular hemolysis and when haemoglobin turnover is increased such as during haemolytic anaemia, transfusion reactions and malaria.

METHOD PRINCIPLE

The haptoglobin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to haptoglobin concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 40ml
 2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human haptoglobin antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found

to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, fresh serum is recommended. Samples may be stored up to 2 weeks at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

adults	0.26 – 1.85 g/l
newborns	0.05 – 0.48 g/l
> 60 year M	0.35 – 1.64 g/l
> 60 year F	0.40 – 1.75 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Hitachi. Results may vary if a different instrument or a manual procedure is used.

- LOQ (Limit of Quantitation):** 0.002 g/l

- Linearity:** up to 3.86 g/l.

- Interferences:**

Hemoglobin up to 0.06 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 1669.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.49	0.007	1.44
level 2	1.31	0.008	0.63
level 3	2.14	0.017	0.77

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.62	0.018	2.84
level 2	1.48	0.037	2.48
level 3	2.27	0.050	2.21

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$$y = 1.03x - 0.0 \text{ g/l};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

No.	41	Prim.Wave.	340
Test	HAPTO	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	250
Unit	g/l	R2 Vol.	50
Decimals	2	Line. Limit	
Incubation	25	Antigen Check	
Reaction	-1 25	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample Vol.	45	Full Name	HAPTO
Dilution	5	Print No.	41

Calibration

Rule	Spline
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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