## ACCENT-300 ALAT

Nr kat. 7-316 (PL)

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Cat. No **7-316** (EN)

### INTENDED USE


Diagnostic kit for determination of alanine aminotransferase activity intended to use in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Alanine aminotransferase (ALAT, ALT, GPT) is an enzyme participated in amino acids metabolism. ALAT is present in all tissues but the highest level is found in liver and kidney cells. Damage of hepatocytes or kidney cells causes significant release of ALAT into the circulation. Measurement of ALT activity in serum is valuable in the diagnosis of liver diseases: jaundice, mononucleosis or hepatic cirrhosis.

### METHOD PRINCIPLE

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at  $\lambda=340$  nm is directly proportional to alanine aminotransferase activity.

### REAGENTS

#### Package

1-Reagent 4 x 48 ml  
2-Reagent 1 x 48 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

Tris (pH 7.5) 100 mmol/l  
L-alanine 500 mmol/l  
LDH > 36.7  $\mu\text{kat/l}$   
2-oxoglutarate 15 mmol/l  
NADH 0.18 mmol/l  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, heparinized or EDTA plasma free from hemolysis. Hemolysis should be avoided, since ALAT activity in erythrocytes is 3 to 5 times higher than in normal serum.

Do not freeze the samples. ALAT activity remains stable in specimen up to 3 days at 15-25°C or up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

### REFERENCE VALUES <sup>6</sup>

serum / plasma	37°C	
women	up to 31 U/l	up to 0.517 $\mu\text{kat/l}$
men	up to 41 U/l	up to 0.683 $\mu\text{kat/l}$

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 8.2 U/l (0.137  $\mu\text{kat/l}$ ).
- Linearity:** up to 1100 U/l (18.3  $\mu\text{kat/l}$ ).
- Specificity / Interferences**

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	41.45	1.96	4.73
level 2	146.20	1.60	1.09
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	41.99	1.65	3.94
level 2	148.05	3.03	2.05

#### Method comparison

A comparison between ALAT values determined at ACCENT-300 (y) and at COBAS INTEGRA 400 (x) using 69 samples gave following results:

$$y = 1.1805x - 4.2315 \text{ U/l};$$

$$R = 0.9998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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