## ACCENT-300 CK

Nr kat. **7-320** (PL) ADP 2 mmol/l

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Cat. No **7-320** (EN)

### INTENDED USE

Diagnostic kit for determination of creatine kinase activity, used in automatic analyser ACCENT-300.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased level of CK is usually the result of muscle injury, myocardial or pulmonary infarction.

### METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC).

The rate of absorbance changes at  $\lambda=340$  nm is directly proportional to creatine kinase activity.

### REAGENTS

<b>Package</b>	
1-Reagent	2 x 40 ml
2-Reagent	1 x 17 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 10 weeks.

### Concentrations in the test

<b>1-Reagent</b>	
imidazole buffer	100 mmol/l
glucose	20 mmol/l
N-acetylcysteine	20 mmol/l
magnesium acetate	10 mmol/l
EDTA	2 mmol/l
NADP	2 mmol/l
ADP	2 mmol/l
AMP	5 mmol/l
HK	> 2.5 U/ml
<b>2-Reagent</b>	
diadenosinepentaphosphate	10 $\mu$ mol/l
glucose-6-phosphate-dehydrogenase (G6P-DH)	> 1.5 U/ml
creatine phosphate preservative	30 mmol/l
ACCENT-300 CK	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Do not use reagents past the expiry date.
- Do not interchange caps among reagents.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains imidazole.

### Danger

H360 May damage fertility or the unborn child.  
P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P308+P313 IF exposed or concerned:

Get medical advice/attention.

P405 Store locked up.

P501 Dispose of the contents/containers in accordance with the current legislation on waste treatment.

### SPECIMEN

Serum, free from hemolysis.  
CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light.  
Specimens can be stored up to 4-8 hours at 15-25 °C or 1-2 days at 2-8°C or 1 month at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
Deionized water is recommended as a reagent blank.

### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: CK – CREATININE, CK – MG, CK-MB – CK. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT\_300\_CARRYOVER.

### REFERENCE VALUES <sup>8</sup>

serum	37°C	
female	< 167 U/l	< 2.78 $\mu$ kat/l
male	< 190 U/l	< 3.17 $\mu$ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 10 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 31.1 U/l (0.52  $\mu$ kat/l).
- Linearity:** up to 2000 U/l (33.3  $\mu$ kat/l).

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor

#### Specificity / Interferences

Haemoglobin up to 0.156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	129.44	1.75	1.35
level 2	543.87	2.93	0.54
Reproducibility (day to day) n = 80	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	144.5	2.40	1.66
level 2	514.6	18.98	3.69

#### Method comparison

A comparison between CK values determined at **ACCENT-300** (y) and at **COBAS INTEGRA 400** (x) using 41 samples gave following results:  
 $y = 1.0204x - 1.5393$  U/l;  
 $R = 0.9995$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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