ACCENT-300 CREATININE

Nr kat. 7-333

(PL)

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Cat. No 7-333 (EN)

INTENDED USE

Diagnostic kit for determination of creatinine concentration, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

METHOD PRINCIPLE

Modified Jaffe's method, without deproteinization. In alkaline solution picrate reacts with creatinine to form a yellow-red 2,4,6-trinitrocylohexadienate. The colour intensity is proportional to the creatinine concentration.

REAGENTS

Package

1-Reagent 2 x 40 ml
2-Reagent 1 x 20.5 ml

The reagents when stored at 15-25°C are stable up to expiry date printed on the package. The reagents are stable for 8 weeks on board the analyser at 2-10°C.

Concentrations in the test

sodium hydroxide ≤ 450 mmol/l
carbonate buffer ≤ 150 mmol/l
picric acid ≤ 38,8 mmol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance of the working reagent is less than 0.225 (read against distilled water, wavelength $\lambda = 500$ nm).
- 1-Reagent and 2-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains sodium hydroxide.
2-REAGENT contains picric acid.

Danger

H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305 +P351 +P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

SPECIMEN

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine collected without preservatives.⁹

Urine preparation: Before analysis urine sample should be diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Mix well probes before measurement.

Specimen can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: CREATININE - BIL TOTAL, CREATININE - BIL DIRECT, CK - CREATININE, CK-MB - CREATININE, URINE PROTEINS - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007-ACCENT-300_CARRYOVER.

REFERENCE VALUES ⁷

serum / plasma	mg/dl	μmol/l
female	0.6 – 1.1	53 – 97
male	0.7 – 1.3	62 – 115
24-hours urine	mg/kg/24h	μmol/kg/24h
female	11 – 20	97 – 177
male	14 – 26	124 – 230

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every week, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.44 mg/dl (38.9 μmol/l).
- Linearity:** up to 20 mg/dl (1768 μmol/l).

If creatinine concentration exceeds 20 mg/dl, dilute the sample with 0.9% NaCl in the ratio of 1 to 2 and repeat the assay. Multiply the result by 3.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 500 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.58	0.07	4.27
level 2	3.71	0.04	1.04
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.26	0.11	8.92
level 2	3.92	0.16	4.00

Method comparison

A comparison between creatinine values determined at ACCENT-300 (y) and at COBAS INTEGRA 400 (x) using 29 samples gave following results:

$$y = 0.9661x + 0.0559 \text{ mg/dl};$$

$$R = 0.9915 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

No.	8	Prim.Wave.	510
Test	CREA	Sec.Wave.	700
Method	Fixed Time	Sample Vol.	25
Direction	Ascend	R1 Vol.	200
Unit	mg/dl	R2 Vol.	50
Decimals	2	Line. Limit	20
Incubation	3	Antigen Check	
Reaction	1 5	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	0.44
Upper	2.5	Upper	20
Sample Vol.	45	Full Name	Creatinine
Dilution	5	Print No.	8

Calibration

Rule	Two Point Linear
K Factor	0
Replicates	3
Interval	7
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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