## ACCENT-300 ALPHA 1-MICROGLOBULIN

Nr kat. 7-335 (PL)

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Cat. No 7-335 (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$  1-microglobulin concentration used in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$  1-microglobulin ( $\alpha$ Mi) is a low molecular weight glycoprotein (24-33 kD) which was initially isolated from the urine of patients with renal tubular disorders in 1975. It is mainly synthesized in the liver and is widely distributed in various body fluids.

The measurement of  $\alpha$ Mi in serum and urine has been considered to be useful for the diagnosis of functional renal disorders, the assessment of the progress and prognosis of diseases.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between  $\alpha$ Mi in a sample and anti- $\alpha$ Mi antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of  $\alpha$ Mi in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

#### Package

1-Reagent 1 x 32 ml  
2-Reagent 1 x 32 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-10°C. The reagents are stable for 6 weeks on board the analyser at 2-10°C.

#### Concentrations in the test

suspension of latex particles sensitized with (rabbit) anti- $\alpha$ Mi antibodies 0.25 w/v%  
glycine buffer solution  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma, urine.  
If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless, it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES <sup>3</sup>

serum, plasma	10.0 – 30.0 mg/l
urine	1.0 – 5.0 mg/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use control serum for determination of  $\alpha$ Mi with each batch of samples, eg.: ROCHE or BIORAD.

For the calibration of automatic analysers systems the CORMAY ALPHA 1-MGLOB CALIBRATORS (S) (Cat. No 4-286) for serum samples and the CORMAY ALPHA 1-MGLOB CALIBRATORS (U) (Cat. No 4-285) for urine samples is recommended.

The calibration curve should be prepared every time the test is performed, with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analyser HITACHI 917. Results may vary if a different instrument is used.

#### Analytical range:

1.5 – 200 mg/l (serum)  
0.3 – 50.0 mg/l (urine)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.5 g/dl, NH<sub>4</sub>Cl do 400 mg/dl, bilirubin up to 31 mg/dl do not interfere with the test in urine.

#### Precision

Repeatability (run to run) n = 20	Mean [mg/l]	SD [mg/l]	CV [%]
level 1	0.5	0.0	3.97
level 2	1.6	0.0	1.81
level 3	13.9	0.1	0.42

#### Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 55 serum samples gave following results:  
 $y = 1.00 x + 2.83$  mg/l;  
R = 1.00 (R – correlation coefficient)

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 55 urine samples gave following results:

$y = 1.00 x - 0.52$  mg/l;  
R = 1.00 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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#### PROGRAM NA ANALIZATOR – surowica

##### Parameters

No.	51	Prim.Wave.	578
Test	AMGs	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	200
Unit	mg/l	R2 Vol.	200
Decimals	1	Line. Limit	
Incubation	5	Antigen Check	
Reaction	5   23	Substrat	0
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower	0	Lower	0
Upper	0	Upper	0
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample Vol.	45	Full Name	AMGs
Dilution	5	Print No.	51

##### Calibration

Rule	Logistic 5P
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0

#### PROGRAM NA ANALIZATOR – mocz

##### Parameters

No.	52	Prim.Wave.	578
Test	AMGu	Sec.Wave.	700
Method	Endpoint	Sample Vol.	9
Direction	Ascend	R1 Vol.	180
Unit	mg/l	R2 Vol.	180
Decimals	2	Line. Limit	
Incubation	5	Antigen Check	
Reaction	5   23	Substrat	0
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower	0	Lower	0
Upper	0	Upper	0
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample Vol.	45	Full Name	AMGu
Dilution	5	Print No.	52

##### Calibration

Rule	Logistic 5P
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0