## ACCENT-300 BIL TOTAL

Nr kat. 7-345 (PL)

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### LITERATURA

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

Data wydania: 10.2020

## ACCENT-300 BIL TOTAL

Cat. No **7-345** (EN)

### INTENDED USE

Diagnostic kit for determination of total bilirubin concentration, intended to use in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Bilirubin is a yellow pigment – product of heme degradation. For clinical purposes, bilirubin is expressed as two fractions: conjugated and unconjugated. In hepatocytes bilirubin is enzymatically conjugated with glucuronic acid residues. This form is called direct or conjugated. Bilirubin without glucuronic acid modification is bound to albumin and is termed unconjugated or indirect. Indirect bilirubin is calculated as the difference between total and direct bilirubin.

Serum bilirubin measurement is widely used as a screening test for liver functions. Hiperbilirubinemia is usually the result of jaundice (mechanical, hemolytic), Dubin-Jonson syndrome, Gilbert's syndrome, Crigler-Najjar syndrome, bile ducts disease.

### METHOD PRINCIPLE

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in a acidic solution, total bilirubin (both conjugated – direct, and unconjugated bilirubin) is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the total bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

### REAGENTS

#### Package

1-Reagent 5 x 40 ml  
 2-Reagent 2 x 25.5 ml

The reagents when stored at 10-25°C are stable up to expiry date printed on the package. The reagents are stable for 7 weeks on board the analyser at 2-10°C.

### Concentrations in the test

#### 1-Reagent

citrate buffer (pH 2.8) 90 mmol/l  
 detergent

#### 2-Reagent

phosphate buffer (pH 7.0) 4.6 mmol/l  
 sodium metavanadate 3.0 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Do not freeze reagents.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum free from haemolysis.  
 Serum should be separated from red blood cells as soon as possible after blood collection. Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended. It is recommended to follow CLSI procedures regarding specimen collecting and handling.  
 Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark at 2-8°C, at the most 3 days.  
 Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 For reagent blank deionized water is recommended.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: CREATININE – BIL TOTAL. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

Age	mg/dl	µmol/l
0 - 1 days	< 8	< 137
1 - 2 days	< 12	< 205
3 - 5 days	< 16	< 274
5 days - 60 years	0.3 - 1.2	5 - 21
60 - 90 years	0.2 - 1.1	3 - 19
> 90 years	0.2 - 0.9	3 - 15

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176), LEVEL 2 (Cat. No 5-175; 5-177) are recommended. As a 0 calibrator deionised water should be used.

The calibration curve should be prepared every 7 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-300 and Hitachi 911. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.15 mg/dl (2.57 µmol/l).

- Linearity:** up to 43 mg/dl (735 µmol/l).

- Specificity / Interferences**

Haemoglobin up to 0.25 g/dl, ascorbic acid up to 500 mg/l and triglycerides up to 250 mg/dl do not interfere with the test.

- Precision**

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[µg/dl]	[%]
level 1	0.93	0.01	1.23
level 2	4.11	0.01	0.35
Reproducibility (day to day) n = 80	Mean	SD	CV
	[µmol/l]	[µmol/l]	[%]
level 1	25.2	0.99	3.9
level 2	77.0	0.91	1.2
level 3	79.7	1.10	1.4

- Method comparison**

A comparison between total bilirubin values determined at ACCENT-300 (y) and at Olympus AU400 (x) using 39 samples gave following results:

$$y = 1.0361 x - 0.0217 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2);116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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1. Tokuda K, Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
2. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
3. Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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#### Calibration

Rule	Multipoint Linear
K Factor	
Replicates	3
Interval	49
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0