ACCENT-300 BILE ACIDS

Cat. No 7-361 (EN)

INTENDED USE

Diagnostic kit for determination of total bile acids concentration, used in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.


INTRODUCTION

Bile acids are the main product of degradation of endogenous cholesterol formed in the liver. Total bile acids are metabolized in the liver and are a valuable indicator of normal or abnormal liver function. Serum total bile acids are increased in patients with viral hepatitis, liver cirrhosis and liver cancer.

METHOD PRINCIPLE

Enzymatic method with 3- α -hydroxysteroid dehydrogenase (3- α -HSD).

Bile acids under the influence of 3-hydroxysteroid dehydrogenase (3- α -HSD) in the presence of thio-NAD are converted to 3-ketosteroids and thio-NADH. The reaction is reversible and 3- α -HSD can convert 3-ketosteroids and NADH to bile acids and NAD.

The rate of thio-NADH formation can be monitored at 405 nm and is proportional to the bile acids activity.

REAGENTS

Package
 1-Reagent 2 x 16 ml
 2-Reagent 2 x 4 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

Concentrations in the test

1-Reagent
 Thio-NAD > 0.1 mmol
 Buffer
2-Reagent
 3- α -HSD > 2 kU/l
 NADH > 0.1 mmol
 Buffer

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Yellow or yellow-brown color of the reagent does not affect the reagents performance.

- Reagents from different lots must not be interchanged.
- Samples from patients treated with ursodeoxycholic acid (UDCA) are not suitable for the determination of total bile acid concentrations.

SPECIMEN

Serum.
 Total bile acids concentration is increased after meals, therefore samples should be collected under fasting conditions. Serum samples are stable for a 7 days at 4 °C or for 3 month at -20 °C.
 Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank deionized water is recommended.

REFERENCE VALUES³

serum	2.5 – 6.8 $\mu\text{mol/l}$ (1.25 – 3.4 $\mu\text{g/ml}$)
-------	--

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY BILE ACIDS CONTROLS (Cat. No 5-149).

For the calibration of automatic analysers systems the CORMAY BILE ACIDS CALIBRATOR (Cat. No 3-125) is recommended. As a 0 calibrator deionised water should be used.

The calibration curve should be prepared every 2 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-300 and Hitachi 717. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 2.0 $\mu\text{mol/l}$ (1.0 $\mu\text{g/ml}$).
- Linearity:** up to 150 $\mu\text{mol/l}$ (75 $\mu\text{g/ml}$).

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 50 mg/dl, ascorbic acid up to 50 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean [$\mu\text{mol/l}$]	SD [$\mu\text{mol/l}$]	CV [%]
level 1	31.00	0.71	2.30
level 2	107.49	1.18	1.10

Reproducibility (day to day) n = 20	Mean [$\mu\text{mol/l}$]	SD [$\mu\text{mol/l}$]	CV [%]
level 1	8.12	0.24	2.90
level 2	23.0	0.61	2.60

Method comparison

A comparison between bile acids values determined at ACCENT-300 (y) and at Hitachi 912 (x) using 29 samples gave following results:

$$y = 1.0005x - 1.1559 \mu\text{mol/l};$$

$$R = 0.9994 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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