A-400 GLUCOSE

Nr kat. 7-401 (PL)

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LITERATURA

- Barham P., Trinder P.: An improved colour reagent for the determination of blood glucose by the oxidase system; Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2017, Diabetologia Praktyczna, tom 3, suplement A, 2017.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

Data wydania: 10.2020.

A-400 GLUCOSE

Cat. No **7-401** (EN)

INTENDED USE

Diagnostic kit for determination of glucose concentration, may be used in automatic analysers BS-400 and BS-480.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Glucose is a simple six-carbon sugar. Oxidative metabolism of glucose provides the energy for most cellular processes. Glucose level in the blood is tightly controlled by several hormones. Elevated glucose level is the classic sign of diabetes mellitus. Glucose level abnormalities (hyper- or hypoglycemia) might be caused also by pancreas tumors and diseases of liver, thyroid gland or adrenal glands.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glucose oxidase.

The colour intensity is proportional to the glucose concentration.

REAGENTS

Package
1-Reagent 4 x 40 ml

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

Concentrations in the test

phosphate buffer (pH 7.0)	< 240 mmol/l
phenol	< 6 mmol/l
glucose oxidase (GOD)	< 480 µkat/l
peroxidase (POD)	< 44 µkat/l
4-aminoantipyrine (4-AA)	< 0,9 mmol/l
stabilizers and preservatives	

Warnings and notes

- Protect from direct sunlight and avoid contamination!

SPECIMEN

EDTA or heparinized plasma in tubes containing sodium fluoride or sodium iodoacetate additive; serum, free from hemolysis; cerebrospinal fluid.

Plasma / Serum. Serum and plasma specimens should be separated from cells within 30 minutes after collection. Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level. Serum and plasma can be stored up to 2 days at 4°C.³ Plasma is the specimen recommended for the glucose determination in the blood.⁵

Cerebrospinal fluid. Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.⁴

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

For reagent blank deionized water is recommended.

Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE – CREA ENZYMATYC, GLUCOSE – PHOSPHORUS II GEN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER.

REFERENCE VALUES

	mg/dl	mmol/l
plasma, serum ^{5,6,7}	70 – 99	3.9 – 5.5
cerebrospinal fluid ⁸	40 – 70	2.2 – 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analyser system BS-400 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

For the calibration of automatic analyser system BS-480 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

5.56 mg/dl (0.31 mmol/l) – BS-400

4.42 mg/dl (0.24 mmol/l) – BS-480

Linearity:

up to 500 mg/dl (27.5 mmol/l) – BS-400

up to 530 mg/dl (29.2 mmol/l) – BS-480

For higher concentrations, dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 2.50 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=20	level 1	88.91	0.26	0.29
	level 2	337.08	1.41	0.42
BS-480 n=10	level 1	90.58	0.28	0.31
	level 2	284.45	1.19	0.42
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=56	level 1	107.19	1.00	0.93
	level 2	274.15	4.16	1.52
BS-480 n=20	level 1	89.57	1.16	1.30
	level 2	276.99	2.20	0.79

Method comparison

A comparison between glucose values determined at **BS-400** (y) and at **OLYMPUS AU400** (x) using 40 samples gave following results:

$$y = 0.9574 x + 8.3695 \text{ mg/dl;}$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **BS-480** (y) and at **COBAS INTEGRA 400 Plus** (x) using 40 samples gave following results:

$$y = 0.9964 x - 0.0475 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2017, Diabetologia Kliniczna, tom 3, suplement A, 2017.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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- Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2017, Diabetologia Kliniczna, tom 3, suplement A, 2017.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACCC Press (2000).

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• Calibration

Calibration		Judgment Criteria		Blank Abs.	
Rule	Two-point Linear	Sensitivity		Error Limit	
Replicate	3	Factor Diff.		Corr. Coeff.	
K		SD			
• QC				Auto QC	
Rules		Cum. Sum Check		Interval	
Westgard Multi-rule		1.0 - 2.7			
v	1-2S	v	4-1S		
v	1-3S	v	10-X		
v	2-2S	v			
		• 1.0 - 3.0			
		0.5 - 5.1			

• BS-480

Chem	GLUC	No.	001	Sample Type	SERUM/CSF
Chemistry	GLUCOSE	Print name	GLUC		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	505	Sec Wave	800		
Unit	mg/dl	Decimal	0.1		
Blank Time	11 12	Reaction Time	63 65		
Standard	Sample Vol 2.5 μL	Aspirated 20 μL	Diluent 180 μL	Reagent Vol 200 μL	Diluent μL
Decreased	2.5 μL	20 μL	180 μL	R2	μL
Increased	μL	μL	μL	R3	μL
	μL	μL	μL	R4	μL
	Sample Blank	v	Auto Retun		
Linearity Range (Standard)	4.42	530	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-33000	33000
R1 Blank Abs	-33000	33000	Uncapping Time	77	Day(s)
Blank Response	-33000	33000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension		
	Prozone Check	o	Rate Check		
Q1	0	Q2	0	• Antigen Addition	
PC	0	ABS	0	Q3	0
				Q4	0

Calibration Settings

Math Model	Multi-point Linear	Auto Calibration	
Factor		<input type="checkbox"/>	Bottle Changed
Replicates	3	<input type="checkbox"/>	Lot Changed
		<input type="checkbox"/>	Cal Time
Acceptance Limits			
Cal Time	1848	Hour	
Slope Diff		SD	
Sensitivity		Repeatability	
Deter Coeff			

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