## A-400 UREA

Nr kat. **7-406** (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of urea concentration, intended to use in automatic analysers BS-400 and BS-480.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Urea is a product of amino acids catabolism. It is produced in liver and excreted in urine. Urea in the blood is reported as the blood urea nitrogen (BUN). Increased urea concentration in the serum, called uremia, is observed due to dehydration, renal failure, high-protein diet, increased protein catabolism caused by tissue injury or massive bleeding into the alimentary tract. The reason of reduced urea level could be overhydration, low-protein diet or starvation and severe liver disease.

### METHOD PRINCIPLE

Kinetic, enzymatic method with urease and glutamate dehydrogenase.

The rate of absorbance changing at  $\lambda=340$  nm is proportional to the urea concentration.

### REAGENTS

#### Package

1-Reagent 3 x 40 ml  
2-Reagent 1 x 31 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

Tris (pH 7.8) 96 mmol/l  
ADP 0.6 mmol/l  
urease 266.7  $\mu\text{kat/l}$   
GLDH 16  $\mu\text{kat/l}$   
NADH 0.26 mmol/l  
2-oxoglutarate 9 mmol/l  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN<sup>9,10,11</sup>

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine.

Do not use heparine ammonium salt and fluoride as anticoagulants.

Specimen can be stored up to 7 days at 2-8°C.

**Urine preparation:** Samples with visible turbidity or the presence of precipitates should be pre-centrifuged.

Urine sample should be mixed well before analysis, diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Bacterial growth in the specimen may cause erroneously elevated results.

24-hours urine samples should be adjusted to pH < 7 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionized water is recommended as a reagent blank.

### REFERENCE VALUES<sup>8</sup>

serum / plasma	mg/dl	mmol/l
	< 50	< 8.3
24-hours urine	g/24h	mmol/24h
	20 – 35	300 – 550

1 mg of urea corresponds to 0.467 mg of urea nitrogen (BUN)  
It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

**For the calibration of automatic analyser system BS-400** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

**For the calibration of automatic analyser system BS-480** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

### Sensitivity:

4.5 mg/dl (0.749 mmol/l) – BS-400  
3.3 mg/dl (0.549 mmol/l) – BS-480

### Linearity:

up to 350 mg/dl (58.3 mmol/l) – BS-400  
up to 430 mg/dl (71.6 mmol/l) – BS-480

### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=10	level 1	31.36	0.53	1.69
	level 2	101.59	0.62	0.61
BS-480 n=10	level 1	32.29	0.53	1.65
	level 2	102.03	0.86	0.84
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=20	level 1	34.95	0.35	0.99
	level 2	97.78	1.24	1.27
BS-480 n=20	level 1	32.70	0.73	2.24
	level 2	100.10	2.06	2.06

### Method comparison

A comparison between urea concentration determined at **BS-400** (y) and at **ADVIA 1650** (x) using 25 samples gave following results:

$$y = 0.953x + 0.4735 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between urea concentration determined at **BS-480** (y) and at **ADVIA 1650** (x) using 57 samples gave following results:

$$y = 0.9236x + 1.1199 \text{ mg/dl};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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• **Calibration**

<b>Calibration</b>	<b>Judgment Criteria</b>	Blank Abs.	<input type="text"/>
Rule <input type="text" value="Two-point Linear"/>	Sensitivity <input type="text"/>	Error Limit	<input type="text"/>
Replicate <input type="text" value="3"/>	Factor Diff. <input type="text"/>	Corr. Coeff.	<input type="text"/>
K <input type="text"/>	SD <input type="text"/>		
<b>• QC</b>			
<b>Rules</b>		<b>Auto QC</b>	
Westgard Multi-rule	Cum. Sum Check <input type="checkbox"/>	Interval <input type="text"/>	
<input type="checkbox"/> 1-2S	<input type="checkbox"/> R-4S		
<input type="checkbox"/> 1-3S	<input type="checkbox"/> 4-1S		
<input type="checkbox"/> 2-2S	<input type="checkbox"/> 10-X		
	• 1.0 - 2.7		
	• 1.0 - 3.0		
	• 0.5 - 5.1		

• **BS-480**

Chem <input type="text" value="UREA"/>	No. <input type="text" value="010"/>	Sample Type <input type="text" value="SERUM/URINE"/>
Chemistry <input type="text" value="UREA"/>	Print name <input type="text" value="UREA"/>	
Reaction Type <input type="text" value="Kinetic"/>	Reaction Direction <input type="text" value="Decrease"/>	
Pri Wave <input type="text" value="340"/>	Sec Wave <input type="text" value="450"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text" value="38"/> <input type="text" value="48"/>	Reaction Time <input type="text" value="55"/> <input type="text" value="60"/>	
Standard <input type="text" value="2"/> <input type="text" value="2"/> <input type="text" value="2"/>	Sample Vol <input type="text" value="2"/> <input type="text" value="2"/> <input type="text" value="2"/>	Aspirated <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/>
Decreased <input type="text" value="2"/> <input type="text" value="2"/> <input type="text" value="2"/>	Diluent <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/>	R1 <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/>
Increased <input type="text" value="2"/> <input type="text" value="2"/> <input type="text" value="2"/>	R2 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	R2 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun	R3 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	R3 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
	R4 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	R4 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
Linearity Range (Standard) <input type="text" value="3.3"/> <input type="text" value="430"/>	Linearity Limit <input type="text" value="0.2"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	
R1 Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time <input type="text" value="84"/> Day(s)	
Blank Response <input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	• Antigen Addition <input type="checkbox"/>	
Q1 <input type="text" value="0"/> Q2 <input type="text" value="0"/>	Q3 <input type="text" value="0"/> Q4 <input type="text" value="0"/>	
PC <input type="text" value="0"/> ABS <input type="text" value="0"/>		

<b>Calibration Settings</b>	<b>Auto Calibration</b>
Math Model <input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="3"/>	<input type="checkbox"/> Cal Time
<b>Acceptance Limits</b>	
Cal Time <input type="text" value="840"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text"/>
Deter Coeff <input type="text"/>	

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