A-400 UA

Nr kat. 7-408 (PL)

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Cat. No **7-408** (EN)

INTENDED USE



Diagnostic kit for determination of uric acid concentration, used in automatic analysers BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

REAGENTS

Package

1-Reagent 2 x 40 ml
 2-Reagent 1 x 21.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

Concentrations in the test

buffer PIPES (pH 7.0) 100 mmol/l
 4-aminoantipyrine 0.78 mmol/l
 ADPS 0.67 mmol/l
 ferricyanide potassium 3.8 µmol/l
 peroxidase (POD) > 38.34 µkat/l
 uricase > 1.65 µkat/l
 preservative

Warnings and notes

- Protect from direct sunlight and contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning.

H315 Causes skin irritation.
 H319 Causes serious eye irritation.
 P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum, heparinized plasma free from hemolysis. Do not use EDTA and fluoride as anticoagulants. Specimen can be stored 3-5 days at 2-8°C or 6 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use. For reagent blank deionized water is recommended.

Actions required:

When performing assays in the analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT II GEN – UA, URINE PROTEINS II GEN – UA. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER

REFERENCE VALUES ⁵

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. As a 0 calibrator deionized water should be used. The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:**
 0.15 mg/dl (8.92 µmol/l) – BS-400
 0.08 mg/dl (4.76 µmol/l) – BS-480
- Linearity:**
 up to 50 mg/dl (2974 µmol/l) – BS-400
 up to 32 mg/dl (1903 µmol/l) – BS-480

Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 31 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=10	level 1	4.70	0.03	0.66
	level 2	9.92	0.06	0.61
BS-480 n=10	level 1	4.82	0.02	0.48
	level 2	9.58	0.03	0.30
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=56	level 1	5.61	0.08	1.51
	level 2	9.48	0.52	5.43
BS-480 n=20	level 1	4.75	0.03	0.62
	level 2	9.47	0.06	0.61

Method comparison

A comparison between uric acid concentration at **BS-400** (y) and at **Cobas Integra 400 Plus** (x) using 82 samples gave following results:
 $y = 0.9414x + 0.2702$ mg/dl;
 $R = 0.996$ (R – correlation coefficient)

A comparison between uric acid concentration at **BS-480** (y) and at **Cobas Integra 400 Plus** (x) using 45 samples gave following results:
 $y = 1.005x + 0.0548$ mg/dl;
 $R = 0.992$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATORY / APPLICATION for / АДАПТАЦИЯ для:

• **BS-400**

• Basic	
Test information	Reagent Volume
No. <input type="text" value="9"/>	R1 <input type="text" value="160"/>
Test <input type="text" value="UA"/>	R2 <input type="text" value="40"/>
Full Name <input type="text" value="Uric Acid"/>	R3 <input type="text"/>
Std. No. <input type="text" value="9"/>	R4 <input type="text"/>
Reaction Parameters	Sample Volume
Reac. Type <input type="text" value="Endpoint"/>	Standard <input type="text" value="4"/> <input type="text" value="15"/> <input type="text" value="10"/>
Pri. Wave <input type="text" value="546"/>	Increased <input type="text" value="8"/> <input type="text" value="15"/> <input type="text" value="10"/>
Sec. Wave <input type="text" value="800"/>	Decreased <input type="text" value="2"/> <input type="text" value="15"/> <input type="text" value="10"/>
Judgment Criteria	Result Setup
Absorbance <input type="text" value="0"/> <input type="text" value="0"/>	Decimal <input type="text" value="0.01"/> Slope <input type="text" value="1"/>
Incr. Test <input type="text" value="0"/>	Unit <input type="text" value="mg/dl"/> Inter <input type="text" value="0"/>
Decre. Test <input type="text" value="0"/>	
Lin. Range <input type="text" value="0.15"/> <input type="text" value="50"/>	
Lin. Limit <input type="text"/>	
Subs. Limit <input type="text"/>	
Q1 <input type="text" value="0"/>	<input type="checkbox"/> Prozone <input type="checkbox"/> Rate <input type="text" value="0"/>
PC <input type="text" value="0"/>	Q2 <input type="text" value="0"/> Q3 <input type="text" value="0"/> Q4 <input type="text" value="0"/>
	ABS <input type="text" value="0"/>
• Calibration	Judgment Criteria
Rule <input type="text" value="Multi-point Linear"/>	Sensitivity <input type="text"/>
Replicate <input type="text" value="3"/>	Factor Diff. <input type="text"/>
K <input type="text"/>	SD <input type="text"/>
• QC	Blank Abs. <input type="text"/>
Rules	Error Limit <input type="text"/>
Westgard Multi-rule	Corr. Coeff. <input type="text"/>
<input type="checkbox"/> 1-2S <input type="checkbox"/> R-4S	Cum. Sum Check <input type="text"/>
<input type="checkbox"/> 1-3S <input type="checkbox"/> 4-1S	<input type="checkbox"/> 1.0 - 2.7
<input type="checkbox"/> 2-2S <input type="checkbox"/> 10-X	<input type="checkbox"/> 1.0 - 3.0
	<input type="checkbox"/> 0.5 - 5.1
	Auto QC
	Interval <input type="text"/>

• **BS-480**

Chem <input type="text" value="UA"/>	No. <input type="text" value="009"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="Uric Acid"/>	Print name <input type="text" value="UA"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="Increase"/>	
Pri Wave <input type="text" value="546"/>	Sec Wave <input type="text" value="800"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="48"/> <input type="text" value="49"/>	Reaction Time <input type="text" value="71"/> <input type="text" value="73"/>	
Standard <input type="text" value="4"/> <input type="text" value="4"/> <input type="text" value="20"/> <input type="text" value="180"/>	Reagent Vol <input type="text" value="160"/> <input type="text" value="40"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="4"/> <input type="text" value="20"/> <input type="text" value="180"/>	R1 <input type="text" value="160"/> <input type="text" value="40"/>	R2 <input type="text" value="40"/>
Increased <input type="text" value="4"/> <input type="text" value="20"/> <input type="text" value="180"/>	R3 <input type="text" value="40"/>	R4 <input type="text" value="40"/>
<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun		
Linearity Range (Standard) <input type="text" value="0.08"/> <input type="text" value="32"/>	Linearity Limit <input type="text"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	
R1 Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time <input type="text" value="77"/> Day(s)	
Blank Response <input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	• Antigen Addition	
O1 <input type="text" value="0"/> O2 <input type="text" value="0"/>	O3 <input type="text" value="0"/> O4 <input type="text" value="0"/>	
PC <input type="text" value="0"/> ABS <input type="text" value="0"/>		

Calibration Settings	Auto Calibration
Math Model <input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	Replicates <input type="text" value="3"/>
	<input type="checkbox"/> Lot Changed
	<input type="checkbox"/> Cal Time
Acceptance Limits	
Cal Time <input type="text" value="1848"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text"/>
Deter Coeff <input type="text"/>	

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).